



IMAGING AND ANALYSIS SOFTWARE

Getting Started with ChemiDoc™ and ChemiDoc MP Sample Trays

Quick Guide

The following trays, accessories, and applications are supported for use with the imaging systems. The Blot/UV/Stain-Free Tray ships with the product; all other trays are optional accessories.

Tray Type or Accessory	Description
Blot/UV/Stain-Free Tray (catalog #12003028)	<p>Use with the following applications:</p> <ul style="list-style-type: none">Fluorescence blotsChemiluminescenceStain-free gels and blotsRhodamineCoomassie blue stainSYPRO StainsEthidium bromide (EtBr)SYBR® StainsOriole™ Fluorescent Gel StainGelRed <p>Caution: You must use the accessory ChemiDoc UV Shield with the Blot/UV/Stain-Free Tray during any band excision procedure. This UV shield reduces exposure to UV rays.</p> <p>Caution: Before you perform band excision, you must protect yourself from UV radiation by wearing protective gear, including UV protective safety glasses, a mask, lab coat, and gloves.</p>
White Tray (catalog #12003026)	<p>Use with protein applications that call for colorimetric stains:</p> <ul style="list-style-type: none">Coomassie blue stainCopper stainSilver stainZinc stainFast Blast™ DNA Stain
Blue Tray (catalog #12003027)	<p>Use with all nucleic acid applications that use blue excitable safe stains:</p> <ul style="list-style-type: none">SYBR® GreenSYBR® SafeSYBR® GoldGelGreenFluoresceinGelStar <p>Note: The Blue Tray prevents UV damage to nucleic acids during the imaging process and protects the user from exposure to UV rays during band excision. View Goggles (catalog #10042285) are included.</p>
ChemiDoc UV Shield (catalog #12003914)	Must be used with the Blot/UV/Stain-Free Tray when performing band excision.

Caution: Chemiluminescent substrates are caustic and will etch all forms of glass over time. Use plastic film or sheet protectors to avoid applying the substrate directly to the glass on the tray. Wash trays immediately after use if substrate was applied.

Caution: Do not clean the trays in a dishwasher; doing so can damage the trays. Avoid using hot water.

Using the Trays

1 Clean the Tray

Clean the tray with standard laboratory detergents or mild solvents, such as ethanol or methanol. Use lint free tissue

to wipe it dry. Dust particles or lint on a sample tray can glow under UV illumination.

2 Dry the Tray

Dust particles or lint on a sample tray can glow under UV illumination.

Wipe the tray dry using a lint free tissue.

3 Turn on the System and Insert the Tray

Turn on the imager. Insert the tray into the imager and close the door.

BIO-RAD

Coomassie is a trademark of BASF Aktiengesellschaft. GelRed and GelGreen are trademarks of Biotium, Inc. SYBR® and SYPRO are trademarks of Life Technologies Corporation. GelStar is a trademark of FMC Corporation.



**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site bio-rad.com **USA** 1 800 424 6723 **Australie** 61 2 9914 2800 **Austrie** 43 1 877 89 01 177 **Belgium** 32 (0)3 710 53 00 **Brazil** 55 11 3065 7550
Canada 1 905 364 3435 **China** 86 21 6169 8500 **Czech Republic** 420 241 430 532 **Denmark** 45 44 52 10 00 **Finland** 358 09 804 22 00
France 33 01 47 95 69 65 **Germany** 49 89 31 884 0 **Hong Kong** 852 2789 3300 **Hungary** 36 1 459 6100 **India** 91 124 4029300
Israel 972 03 963 6050 **Italy** 39 02 216091 **Japan** 81 3 6361 7000 **Korea** 82 2 3473 4460 **Mexico** 52 555 488 7670 **The Netherlands** 31 (0)318 540 666
New Zealand 64 9 415 2280 **Norway** 47 23 38 41 30 **Poland** 48 22 331 99 99 **Portugal** 351 21 472 7700 **Russia** 7 495 721 14 04
Singapore 65 6415 3188 **South Africa** 27 (0) 861 246 723 **Spain** 34 91 590 5200 **Sweden** 46 08 555 12700 **Switzerland** 41 026674 55 05
Taiwan 886 2 2578 7189 **Thailand** 66 662 651 8311 **United Arab Emirates** 971 4 8187300 **United Kingdom** 44 020 8328 2000

