

Notice

The information contained in this document is subject to change without notice.

APPLIED BIOSYSTEMS MAKES NO WARRANTY OF ANY KIND WITH REGARD TO THIS MATERIAL, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. Perkin-Elmer shall not be liable for errors contained herein or for incidental consequential damages in connection with the furnishing, performance, or use of this material.

Copyright Information

Reproduction or publication of this document in any form or format is prohibited without written permission of The PE Corporation

Copyright © 1994 The PE Corporation

All rights reserved.

Printed in the United States of America.

Software programs are protected by copyright. It is unlawful to duplicate these programs in any manner other than for the user's personal use. Specifically, it is unlawful to use or duplicate these programs other than for use with the purchaser's original computer.

Release Information

Manual Part No.	Manual Version
0993-8984	B

Trademarks

Applied Biosystems is a registered trademark of the PE Corporation

Epson is a trademark of Seiko Epson Corporation.

AmpliTaq, AmpliType, GeneAmp, GeneAmplimer, AmpliWax, EnviroAmp, MicroAmp, and Twelve Paq are registered trademarks/trademarks of Roche Molecular Systems, Inc. and F. Hoffman-LaRoche A.G.

The GeneAmp PCR Process, Taq DNA Polymerase, and AmpliTaq DNA Polymerase are covered by U.S. patents owned by Roche Molecular Systems, Inc., and issued and pending patents owned by F. Hoffman-La Roche A.G.

Registered trademarks, etc. used in this document, even when not specifically marked as such, are protected by law.




DNA Thermal Cycler 480

Quick Reference Guide

There are three main functions of the software that you can access after you turn on the DNA Thermal Cycler 480.


1 Access a Default File

Access a default file to customize, store, or run a PCR file and method.




2 Get to the File Directory

Review a directory of all or some of your PCR files.



3 Get to the Configuration and Diagnostic Files


Perform two configuration tasks and run seven diagnostic tests.



The Incubation Files


Soak File #1

Soak files hold a single temperature until you abort the run.



Time Delay File #2


Time Delay files hold a single temperature for a specified period of time.



The Cycle Files


Thermo-Cycle File #3

Thermo-Cycle files contain programmed thermal transitions and thermal plateaus for PCR cycling. You can also select the Auto Segment Extension feature.



Step-Cycle File #4

Step-Cycle files contain programmed thermal plateaus for PCR cycling. You can also select the Auto Segment Extension feature.



How to Customize a Soak File, File #1

SOAK FILE STEP to edit file #1	<i>This is the first display in the Time Delay file. Press STEP to edit the first parameter.</i>
Soak temperature (-5 to 100) 25 ° C	<i>Enter the desired soak temperature.</i> <ul style="list-style-type: none">• Press STEP to accept a default value and move to the next display. Or, to change a value and move to the next display, use the numeric keys then press ENTER.
End of file RUN-STORE-PRINT	

How to Customize a Time Delay File, File #2

TIME DELAY FILE STEP to edit file #2	<i>This is the first display in the Time Delay file. Press STEP to edit the first parameter.</i>
Delay Time (0-999) 7min 0sec	<i>Enter the desired delay time and temperature.</i> <ul style="list-style-type: none">• Press STEP to accept a default value and move to the next display or parameter. Or, to change a value and move to the next display or parameter, use the numeric keys then press ENTER.• Press HELP to move to the previous parameter.• Press STOP to move to the last display in the file.
Delay Temperature (-5 to 100) 72 ° C	
Beep while Delay NO-YES	<i>Choose YES to turn on a beeper, which will sound while the file is running, by pressing NO, then ENTER. Or, press ENTER to choose NO.</i>
Link to stored file (0-99, 0=SHUT-OFF) 0	<i>Enter the file number of the file you wish to link to in a method, then ENTER, or press STEP to accept the default link to shut-off.</i>
End of file RUN-STORE-PRINT	<i>Press NO to select one of the options, then ENTER. Press STOP to return to the first display in the file.</i>

How to Customize a Thermo-Cycle File, File #3

THERMO-CYCLE FILE
STEP to edit file #3

This is the first display in the Thermo-Cycle file.
Press **STEP** to edit the first parameter.

Seg 1 Target Temp
(-5 to 100) 94°C

Seg 1 Segment Time
(0-999) 0min 1sec

Enter the desired temperatures and times for Segments #1, 2, 3, 4, 5, and 6. Odd-numbered segments are thermal transitions and even-numbered segments are thermal plateaus.

- Press **STEP** to accept a default value and move to the next display or parameter.
Or, to change a value and move to the next display or parameter, use the numeric keys then press **ENTER**.
- Press **HELP** to move to the previous parameter.
- Press **STOP** to move to the last display in the file.



Seg 7 Target Temp
(-5 to 100) 0°C

Seg 7 Segment Time
(0-999) 0min 0sec

Segment #7 is the terminating segment in the file and must have a time and temperature of zero. The terminating segment marks the end of one cycle.

- Press **STEP** to accept a default value and move to the next display or parameter.
Or, to change a value and move to the next display or parameter, use the numeric keys then press **ENTER**.
- Press **HELP** move to the previous parameter

Auto Seg. Extension
NO-YES

Choose YES to automatically extend any segment by pressing NO, then **ENTER**. Then enter the number of the segment you wish to extend (usually 6) and the extension time per cycle in seconds; Or, press **ENTER** to choose NO.

Cycle Count
(1-99) 25

Enter the number of cycles you wish to run. Use the numeric keys and **ENTER** to change the value, or press **STEP** to accept the default value.

Link to stored file
(0-99, 0=SHUT-OFF) 0

Enter the file number of the file you wish to link to in a method, then **ENTER**, or press **STEP** to accept the default link to shut-off.

End of file
RUN-STORE-PRINT

Press **NO** to select one of the options, then **ENTER**. Press **STOP** to return to the first display in the file.

How to Customize a Step-Cycle File, File #4

STEP-CYCLE FILE
STEP to edit file #4

This is the first display in the Step-Cycle file.
 Press **STEP** to edit the first parameter.

Seg 1 Target Temp
 (-5 to 100) 94°C

Seg 1 Segment Time
 (0-999) 1min 0sec

Enter the desired temperatures and times for Segments #1, 2, and 3, which are thermal plateaus:

- Press **STEP** to accept a default value and move to the next display or parameter. Or, to change a value and move to the next display or parameter, use the numeric keys then press **ENTER**.
- Press **HELP** to move to the previous parameter
- Press **STOP** to move to the last display in the file.



Seg 4 Target Temp
 (-5 to 100) 0°C

Seg 4 Segment Time
 (0-999) 0min 0sec

Segment #4 is the terminating segment in the file and must have a time and temperature of zero. The terminating segment marks the end of one cycle.

- Press **STEP** to accept a default value and move to the next display or parameter. Or, to change a value and move to the next display or parameter, use the numeric keys then press **ENTER**.
- Press **HELP** move to the previous parameter

Auto Seg. Extension
NO-YES

Choose **YES** to automatically extend any segment by pressing **NO**, then **ENTER**. Then enter the number of the segment you wish to extend, and the extension time per cycle in seconds. Or, press **ENTER** to choose **NO**.

Cycle Count
 (1-99) 25

Enter the number of cycles you wish to run. Use the numeric keys and **ENTER** to change the value, or press **STEP** to accept the default value.

Link to stored file
 (0-99, 0=SHUT-OFF) 0

Enter the file number of the file you wish to link to in a method, then **ENTER**, or press **STEP** to accept the default link to shut-off.

End of file
RUN-STORE-PRINT

Press **NO** to select one of the options, then **ENTER**. Press **STOP** to return to the first display in the file.

How to Store a Customized User File

Use these instructions to store your customized user files.

End of file RUN-STORE-PRINT	When the last display in your customized file appears, press NO to select the STORE option, then press ENTER .
File Store Enter User #1031	Assign a one- to four-digit User Number (0 to 9999) to your customized file by using the numeric keys and ENTER . <ul style="list-style-type: none">Assign the same user number to a group of files.We suggest using your birthday, in the month-day format, for your user number.
File Store Enter File #7	
Yes to store User #1031 File #7	Assign a one- or two-digit File Number (7 to 99) to your customized file by using the numeric keys and ENTER . <ul style="list-style-type: none">Use this number for direct access to your customized file anytime in the future. In this example, you would press FILE, 7, then ENTER.
	Press YES to store your file. In this example, we assigned User #1031 and File #7 to our customized file.

How to Use the File Directory

Use the File Directory to view all or some of the stored files.

File Directory ALL-BY TYPE-BY USER	Press NO to select the kind of directory you wish to view, then press ENTER .
	This option lets you view all user files with the same user number.
	This option lets you view all files of a specific type. For example, all Soak files, all Time Delay files, etc.
	This option lets you view all stored files in numerical order (by file number)


How to Run a File

Use these instructions to run any file. You can run a file even if you haven't stored it yet.

- | | | |
|--|-----------|---|
| <ol style="list-style-type: none">1. Access the desired file by pressing FILE, the desired file number, then ENTER.2. Press the START key.3. If you have configured your instrument for use with a printer, use the NO key to select a print mode, then press ENTER.4. The run begins. | <p>or</p> | <ol style="list-style-type: none">1. From the last display in a file, press ENTER when the cursor is positioned under the RUN option.2. If you have configured your instrument for use with a printer, use the NO key to select a print mode (OFF, CYCLE, or SEGMENT), then press ENTER.3. The run begins. |
|--|-----------|---|

Run-Time Displays for Cycle Files

While a file is running, a display with information about the run appears. The following example shows the run-time displays for a Thermo-Cycle and a Step-Cycle file run.




How to Create a Method

A method is a series of files linked together in a specific sequence. To create a method:

1. Edit or create the first file in the sequence.
Link this file to the file number of the second file in the sequence, using the "Link to stored file" parameter..
2. Store the file. Remember to assign the correct file number.
3. Edit or create the second file in the sequence.
Link this file to the file number of the third file in the sequence.
4. Edit or create the third file in the sequence.
Link the file to shut-off if it is the last file in the method, unless the last file in the method is a Soak file, in which case you will have to manually abort the run.
5. Store the file. Remember to assign the correct file number.

How to Run a Method

1. Run the first file in the method.
2. When that file has been run, the second file in the method will run automatically.
3. When the second file has been run, the third file in the method will run automatically.
4. When the last file in the method has been run, temperature control will end and the block temperature will drift to ambient.



In the example method on the right, File #60, a pre-PCR cycle incubation Time Delay file, was linked to File #61, a Thermo-Cycle file. This file was linked to File #62, a post-PCR cycle incubation file. Soak files run indefinitely until aborted by the user.

The Configuration File

This file contains the two configuration parameters.

Select Function CONFIG-DIAGNOSTIC	<i>Press ENTER to access the Configuration file.</i>
User Configuration STEP to edit	<i>Press STEP to edit the first configuration parameter.</i>
Pause Time-out limit (9-999) 99sec	<p data-bbox="613 457 1260 516"><i>The pause time-out limit is the time at which a manually paused run will automatically abort.</i></p> <ul data-bbox="613 537 1260 695" style="list-style-type: none"><li data-bbox="613 537 1260 646">• <i>Press STEP to accept a default value and move to the next display. Or, to change a value and move to the next display, use the numeric keys then press ENTER.</i><li data-bbox="613 657 1260 695">• <i>Press HELP to move to the previous display.</i>
Printer Mode OFF-ON	<p data-bbox="613 726 1227 785"><i>Press NO to select the ON option then press ENTER to configure your instrument for printer output.</i></p> <p data-bbox="613 785 1227 842"><i>Or, press STEP to select the OFF option if you are not using a printer.</i></p>
End of config. file STEP to continue	<p data-bbox="613 877 1227 909"><i>Press STEP to re-edit the first configuration parameter.</i></p> <p data-bbox="613 909 1227 936"><i>Press STOP to exit the file.</i></p>

The Diagnostic File

This file contains six diagnostic tests.

Select Function
CONFIG-DIAGNOSTIC

Press **NO** to select *DIAGNOSTIC*, then press **ENTER** to access the tests in the Diagnostic file.

Diagnostic Tests
Enter test # (1-6)

Use the numeric keys and **ENTER** to select the test you wish to run:

Test #1 - Display/Keypad Test

Test #2 - Delete User Files (This test deletes ALL user files.)

Test #3 - Heater Test

Test #4 - Chiller Test

Test #5 - Overshoot Test

Test #6 - Undershoot Test

Using the Keypad

<i>To Do This</i>	<i>Press This Key</i>
Start a run from the first display in the file	Start
View the next display or parameter in a file Accept a default value or condition and move to the next parameter in a file	Step
While editing a file, view the previous display parameter While running a file, toggle between the two run-time displays	Help
While editing a file, move to the last display in the file	Stop
Pause a run in progress	Stop
Abort a run in progress	Stop Stop
Access a file (for example, File #4)	File 4 Enter
Get to the File Directory	File Enter
Get to the Configuration and Diagnostic Files	File Yes
Delete an individual file from the first display in that file	File Clear
Enter numerical values	0 to 9 , Enter
Erase numerical values not yet entered	Clear
Move the screen cursor to the next option in a hyphenated list on a display	No
Accept a selection	Enter

Technical Support

To Reach Us on the Web

Applied Biosystems web site address is:

<http://www.appliedbiosystems.com/techsupport>

We strongly encourage you to visit our web site for answers to frequently asked questions, and to learn more about our products. You can also order technical documents and/or an index of available documents and have them faxed or e-mailed to you through our site (see the "Documents on Demand" section below).

Hours for Telephone Technical Support

In the United States and Canada, technical support is available at the following times.

Product	Hours
Chemiluminescence	9:00 a.m. to 5:00 p.m. Eastern Time
LC/MS	9:00 a.m. to 5:00 p.m. Pacific Time
All Other Products	5:30 a.m. to 5:00 p.m. Pacific Time

See the "Regional Offices Sales and Service" section below for how to contact local service representatives outside of the United States and Canada.

To Reach Us by Telephone or Fax in North America

Call Technical Support at 1-800-831-6844, and select the appropriate option (below) for support on the product of your choice at any time during the call. (To open a service call for other support needs, or in case of an emergency, press 1 after dialing 1-800-831-6844.)

For Support On This Product	Dial 1-800-831-6844, and...	
ABI PRISM® 3700 DNA Analyzer	Press	FAX
	8	650-638-5981
ABI PRISM® 3100 Genetic Analyzer	Press	FAX
	26	650-638-5891
DNA Synthesis	Press	FAX
	21	650-638-5981

For Support On This Product	Dial 1-800-831-6844, and...	
Fluorescent DNA Sequencing	Press 22	FAX 650-638-5891
Fluorescent Fragment Analysis (includes GeneScan® applications)	Press 23	FAX 650-638-5891
Integrated Thermal Cyclers	Press 24	FAX 650-638-5891
Bioinformatics (includes BioLIMS™, BioMerge™, and SQL GT™ applications)	Press 25	FAX 505-982-7690
PCR and Sequence Detection	Press 5 , or call 1-800-762-4001, and press 1 for PCR, or 2 for Sequence Detection	FAX 240-453-4613
FMAT	Telephone 1-800-899-5858, and press 1, then press 6	FAX 508-383-7855
Peptide and Organic Synthesis	Press 31	FAX 650-638-5981
Protein Sequencing	Press 32	FAX 650-638-5981
Chemiluminescence	Telephone 1-800-542-2369 (U.S. only), or 1-781-271-0045 (Tropix)	FAX 781-275-8581 (Tropix) 9:00 a.m. to 5:00 p.m. ET

For Support On This Product	Dial 1-800-831-6844, and...	
LC/MS	Telephone 1-800-952-4716	FAX 650-638-6223 9:00 a.m. to 5:00 p.m. PT

Documents on Demand

Free 24-hour access to Applied Biosystems technical documents, including MSDSs, is available by fax or e-mail.

You can access Documents on Demand through the internet or by telephone:

If you want to order...	Then...
through the internet	Use http://www.appliedbiosystems.com/techsupport You can search for documents to order using keywords. Up to five documents can be faxed or e-mailed to you by title.
by phone from the United States or Canada	a. Call 1-800-487-6809 from a touch-tone phone. Have your fax number ready. b. Press 1 to order an index of available documents and have it faxed to you. Each document in the index has an ID number. (Use this as your order number in step "d" below.) c. Call 1-800-487-6809 from a touch-tone phone a second time. d. Press 2 to order up to five documents and have them faxed to you.
by phone from outside the United States or Canada	a. Dial your international access code, then 1-858-712-0317, from a touch-tone phone. Have your complete fax number and country code ready (011 precedes the country code). b. Press 1 to order an index of available documents and have it faxed to you. Each document in the index has an ID number. (Use this as your order number in step "d" below.) c. Call 1-858-712-0317 from a touch-tone phone a second time. d. Press 2 to order up to five documents and have them faxed to you.

To Reach Us by E-Mail Contact technical support by e-mail for help in the following product areas.

For this product area	Use this e-mail address
Chemiluminescence	info@appliedbiosystems.com
Genetic Analysis	galab@appliedbiosystems.com
LC/MS	apisupport@sciex.com
PCR and Sequence Detection	pcrlab@appliedbiosystems.com
Protein Sequencing, Peptide and DNA Synthesis	corelab@appliedbiosystems.com

Regional Offices Sales and Service If you are outside the United States and Canada, you should contact your local Applied Biosystems service representative.

The Americas	
United States Applied Biosystems 850 Lincoln Centre Drive Foster City, California 94404 Tel: (650) 570-6667 (800) 345-5224 Fax: (650) 572-2743	Latin America (Del.A. Obregon, Mexico) Tel:(305) 670-4350 Fax: (305) 670-4349

Europe	
Austria (Wien) Tel: 43 (0)1 867 35 75 0 Fax: 43 (0)1 867 35 75 11	Hungary (Budapest) Tel: 36 (0)1 270 8398 Fax: 36 (0)1 270 8288
Belgium Tel: 32 (0)2 712 5555 Fax: 32 (0)2 712 5516	Italy (Milano) Tel: 39 (0)39 83891 Fax: 39 (0)39 838 9492
Czech Republic and Slovakia (Praha) Tel: 420 2 61 222 164 Fax: 420 2 61 222 168	The Netherlands (Nieuwerkerk a/d IJssel) Tel: 31 (0)180 331400 Fax: 31 (0)180 331409
Denmark (Naerum) Tel: 45 45 58 60 00 Fax: 45 45 58 60 01	Norway (Oslo) Tel: 47 23 12 06 05 Fax: 47 23 12 05 75

Europe <i>(continued)</i>	
Finland (Espoo) Tel: 358 (0)9 251 24 250 Fax: 358 (0)9 251 24 243	Poland, Lithuania, Latvia, and Estonia (Warszawa) Tel: 48 (22) 866 40 10 Fax: 48 (22) 866 40 20
France (Paris) Tel: 33 (0)1 69 59 85 85 Fax: 33 (0)1 69 59 85 00	Portugal (Lisboa) Tel: 351 (0)22 605 33 14 Fax: 351 (0)22 605 33 15
Germany (Weiterstadt) Tel: 49 (0) 6150 101 0 Fax: 49 (0) 6150 101 101	Russia (Moskva) Tel: 7 095 935 8888 Fax: 7 095 564 8787
Spain (Tres Cantos) Tel: 34 (0)91 806 1210 Fax: 34 (0)91 806 1206	South Africa (Johannesburg) Tel: 27 11 478 0411 Fax: 27 11 478 0349
Sweden (Stockholm) Tel: 46 (0)8 619 4400 Fax: 46 (0)8 619 4401	United Kingdom (Warrington, Cheshire) Tel: 44 (0)1925 825650 Fax: 44 (0)1925 282502
Switzerland (Rotkreuz) Tel: 41 (0)41 799 7777 Fax: 41 (0)41 790 0676	South East Europe (Zagreb, Croatia) Tel: 385 1 34 91 927 Fax: 385 1 34 91 840
Middle Eastern Countries and North Africa (Monza, Italia) Tel: 39 (0)39 8389 481 Fax: 39 (0)39 8389 493	Africa (English Speaking) and West Asia (Fairlands, South Africa) Tel: 27 11 478 0411 Fax: 27 11 478 0349
All Other Countries Not Listed (Warrington, UK) Tel: 44 (0)1925 282481 Fax: 44 (0)1925 282509	

Japan
Japan (Hatchobori, Chuo-Ku, Tokyo) Tel: 81 3 5566 6100 Fax: 81 3 5566 6501

Eastern Asia, China, Oceania	
Australia (Scoresby, Victoria) Tel: 61 3 9730 8600 Fax: 61 3 9730 8799	Malaysia (Petaling Jaya) Tel: 60 3 758 8268 Fax: 60 3 754 9043
China (Beijing) Tel: 86 10 6238 1156 Fax: 86 10 6238 1162	Singapore Tel: 65 896 2168 Fax: 65 896 2147
Hong Kong Tel: 852 2756 6928 Fax: 852 2756 6968	Taiwan (Taipei Hsien) Tel: 886 2 2698 3505 Fax: 886 2 2698 3405
Korea (Seoul) Tel: 82 2 593 6470/6471 Fax: 82 2 593 6472	Thailand (Bangkok) Tel: 66 2 719 6405 Fax: 66 2 319 9788

Worldwide Sales Offices

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches into 150 countries on six continents. For international office locations, please call our local office or refer to our web site at www.appliedbiosystems.com.

Headquarters

850 Lincoln Centre Drive
Foster City, CA 94404 USA
Phone: +1 650.638.5800
Toll Free: +1 800.345.5224
Fax: +1 650.638.5884

Technical Support

For technical support:
Toll Free: +1 800.831.6844 ext 23
Fax: +1 650.638.5891

www.appliedbiosystems.com



PE Corporation is committed to providing the world's leading technology and information for life scientists. PE Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Printed in the USA, 10/2000
Part Number 0993-8984B