



Belimed
Infection Control

This manual is valid for:

Product	Steam sterilizer MST 5-5-9
Product version	TOP 5000
Local regulations	American Society of Mechanical Engineers (ASME)

Editor	Belimed Sauter AG, CH-8583 Sulgen
Date of issue	July 2013



Table of contents

USER MANUAL.....	11
1 INTRODUCTION.....	13
1.1 Become familiar with the system.....	13
1.2 Intended use.....	13
1.3 Summary of safety precautions.....	14
2 INSTALLATION.....	17
2.1 Sterilizer sizes.....	17
2.2 Required building system utilities.....	18
2.2.1 Utility connections.....	18
2.2.2 Utility connection positioning and dimensions.....	19
2.2.3 Utility connections.....	20
2.3 Installation instructions.....	20
2.3.1 Introduction.....	20
2.3.2 Space considerations.....	21
2.3.3 Mounting details.....	21
2.3.4 Utility service requirements.....	21
2.3.5 Preparation.....	22
3 LABELING.....	24
3.1 Type and Model Designation.....	24
3.2 Marking.....	24
3.3 Warning labels.....	25
4 GENERAL DESCRIPTION OF THE COMPONENTS AND FUNCTIONS.....	27
4.1 General view.....	27
4.2 Indicators and gauges.....	28
4.3 Operating End Control Panel.....	29
4.4 Non-Operating End Control Panel (option).....	30

4.5	Printer.....	31
4.6	Printouts (example).....	31
4.7	Door operation.....	32
4.7.1	Manual doors	32
4.7.2	Motorized doors	32
4.8	General precautions for the sterilization procedure	33
5	OPERATING THE STERILIZER.....	34
5.1	Turn on the sterilizer	34
5.2	Log on	35
5.3	Loading the sterilizer	37
5.3.1	Loading sterilizer	37
5.4	Select a cycle.....	39
5.4.1	Factory-set cycles and cycle values	39
5.4.2	Select a cycle.....	39
5.5	Start a cycle	39
5.6	In-cycle.....	40
5.7	Cycle failed	40
5.8	Unloading the sterilizer	41
5.8.1	Unloading sterilizer.....	41
5.9	Cycle documentation	41
5.10	Log off	42
5.11	Turn off the sterilizer	42
6	STERILIZER FACTORY CYCLES SETTINGS	43
6.1	270°F Pre-vacuum Cycles	44
6.2	250°F Liquid Cycle.....	45
6.3	270°F Express Cycle.....	47
6.4	270°F Flash Cycles	48

6.5	DART / Bowie-Dick Test cycle	49
6.6	Warm-up & Leak Test cycle	50
7	TECHNIQUES OF STERILIZATION	51
7.1	Guidelines for preparing and sterilizing wrapped packs.....	51
7.1.1	Preparing load:	51
7.1.2	Wrapped fabric packs	51
7.1.3	Wrapped instrument sets	52
7.1.4	Wrapped utensils	53
7.1.5	Loading the shelves:.....	54
7.1.6	Steam sterilization cycles.....	55
7.1.7	Material qualification	55
7.2	Guidelines for preparing and sterilizing liquids	56
7.2.1	Preparing load	56
7.2.2	Liquid cycle	57
7.3	Appropriated use of Flash cycle.....	57
7.4	Appropriated use of Express cycle	57
7.5	Guideline for determining wet packs	58
7.5.1	Introduction	58
7.5.2	Evaluation of wet packs	58
7.5.3	Moisture retention acceptance criteria	58
7.5.4	Solving wet pack problems	58
7.6	Sterility maintenance	59
7.6.1	Cooling	59
7.6.2	Inspection	59
7.6.3	Sterile storage.....	59
7.6.4	Expiration dating	59
8	OPERATING DEVICE TOP5000.....	60
8.1	User administration	60
8.1.1	User's rights	60
8.1.2	Direct login.....	62
8.1.3	Administrate user	63
8.2	Info	64

8.3	Readings	65
8.4	Set-up	66
8.4.1	Early start / night-time shutdown	67
8.4.2	Set clock	67
8.4.3	Language and units.....	68
8.4.4	Steam generator	69
8.4.5	Drain steam generator	69
8.4.6	Configure cycle selection	69
8.4.7	Display contrast	69
8.4.8	Online docu.....	70
8.5	Print last cycle	70
8.6	Alarm messages	70
8.6.1	Critical alarms	71
8.6.2	Alarm list	71
8.7	Cycle intervention	74
8.7.1	Cycle advance	74
8.7.2	Cycle abort.....	74
8.7.3	Cycle continue	75
8.8	Batch data	75
8.8.1	Manual data entry	75
8.8.2	Batch data entry with Documentation system 8535-BC	75
8.8.3	Batch data input via scanner with ICS 8665 Documentation System	76
8.9	Maintenance messages	77
8.10	Maintenance functions	77
8.10.1	Replace door seal	78
8.10.2	Calibrate touch screen	79
8.10.3	Sensors adjustment	79
8.10.4	Machine set-up	79
8.10.5	Maintenance interval.....	83
8.10.6	Printer self-test.....	84
8.10.7	Modify Cycle Parameters and copy program	85
8.10.8	PLC functions.....	86
8.10.9	I/O- Test.....	87
9	ROUTINE MONITORING	89
9.1	Biological monitoring	89

9.2	Testing for pre-vacuum efficiency.....	89
9.2.1	Leak Test.....	89
9.2.2	Bowie-Dick Test.....	89
10	ROUTINE MAINTENANCE.....	90
10.1	Replacement of ribbon or paper roll.....	90
10.2	Cleaning of Chamber and Front panels.....	92
	TECHNICAL MANUAL.....	93
11	GENERAL SAFETY PRECAUTIONS	95
11.1	Summary of warnings and cautions	95
12	MAINTENANCE / PERIODIC TESTS.....	97
12.1	Safety information	97
12.2	Maintenance intervals	98
12.3	Maintenance personnel.....	98
12.4	Maintenance plan	98
12.4.1	Daily Maintenance (before start sterilizing)	98
12.4.2	Weekly Maintenance	99
12.4.3	Maintenance level 1.....	100
12.4.4	Maintenance level 2.....	101
12.4.5	Maintenance level 3.....	102
12.4.6	Yearly safety check.....	102
12.4.7	Periodic pressure vessel check	102
13	CLEANING	103
14	CALIBRATION VERIFICATION.....	104
14.1	Temperature calibration verification.....	104
14.2	Chamber pressure transducer calibration verification	105
14.3	Jacket and steam pressure transducer calibration verification	105
15	SENSORS ADJUSTMENT	106

15.1	Temperature sensors adjustment	106
15.1.1	Offset calibration:	106
15.1.2	Gain calibration:	106
15.2	Pressure sensors adjustment.....	107
15.2.1	Offset calibration:	107
15.2.2	Gain calibration:	107
16	SERVICING INSTRUCTIONS.....	108
16.1	Opening control panel.....	108
16.2	Cleaning and replacement of door seal.....	108
16.2.1	Door seal on operating end.....	108
16.2.2	Door seal on non operating end.....	109
16.3	Adjust door sliding clutch (option motorized doors)	110
16.3.1	Check the closing force.....	110
16.3.2	Door sliding clutch adjustment	110
16.4	Replacing the batteries (accumulators).....	111
16.5	Arrangements to be taken before and after an operation interruption	111
17	TESTS AND INSPECTIONS	112
17.1	Emergency-Stop (double door sterilizers only)	112
17.2	Safety valves.....	112
17.2.1	General:	112
17.2.2	Inspection in fitted condition.....	112
17.2.3	Inspection in removed condition.....	113
17.3	Safety equipment, door lock.....	113
17.3.1	Chamber atmospheric pressure, safety switch.....	113
17.3.2	Door seal pressure.....	114
17.3.3	Door limit switches and steam supply to chamber monitoring.....	114
17.3.4	Door limit switch adjustment	115
18	SETTINGS	116
18.1	Pressure switches, pressure regulators and safety valves	116
19	TESTS AFTER ELECTRICAL REPAIR	117

20	LOGGING THE MAINTENANCE WORK CARRIED OUT	118
21	LIST OF FIGURES	119
22	LIST OF TABLES	121
23	ELECTRICAL WIRING DIAGRAM STERILIZER	123
24	BILL OF MATERIAL TO ELECTRICAL WIRING DIAGRAM STERILIZER.....	161
25	ELECTRICAL WIRING DIAGRAM STEAM GENERATOR V30 (OPTIONAL).....	169
26	ELECTRICAL WIRING DIAGRAM STEAM GENERATOR V30 (OPTIONAL).....	185
27	P & I DIAGRAM STERILIZER WITHOUT COOLING CIRCUIT	187
28	BILL OF MATERIAL TO P & I DIAGRAM STERILIZER WITHOUT COOLING CIRCUIT	189
29	P & I DIAGRAM STERILIZER WITH COOLING CIRCUIT	195
30	BILL OF MATERIAL TO P & I DIAGRAM STERILIZER WITH COOLING CIRCUIT.....	197
31	P & I DIAGRAM ELECTRIC STEAM GENERATOR (OPTIONAL)	203
32	BILL OF MATERIAL TO P & I DIAGRAM ELECTRIC STEAM GENERATOR.....	205

Belimed Steam Sterilizer TOP 5000 USER MANUAL



Model 5-5-9
535x535x965 mm

Limitation of liability and indemnity

In no event, whether as a result of breach of contract, warranty or tort (including negligence and strict liability) shall Belimed Sauter AG or its suppliers be liable for any consequential or incidental damage including, but not limited to loss of profits or revenues, loss of use of the Products or any associated equipment, loss of the Buyer's Products, damage to associated equipment, cost of capital, cost of substitute products, facilities, service or replacement power, downtime costs, caused by such Products or claims of the users for such damages. Buyer and ultimate user hereby agree to indemnify Belimed Sauter AG and to hold Belimed Sauter AG harmless from any and all liability for such consequential and incidental damages. The responsibility of Belimed Sauter AG for damages due to injuries or death or death of employees of the Buyer or ultimate user of the Product, caused by the product, shall be limited to that portion of such damages as might be attributable to the negligence or strict liability of the Belimed Sauter AG Company. The Buyer and ultimate user agree to indemnify Belimed Sauter AG and hold Belimed Sauter AG harmless for any further damages, indemnity or contribution. If Buyer transfers title to or leases the Products sold to any third party, buyer shall obtain from such third party a provision affording Belimed Sauter AG and its suppliers the protection of this article relating to Limitations of Liability and Indemnities.

The Belimed Sauter AG's liability for any claim of any kind (including negligence and strict liability) for any loss or damage arising out of, or resulting from this agreement, or from the performance of breach thereof, or from the Products or Services furnished hereunder, shall in no case exceed the price of the specified Product, system, component or service which gives rise to the claim. Except as to title, any such liability shall terminate one year from the date of installation of any Product or upon the expiration of the warranty period applicable to each type of product covered hereby, whichever time period expires first.

IMPORTANT: A SUMMARY OF THE SAFETY PRECAUTIONS TO BE OBSERVED WHEN OPERATING THIS EQUIPMENT CAN BE FOUND ON PAGE 14 OF THIS MANUAL AND BEHIND THE SERVICE DOOR OF THE UNIT. DO NOT OPERATE THE STERILIZER UNTIL YOU HAVE BECOME FAMILIAR WITH THIS INFORMATION.

Manufacturer

Belimed Sauter AG

Zelgstrasse 8

CH- 8583 Sulgen Switzerland

Phone 011 41 71 644 85 00

Fax 011 41 71 644 86 01

1 Introduction

1.1 Become familiar with the system

This manual contains important information on proper use and maintenance of this sterilizer. All operators and departments heads are urged to carefully review and become familiar with the warnings, cautions and instructions contained herein. This sterilizer is specifically designed to process goods using only the cycles as specified in this manual. If there is any doubt about a specific material or product, contact the manufacturer of the product for the recommended sterilization technique.

1.2 Intended use

The Belimed Steam Sterilizer TOP 5000 is designed to be used for the terminal sterilization of porous and non porous, heat and moisture stable materials in the healthcare facilities.

The sterilizer is factory equipped with standard cycles: their applicability is insured only under the specified conditions. Depending of the chosen cycle, materials as different as textiles, glassware, unwrapped or wrapped instrument trays with single or multiple instruments can be sterilized.

The Belimed Steam Sterilizer TOP 5000 is factory equipped with cycles which has been tested in accordance with AAMI/ANSI ST-8 as well as EN285 standards under defined load conditions. The cycles parameters can be easily modified to perform the sterilization under optimized conditions. A modification of the parameters of the standard cycles is under the responsibility of the user. BELIMED can offer assistance or quote complete validation procedures to the relevant standards.

WARNING

This sterilizer is not designed to process flammable liquids nor liquids in containers that are not designed for sterilization.

CAUTION

Any alteration of the sterilizer which affects its operation will void the warranty and could violate state and local regulations and jeopardize insurance coverage.

CAUTION

Use only original BELIMED replacement parts for maintenance purposes: in all other cases the manufacturer can no more guaranty the specified performances of the equipment nor apply the warranty conditions.

CAUTION

In order to ensure the proper function of the sterilizer, we recommend to perform the leak-rate-test at least once per week, at normal working temperature.

CAUTION

The Bowie-Dick Test Cycle has to be performed daily.

1.3 Summary of safety precautions

The following is a summary of safety precautions to be observed when operating or servicing this unit.

WARNINGS indicate the potential for danger to personnel and CAUTIONS indicate the potential for damage to equipment. The precautions may be repeated where applicable throughout the manual. The following is a list of all safety precautions to be taken. Carefully read them before proceeding to use or service the unit.

WARNINGS

- **RESTRICTION:** It is inappropriate for a health care facility to sterilize liquids for direct patient contact.
- Before entering the chamber (e.g. for cleaning), turn the key switch off (position '0') and keep the key in your pocket as long as you work in the chamber.
- **Danger of Squeeze Hand:** Hands off from chamber door area while closing the door.
- When closing the door, no person is to be inside the chamber. The operator is responsible for the observation of this safety precaution.
- **BURN HAZARD:** The chamber door surface can be hot ($>70^{\circ}\text{C}/158^{\circ}\text{F}$), if the ambient temperature is above $30^{\circ}\text{C}/86^{\circ}\text{F}$.
- **BURN AND SHOCK HAZARD:** Repairs and adjustments should be attempted only by authorized persons fully acquainted with this equipment. Use of inexperienced, unqualified persons to work on the equipment or the installation of unauthorized parts could cause personal injury or result in costly damage!
- Observe the label "Dangerous voltage": turn off main switch before opening
- **BURN HAZARD:** Allow sterilizer, generator (if applicable) and accessories to cool to room temperature before performing any cleaning or maintenance procedures.
- **BURN HAZARD:** Sterilizer and rack/shelves will be HOT after cycle is run. Always wear protective gloves and apron when removing a processed load. Protective gloves and apron should also be worn when reloading sterilizer following previous operation.
- **FALL HAZARD:** To prevent falls keep floors dry by immediately wiping up any liquids or condensation in sterilizer loading or unloading area.
- **EXPLOSION HAZARD:** This sterilizer is not designed to process flammable liquids.
- **BURN HAZARD:** When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedure:

- Use LIQUIDS cycle only. No other cycle is safe for processing liquids.
 - Use only vented closures- do not use screw caps or rubber stoppers with crimped seal.
 - Use only Type I borosilicate glass bottles- do not use ordinary bottles or any container not designed for sterilization.
 - Avoid sudden opening of door at end of cycle. Wait at least 10 minutes before door opening and unloading the sterilizer.
 - Do not allow hot bottles to be jolted. This can cause hot-bottle explosions! Do not move bottles if any boiling or bubbling is present.
 - Allow bottles to cool to touch before attempting to move them from sterilizer shelf to storage area.
-
- **BURN HAZARD:** A steam supply malfunction may cause the sterilizer chamber to fill with scalding water. Do not open chamber door if the unit fails to complete an automatic cycle or if water leaks past the door gasket upon unlocking the door.

 - **STERILITY ASSURANCE HAZARD:** Load sterility may be compromised if the chemical or biological indicator or the Bowie-Dick Test (or DART) indicate a potential problem. If these indicators show a potential problem, refer the situation to a qualified maintenance technician before using the sterilizer further.

 - **STERILITY ASSURANCE HAZARD:** According to AAMI ST-8 a measured leak rate greater 1,3 mbar/minute indicates a problem with the sterilizer. Refer the situation to a qualified maintenance technician before using the sterilizer further.

CAUTIONS

CAUTION: Never use sharp tools to push gasket into groove.

CAUTION: Never use wire brush or steel wool on door and chamber assembly.

CAUTION: Observe the Electrostatic Precautions outlined in the maintenance manual. Always wear a grounding wrist strap when removing or replacing PC boards or ICs.

CAUTION: Solenoid valves are equipped with a special material which can be attacked by oils and grease. When replacing entire valve, wipe threads clean of cuttings oils and use Teflon tape to seal pipe joints.

CAUTION: Handle siphon and bellows assembly gently to avoid damage.

CAUTION: Keep the movement area of the door unobstructed during door opening and closing.

CAUTION: Any alteration of the sterilizer which affects its operation will void the warranty and could violate state and local regulations and jeopardize insurance coverage.

CAUTION: Use only original BELIMED replacement parts for maintenance purposes: in all other cases the manufacturer can no more guaranty the specified performances of the equipment nor apply the warranty conditions.

CAUTION: In order to ensure the proper function of the sterilizer, we recommend to perform the leak-rate-test at least once per week, at normal working temperature.

CAUTION: The Bowie-Dick Test Cycle has to be performed daily.

CAUTION: Sterilize only goods which are declared to be steam-sterilizable by their manufacturer and compatible with the chosen cycle parameters (maximum temperature)

CAUTION: Respect the maximum number of sterilization's given by the manufacturer of the goods to be sterilized. Check their proper function according to the manufacturer's information of the goods (e.g. instruments)

2 Installation

2.1 Sterilizer sizes

The Belimed steam sterilizers are available in the following sizes (Table 2-1):

Model	Configuration	Chamber volume (liters)	Chamber size (H x W x D) (mm)	Overall Dimensions (H x W x D) (mm)	Weights (kg) ¹ Transport / Operation / Inspection
5-5-9 VS1	1- door	275	535x535x965	1900x855x1174	800 650 900
5-5-9 VS2	2- door	275	535x535x915	1900x855x1174	820 670 920

Table 2-1: Sterilizer sizes

Tolerances: Overall Dimensions ± 3 mm

¹ Weight without internal steam generator

2.2 Required building system utilities

Unless otherwise specified in the contract the utilities (refer to Table 2-2 and Table 2-3) must be provided by the customer to the connection points (refer to Figure 2-1 to Figure 2-2)

2.2.1 Utility connections

Symbol	Medium	Pressure in bar g	Customer's supply	Peak	Consumption per batch with normal load
Standard-Supply lines					
SD	Sterilizing Steam, safeguarding 60 psi g	2.5- 4	¾" NPT-f	40 kg/h	9 kg
KW	Tap water for vacuum pump, 15 °C (60°F)	2 - 5	½" NPT-f	1 m³/h	150 l 10 l
VE	Deionized water 15-80°C for internal steam generator only	2 - 5	¼" NPT-f	0.1 m³/h	9.5 l
DL	Compressed air, oil-free	5 - 7	¼" NPT-f	5 Nm³/h	0.25 Nm³
EL1	Electric supply: external steam 3~ 208V 60 Hz (branch circuit breaker)		6 A (10A)	1.6 kW	0.5 kWh
EL2	Electric supply: Option internal steam generator 3~ 480V 60 Hz (branch circuit breaker)		37 A (50A)	30 kW	8 kWh
ZL	Air inflow to service room coming through customer's ventilation duct or from operating room.				
Standard-Waste line					
ALP	Free air exhaust of vacuum pump				
A	Floor drain (waste water open)		1 ¼" NPT-f	max. 15 l/min	
GU	Floor drain				
AL1	Air outflow from service room: heat flow of 1 kW to be dissipated, temperature in service room ~30°C			1.5 kW	
Option waste steam line upwards					
AD1	Waste steam, chamber	≤ 2.75	1" NPT-m	100 kg/h	
AD2	Waste steam, jacket	≤ 2.75	1" NPT-m	100 kg/h	
Option Cooling water circuit					
KWV	Cooling water forward: T1 10 °C, ΔT 15K Cooling energy required	2 - 5	¾ " NPT-f	2 m³/h 8kW	400l 26000kJ
KWR	Cooling water backward		¾ " NPT-f		

Table 2-2: Utilities for chamber size 5-5-9

2.2.2 Utility connection positioning and dimensions




Figure 2-1: Utility connection positions and dimensions

2.2.3 Utility connections

The customer shall provide all supply connection lines with a built in manual stop-cock. All connection assemblies between the sterilizer and supply / exhaust lines (customer's supply) are designed with flexible, pressure-resistant tubes made of 316L material: in that way, mechanical tolerances of ± 25 mm in any direction can be easily compensated which results in short assembly times. Only exception is the steam supply <SD> which is realized in a rigid manner.

When the exhaust line will be neither connected nor led outside the tubing ends in a bend oriented upwards.




Figure 2-2: Utility connections

2.3 Installation instructions

2.3.1 Introduction

The installation of the sterilizer will generally be performed by Belimed technicians. They received a detailed training on all aspects of the installation, running and maintenance of the Belimed Sterilizers. Those details of the installation are part of their specific know how and will not be written down here. If the customer requires to perform the installation himself, Belimed will provide the necessary information and at least insure the supervision of the installation. This is required in order to be able to insure the warranty of the system for the accepted duration.

Each installation drawing pertains to the sterilizing equipment as specified or purchased by the customer. The general notes and recommendations are intended to complete the installation drawings.

2.3.2 Space considerations

Clearance in front of sterilizer, for comfortable loading and unloading operations, should be about twice the length of the transfer carriage unless otherwise specified on installation drawing.

2.3.3 Mounting details

2.3.3.1 Service power

One convenience outlet (120V) is required for power tools next to the sterilizer.

2.3.4 Utility service requirements

2.3.4.1 Terminal Fittings

Unless otherwise specified in the contract piping, shutoff valves and other appurtenances between terminal fittings on the equipment and wall or floor outlets are not furnished by Belimed.

Note: "Supply lines are to be safeguarded and separable in accordance with connection ratings"

2.3.4.2 Pipe sizes

Pipe sizes listed in paragraph 3.2 indicate the equipment termination sizes only. Size of supply piping is dependent on length of pipe run from pressure regulation station for steam line and main water headers to ensure adequate supply service pressure and demand flow at equipment terminals. Effect of coincident draw of multiple unit installations must also be considered.

2.3.4.3 Backflow preventer

If local codes require a reduced pressure principle device on water supply line it shall be provided by the customer.

2.3.4.4 Pressure relief valves

Belimed recommends piping all chamber relief valves to a vented manifold outside the equipment service area .

Caution:

Do not reduce the discharge capacity of the safety relief valve.

Recommended piping practices for relief valve piping can be found in Pressure equipment directive 97/23/EG.

2.3.4.5 Blow down piping

Blow down building steam, water and compressed air supply lines before final connection to equipment.

2.3.4.6 Steam and water pressure

Sterilizer is equipped to operate on pressures listed in Table 2-2 and Table 2-3. If supply pressure exceeds those shown, provide reducing valves. Unless otherwise specified in the contract reducing valves are not furnished by Belimed.

2.3.4.7 Steam quality

Steam should be condensate free and 100% saturated vapor (dry, not superheated) to ensure proper goods drying. It has to fulfill the requirements of AAMI ST 46.

2.3.4.8 Water quality

Water is used for vacuum pump and heat exchangers. Refer to Table 2-3 for recommended water quality.

Condition	Nominal conditions	Maximum conditions
Temperature	4-16°C (40-60°F)	21°C (70°F)
PH	6.5-8	6-8.5
Total hardness as Ca ²⁺	0.2-1.8 mmol/l	2 mmol/l
Total dissolved solids	50-200 mg/l	400 mg/l
Chloride	10-70 mg/l	80 mg/l


Table 2-3: Recommended feed water quality for sterilizers

2.3.4.9 Wiring Terminals

It is the responsibility of the customer to provide proper electrical disconnect device and branch circuit protection in accordance with applicable codes and regulations.

Wiring to the disconnect device is not furnished by Belimed.

Connect the supply line on the main disconnect switch terminals (refer to Figure 2-3)



Wire (cable) connector no
Phase L1 to 1/L1
Phase L2 to 3/L2
Phase L3 to 5/L3
Protective Conductor to PE

Figure 2-3: Disconnect switch

The electrical installation must be carried out by an installer in the fields.

2.3.5 Preparation

2.3.5.1 Standard sterilizers

- Check the sterilizer placement according to installation drawing.
- Measure overall height on the wall opening.
- Mark the center lines of the sterilizer chamber on the floor with a color marker
- For correct positioning refer to installation drawing
- Check the position and carrying out of the floor drains (A and GU)
- Put the stainless steel tub on place (if required).
- Remove transportation screws from all door-counter weights.
- Place the chamber on the stainless steel tub.
- If the sterilizer is installed on a flammable floor, a steel bottom plate, 800mm wide and 2mm thick, must be placed on the floor below the service area, covering the whole service area.

- Adjust the chamber high using the adjusting screws. Refer to the installation drawing.
- Open chamber door
- Remove door locking bolts (below chamber) and insert the bolt 90° turned again.
- Adjust the chamber high with a spirit level along x and y axis by using the 4 adjustment screws on the feeds.
The chamber floor must be leveled to 837mm above ground.
- Close chamber door.
- Tighten the gap between steel tub and sterilizer ground frame with silicone compound.
- Install the 4 pressure sensors.
- Fill 5 liters of water into the vacuum pump circulation tank.

3 Labeling

3.1 Type and Model Designation

The Belimed Steam Sterilizer TOP 5000 model 5-5-9 is available in the standard configuration.

Table 3-1 hereafter resumes size and designation.

Model	Configuration	Chamber volume (liters)	Chamber size (H x W x D) (mm)	Overall Dimensions (H x W x D) (mm)	Weights (kg) ² Transport / Operation / Inspection
5-5-9 VS1	1- door	275	535x535x965	1900x855x1174	650
5-5-9 VS2	2- door	275	535x535x915	1900x855x1174	650

Table 3-1: Size and Designation

3.2 Marking

The sterilizer can be identified through following two metal marking plates

- Chamber/jacket ratings (fig. 4-1) are located just above the chamber door on the operating end.




Figure 3-1: Chamber labeling (according to ASME pressure and fire vessels, example)

² Weight without internal steam generator

- Electrical ratings (Figure 3-2) are given on the label located in the upper part inside the servicing door on the operating end.




Figure 3-2: Labeling Electrical Ratings

Remark to Figure 3-2:

Working pressure: 1 bar (gauge) = 14.5psi / 2.75 bar g = 40 psi g

Voltage: 3~ = 3 phase

3.3 Warning labels

Those labels are placed close to the respective parts which represent a potential hazard:

- DANGER OF SCALDING !




Figure 3-3: Label "DANGER OF SCALDING !"

- DANGER: OVERPRESSURE !




Figure 3-4: Label "Danger overpressure !"

- DANGEROUS VOLTAGE !




Figure 3-5: Label „DANGEROUS VOLTAGE !“

- DANGER OF GET ONE'S HAND CRUSEHD



Figure 3-6: Label „Danger of get one's hand crushed" while closing chamber door

- BURN HAZARD



Figure 3-7: Label „Burn Hazard"

- CAUTION refer to accompanying documents



Figure 3-8: Label „CAUTION - REFER TO MANUAL"

4 General description of the components and functions

4.1 General view

Figure 4-1 shows the main components or control elements which are accessible to the user from the operating end. The number of control elements of the sterilizer is reduced to the absolute minimum in order to maximize the simplicity and security of use.




Figure 4-1: Control elements operating end

Figure 4-2 shows the main components or control elements which are accessible to the user from the non operating end.




Figure 4-2: Control elements non operating end

4.2 Indicators and gauges

The system offers two types of indicators and gauges:

- The process indicators are all located on the display of the operating end control panel and of the non operating end control panel (see section 5.3 and 5.4)
- The service indicators and gauges are all located within the service area (see service and maintenance manual). The pressure gauges in the service area scaled in bar gauge:

The gauges shall show :

Pressure indication	Nominal conditions bar	Maximum conditions bar
Steam supply	2.5-3	2.5-3.5
Chamber	-0.96-2.3	-1..+2.7
Jacket	-0.96-2.3	-1..+2.7
Door seal (door locked)	2.9-3.2	2.8-3.5
Compressed air	5-7	5-7
Tap water	2-5	2-5
Deionized water	2-5	1.5-5

Table 4-1: Pressure indicators in service area




Figure 4-3: Instrumentation

4.3 Operating End Control Panel

Figure 4-4 shows the arrangement of the control elements on the operating end.

The user's control panel with the interactive touch-screen display is described with more details in Figure 4-5 and Figure 4-6 and hereafter.

The chamber temperature indicator can be configured in °F or °C.

The chamber pressure indicator can be configured in mbara, barg, kPaa, psia, psig/in.Hg



Figure 4-4: Operating end control panel

P1 chamber	T1 chamber	status	
14.78 psia	258.1 °F	standby	
more ..			
1 PreVac 270 long Dry		2 PreVac 270 short Dry	
		4 Flash 270 single Ins	
5 Flash 270 multi Ins			cycle start
PreVac 270 long Dry		open door	menu
cycle counter	0		

Figure 4-5: Out-of-cycle display

P1 chamber	T1 chamber	status	door locked
860 mbar	101.8 °C	2. vacuum pulse	
Instruments 134°C		cycle	menu
remain time: 44.4 min		intervention	

Figure 4-6: In-cycle display

4.4 Non-Operating End Control Panel (option)

Sterilizers with 2 doors are equipped with the non-operating end control panel.

Figure 4-7 shows the panel on the non operating end.



Figure 4-7: Non operating end control panel

No.	Item	Name	Function
1	Push button	Open door	Open chamber door at end of cycle
2	Push button	Close door	Close chamber
3	Push button		n/a
4	LED	Open door	Lights up, when the door on the non operating end can be opened
5	LED	In cycle	Lights up while a cycle is running
6	LED	Door locked	Lights up, when doors are locked
7	LED	Alarm	Lights up while an alarm is active
8	Bar Graph	Cycle progress	Shows the progress of a cycle (0...100%)
9	Display	Temperature	Chamber temperature Indicator or remain process time
10	Display	Pressure	Chamber pressure Indicator
11	Push button	Emergency stop	Emergency stop button

4.5 Printer

The cycle documentation printer is a panel mount printer (refer to Figure 4-8).



General Specifications:

Technology: impact dot matrix

Columns: 42 columns

Paper with: 57.5mm

Paper roll diameter 38mm

For replacing paper roll or ribbon refer to chapter 11

Figure 4-8: Panel mount printer

4.6 Printouts (example)

```
BELIMED CYCLE DOCUMENTATION
Hospital      : Charity
Department   : CSSD
Machine-Type  : 12-6-12HS2   No: 045556
Run number    : 345679
Operator      : Brown
Cycle         : 4: Leak Test
Version       : PR: 13.03.2001 SW: 2.2
Cycle start   : 24.03.2001 / 13:35
```

```
Set-Values   : Leak rate <1.3 mbar/min
               Test time 15.0 min
```

Time	Phase	Press.	Temp.
m:s		mbara	°C
0:00	Vacuum	995	28.5
2:20	Stabilize	70	25.5
7:20	Air Leak Test	72	27.9
22:20	Air Break	79	28.3
24:03	Complete	980	29.1

```
Leak rate 0.5 mbar/min
```

CYCLE PASSED

```
Cycle approved: Y / N   Date:
```

Signature:

Figure 4-9: Printout of a Leak Test cycle

4.7 Door operation

4.7.1 Manual doors

The door can be opened by pressing the control button “open door”. The door will be unlocked by sucking back the door seal. Five seconds after pressing the “open door” button the door can be opened.

Close the door manually. When the door is closed the door will be locked automatically.

4.7.2 Motorized doors

If the doors of the Belimed Steam Sterilizer TOP 5000 are motorized: they can be activated by pressing the control buttons “open door” and “close door” on the control panel of the operating or non-operating end: the corresponding door will be activated when the sterilizer status allows the opening or closing procedure. The door motor is equipped with a safety clutch, which stops the door if the closing force is >150N.

F1 chamber	T1 chamber	status	
14.78 psia	258.1 °F	standby	
more ..			
1 PreVac 270 long Dry		2 PreVac 270 short Dry	
		4 Flash 270 single Ins	
5 Flash 270 multi Ins			cycle start
PreVac 270 long Dry		open door	menu
cycle counter	0		

Figure 4-10: Door operation buttons operating end (manual door)



Figure 4-11: Door operation buttons non operation end

4.8 General precautions for the sterilization procedure

The sterilization performances in terms of cycle duration, level of dryness and water/steam consumption can only be achieved when the loading prescriptions are respected. Following aspects are to be followed:

- The goods to be sterilized must be placed on the loading racks especially designed therefore: avoid placing the goods directly onto the chamber floor.
- Each shelf of the loading rack must be loaded individually: do not stack trays.
- When wire baskets are stacked, pay attention that the goods of the underlying basket shall not be compressed. The package must be able to expand during the sterilization process.
- Heavy sterilization items shall always be placed in the lower part of the unit.
- Overfilled baskets and trays can cause problems. Loading of baskets has to be achieved in such a way that a hand can easily fit between the packs without compressing them.
- Foil/ paper packaging has to be stacked in an oblique way to each other, the paper side located downwards.
- Textiles have to be placed vertically whenever possible.
- Instruments shall be equally distributed among the tray
- Loading has to be performed with similar instruments or goods to be processed.
- Goods to be sterilized shall be packaged in a loose way. Sufficient room shall remain free in order that the steam can reach the surface of the goods to be sterilized without hindrance.
- Keep at least 5 cm (2") distance - horizontally and vertically - between the goods to be sterilized and the chamber walls.

5 Operating the sterilizer

5.1 Turn on the sterilizer

The **Sterilizer main power disconnect switch** is located on the right side above the sterilizer. This switch disconnects power to the control. Under normal operation, this switch is left in the **On** position at all times.

The **Steam shut-off valve** (refer to Figure 5-1) is located behind the access door on the top. Ensure this valve is in the open position before trying to operate the sterilizer.

The **Tap water shut-off valve** is located behind the lower maintenance door. Ensure this valve is open.

The **Deionized water shut-off valve** (option internal steam generator) is located behind the lower maintenance door on the left side. Ensure this valve is in the open position before trying to operate the sterilizer.

The **Key Switch** (refer to Figure 5-1) is located on the Operating End Control Panel. Turn the key in the **I** position (on)




Figure 5-1: Utility locations

5.2 Log on

To operate the sterilizer a user must be logged on. The operator has to enter his username and password (refer to Figure 5-2 to Figure 5-5).

Note:

The Belimed service personnel can configure the system that a 'common' user is automatically logged on. The 'common' user can select and start a cycle and can operate the door. User 'System' is entered in the batch log. In this configuration, it is not possible to trace back to the user who started the program

Press the 'menu' button in the out-of-cycle display (refer to Figure 5-2)

P1 chamber	T1 chamber	status	
14.78 psia	258.1 °F	standby	
more ..			
1 PreVac 270 long Dry		2 PreVac 270 short Dry	
		4 Flash 270 single Ins	
5 Flash 270 multi Ins			cycle start
PreVac 270 long Dry		open door	menu
cycle counter	0		

Figure 5-2: Out-of-cycle display

Then press the 'log on' button in the 'menu' display (refer to Figure 5-3)

P1 chamber	T1 chamber	status	door locked
990 mbar	35.7 °C	standby	
menu			
log on		info	readings
set-up		print last cycle	alarms
cycle intervention		batch data	user administration
←		maintenance	

Figure 5-3: Menu display

The log on display shows the entry-field for the password (refer to Figure 5-4).




Figure 5-4: Log on display

To enter a password, press the field below 'password'.

The alphanumeric keyboard is shown (refer to Figure 5-5):




Figure 5-5: Alphanumeric keyboard

Enter the password and press the ↵ button. After entering the Password, press the 'Login' button (Figure 5-4)

You will see message 'Invalid Password' if you have entered an incorrect Password.

If the Password is correct, the display switches to the Menu display.

5.3 Loading the sterilizer

WARNING - BURN HAZARD

Sterilizer and rack/shelves will be HOT after cycle is run. Always wear protective gloves and apron when removing a processed load. Protective gloves and apron should also be worn when reloading sterilizer following previous operation.

WARNING - FALL HAZARD

To prevent falls keep floors dry by immediately wiping up any liquids or condensation in sterilizer loading or unloading area.

Prepare the load as described in chapter 8 Techniques of sterilization.

CAUTION

Keep the movement area of the door unobstructed during door opening and closing.

5.3.1 Loading sterilizer

If you want to enter batch data manually select the 'Batch data' menu (refer to chapter 8.8.1).

If the optional documentation system 8535-BC is installed, enter batch data with the barcode scanner first before loading the sterilizer (refer to chapter 8.8.2).

- Open the chamber door
option motorized doors: by pressing the 'open door' button or actuating the foot pedal.
- Pull-out a shelf (refer to Figure 5-6)
- Put the load on the shelf
- Slowly push-in the shelf
- Close the chamber door
option motorized doors: by pressing the 'close door' button, until door is closed.



Figure 5-6: Pull-out shelf

5.3.1.1 Remove the upper shelf

If an item is higher than 300mm, you can remove the upper shelf.

To remove a shelf:

- Pull-out the shelf to the stop
- Lift-on the latch on the right side (refer to Figure 5-7)
- Pull-out and remove the shelf, use both hands (refer to Figure 5-8)




Figure 5-7: Pull-out shelf



Figure 5-8: Pull-out shelf

5.3.1.2 Remove the lower shelf

For cleaning the chamber and the strainer remove the lower shelf in the same way.

5.4 Select a cycle

5.4.1 Factory-set cycles and cycle values

The predefined programs of the Belimed steam sterilizers TOP5000 are described in Chapter 6.

5.4.2 Select a cycle

Ensure that the load is suitable for the intended sterilization temperature (temperature resistant) before selecting a cycle.

Five cycle selection buttons (P1..P5) are shown on the screen (refer to Figure 5-9). These buttons display the cycle number, the cycle name and the sterilization exposure temperature.

P1 chamber	T1 chamber	status	
14.78 psia	258.1 °F	standby	
more ..			
1 PreVac 270 long Dry	2 PreVac 270 short Dry		
		4 Flash 270 single Ins	
5 Flash 270 multi Ins			cycle start
PreVac 270 long Dry		open door	menu
cycle counter	0		

Figure 5-9: Cycle selection buttons with cycle start button

5.5 Start a cycle

The selected program is started after you press the Start button if:

- a cycle is selected
- chamber doors are closed
- no alarms are pending
- the minimum steam pressure is exceeded

P1 chamber	T1 chamber	status	door locked
13.91 psia	74.8 °F	standby	
more ..			
door not closed			
ok			
Instrumente 134°C			
cycle counter	7	open door	menu

Figure 5-10: Cycle selection with "dialog box"

If the sterilizer is not at start standby, the cause is displayed with a dialog box (Figure 5-10) when you press the Start button.

5.6 In-cycle

During a cycle the following information is displayed:

- Chamber pressure P1
- Chamber temperature T1
- Cycle status (phase)
- Door status
- Cycle graph
- Cycle intervention button
- Menu button




Figure 5-11: In-cycle display

5.7 Cycle failed

If there was an injurious alarm during a cycle, the door on the non operating end is locked.

On the operating end panel the message 'cycle failed' is displayed (refer to Figure 5-12).

If the user must press the 'ok' button. Only the door on the operating end can be opened. Unload the sterilizer.

The load must be wrapped again with dry material before sterilizing.




Figure 5-12: Cycle failed display

5.8 Unloading the sterilizer

WARNING - BURN HAZARD

Sterilizer and rack/shelves will be HOT after cycle is run. Always wear protective gloves and apron when removing a processed load. Protective gloves and apron should also be worn when reloading sterilizer following previous operation.

WARNING - FALL HAZARD

To prevent falls keep floors dry by immediately wiping up any liquids or condensation in sterilizer loading or unloading area.

The door on the non operating end (double-door sterilizers) can be opened if:

- The sterilization cycle is completed (Test cycles: door can be opened on the operating end only)
- The cycle passed (no injurious alarms during cycle)
- The difference between chamber pressure and atmospheric pressure does not exceed ± 80 mbar.

CAUTION:

Keep the movement area of the door unobstructed during door opening and closing.

5.8.1 Unloading sterilizer

- Open the chamber door by pressing the 'open door' button.
- Pull-out the shelf
- Remove the load
- Push-in the shelf
- Close the chamber door by pressing the 'close door' button

5.9 Cycle documentation

At cycle end the cycle documentation is printed on the built-in printer (refer to Figure 5-13).

The Panel Mount printer uses impact dot matrix technology.

For a detailed description of the printer see instructions in chapter 4.5.

For printout examples see in chapter 4.6.



Figure 5-13: Built-in printer

5.10 Log off

After using the sterilizer the user should log off to prevent that unauthorized persons can handle the sterilizer. Press first the 'menu' button on touch screen and then press the 'log off' button in the 'menu' display (refer to Figure 5-14).




Figure 5-14: Menu display

5.11 Turn off the sterilizer

After unloading the last cycle of a day turn off the sterilizer with the key switch (refer to Figure 5-1).

6 Sterilizer Factory Cycles Settings

Belimed steam sterilizers TOP5000 are shipped with the factory-set cycles and cycle values listed in Table 6-1:

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
1	PreVac 270F 4S / 30Dry	270°F	4 minutes	30 minutes	Double-wrapped instrument trays, maximum weight of 17 lbs each. Fabric Packs.
2	PreVac 270F 4S / 5Dry	270°F	4 minutes	6 minutes	Double-wrapped instrument trays, maximum weight of 17 lbs each. Fabric Packs. Instruments for immediate use only.
3	Flash 270F - Single Ins	270°F	3 minutes	1 minute	Unwrapped instrument tray with a single instrument
4	Flash 270F - Multi Ins	270°F	10 minutes	1 minute	Unwrapped instrument tray with non porous multiple instruments (maximum weight of 17lbs)
5	Express 270F	270°F	4 minutes	3 minutes	Single wrapped instrument tray with non porous single instrument
6	Bowie-Dick Test	273°F	3.5 minutes	1 minute	One DART or Bowie-Dick-Test-Pack
7	Leak Test	-	-	-	Empty chamber
8	Warm up & Leak Test	270°F	4 minutes	3 minutes	Empty chamber
9	-	-	-	-	-
10	Liquid 250F	250°F	45 minutes	-	Max 1000 ml of Liquid per bottle. Use only vented bottles, type I, borosilicate glass. Cycle disabled

Table 6-1 Factory set cycles

6.1 270°F Pre-vacuum Cycles

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
1	PreVac 270 long Dry	270°F	4 minutes	30 minutes	Double-wrapped instrument trays, maximum weight of 17 lbs each. Fabric Packs.
2	PreVac 270 short Dry	270°F	4 minutes	6 minutes	Fabric Packs.

Table 6-2: 270°F Pre-vacuum Cycles

The cycle graph provides a visual representation of the Belimed Steam Sterilizer TOP 5000 Pre-vacuum cycles (refer to Figure 6-1):




Figure 6-1: Cycle Graph - 270°F Pre-vacuum cycles

6.2 250°F Liquid Cycle

RESTRICTION:

It is inappropriate for a health care facility to sterilize liquids for contact.

WARNING - EXPLOSION HAZARD:

This sterilizer is not designed to process flammable compounds

WARNING - BURN HAZARD:

When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedures:

Use Liquid cycle only; no other cycle is safe for processing liquids

Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.

Use only Type I borosilicate glass bottles, do not use ordinary glass bottles or any container not designed for sterilization.

Do not allow hot bottles to be jolted; this can cause hot-bottle explosions. Do not move bottles if any boiling or bubbling is present.

This sterilizer is not designed to process flammable liquids nor liquids in containers that are not designed for sterilization.

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
10	LIQUID 250°F	250°F	45 minutes	--	Max 1000ml of Liquid per bottle. Use only vented Type I borosilicate glass bottles.

Table 6-3: 250°F Liquid Cycle

The cycle graph provides a visual representation of the Belimed Steam Sterilizer TOP 5000 LIQUID cycles (refer to Figure 6-2):




Figure 6-2: Cycle Graph - Liquid cycle

6.3 270°F Express Cycle

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
5	Express 270	270°F	4 minutes	3 minutes	Single Wrapped Instrument Tray with non porous single instrument

Table 6-4: 270°F Express Cycle

The cycle graph provides a visual representation of the Belimed Steam Sterilizer TOP 5000 Express cycles (refer to Figure 6-3):




Figure 6-3: Cycle Graph - Express cycle

6.4 270°F Flash Cycles

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
3	Flash 270 single Ins	270°F	3 minutes	1 minute	Unwrapped instrument tray with a single instrument
4	Flash 270 multi ins	270°F	10 minutes	1 minute	Unwrapped instrument tray with non porous multiple instruments (maximum weight of 17lbs)

Table 6-5: 270°F Flash Cycles

The cycle graph provides a visual representation of the Belimed Steam Sterilizer TOP 5000 Flash cycles (refer to Figure 6-4):




Figure 6-4: Cycle Graph - Flash cycles

6.5 DART / Bowie-Dick Test cycle

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
6	DART / Bowie-Dick Test	273°F	3.5 minutes	1 minute	One DART or Bowie-Dick Test Pack

Table 6-6: DART / Bowie-Dick Test cycle

The cycle graph provides a visual representation of the Belimed Steam Sterilizer TOP 5000 DART / Bowie-Dick Test cycle (refer to Figure 6-5):




Figure 6-5: Cycle Graph - DART / Bowie-Dick Test cycle

CAUTION

The Bowie-Dick/ DART Test Cycle has to be performed daily.

6.6 Warm-up & Leak Test cycle

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
8	Warm up & Leak test	270°F	3 minutes	3 minutes	EMPTY CHAMBER

Table 6-7: Warm-up & Leak Test cycle

The cycle graph provides a visual representation of the Belimed Steam Sterilizer TOP 5000 Warm-up & Leak Test cycle (refer to Figure 6-6):




Figure 6-6: Cycle Graph - Warm-up & Leak Test cycle

CAUTION

In order to ensure the proper function of the sterilizer, we recommend to perform the leak-rate-test at least once per week, at normal working temperature.

7 Techniques of sterilization

7.1 Guidelines for preparing and sterilizing wrapped packs

WARNING - BURN HAZARD

BURN HAZARD: Sterilizer and rack/shelves will be HOT after cycle is run. Always wear protective gloves and apron when removing a processed load. Protective gloves and apron should also be worn when reloading sterilizer following previous operation

7.1.1 Preparing load:

- Before sterilization, materials must be thoroughly cleaned.
- Wrappers may be made of 100% cotton, 140 thread count, two-ply fabric, and freshly laundered.
- Reusable wrappers should be laundered between sterilization cycles
- Before wrapping, packaging materials must be dry and should be under room temperature (15°C-25°C) and at a relative humidity ranging from 30%-60% for at least 2 hours.

7.1.2 Wrapped fabric packs

- Place two wrappers on work surface
- Place contents on wrappers. Fabrics should be folded flat, with each succeeding layer placed crosswise to the one below to promote free steam circulation.
- Place internal chemical indicator in the center of the pack
- Wrap contents sequentially in two wrappers (refer to Figure 7-1).
- Secure with sterilizer tape and identify
- The maximum weight must not exceed 8 kg , density factor: not in excess of 0.47kg/dm³, and the size must be limited to 250x500x250mm.
- Single Fabric Test Pack (for validation) should weigh approximately 7 kg and the size must be limited to 220x300x250mm.

Notes:

- Wrapper size should be adequate for the desired method of wrapping. Excessive size wrappers may cause drying problems
- Do not wrap too tightly




Figure 7-1: Double-wrapping: Square fold

7.1.3 Wrapped instrument sets

- Inspect instruments to make sure they are functioning properly, they are clean and dry.
- Open, disassemble or unlock instruments to permit steam contact to all surfaces.
- Use mesh-bottom tray or wire-mesh or equivalent trays.
- Place a fully opened huck towel in the bottom of the tray. This will assist drying.
- Place instruments on towel, distributing them as much as possible. Using an adequate size tray will allow optimal distribution of the instruments.
- Fold the towel excess over the instruments
- Place internal chemical indicator on instruments. Ink side should not come in contact with metal surface.
- Wrap instrument tray sequentially in two wrappers (refer to Figure 7-1).
- Secure with sterilizer tape and identify.

- Recommended weight of each wrapped instrument set is 5-8 kg for minimizing moisture retention. The maximum weight must not exceed 10kg.

Notes:

- Instrument trays must be designed for effective sterilization and drying
- Use a towel which covers the tray with minimum excess overhang.
- Wrapper size should be adequate for the desired method of wrapping. Excessive size wrappers may cause drying problems.

7.1.4 Wrapped utensils

- When placing utensils in a set, separate each clean, dry basin from the one beneath it by a huck towel.
- Open towel to fully cover the metal surface.
- Arrange utensils so that the bottom of each is parallel to the one beneath it. This allows air to escape from the utensils and helps drying.
- Place internal chemical indicator in an area of the package where steam penetration is worst. Ink side should not come in contact with metal surface.
- Wrap utensils sequentially in two wrappers (refer to Figure 7-2).
- Secure with sterilizer tape and identify
- Recommended weight of each wrapped utensil set is 6-8 kg for minimizing moisture retention.

Note:

- Wrapper size should be adequate for the desired method of wrapping. Excessive size wrappers may cause drying problems.




Figure 7-2: Double-wrapping: Envelope fold

7.1.5 Loading the shelves:

- Do not combine fabrics and hard goods (instruments) in one load
- Place fabric packs on edge to promote passage of steam through the pack
- Do not overload shelves
- Do not compress packages
- Provide at least 2 inches between the chamber walls and the packages
- Place packages on sterilizer shelves, do not staple packages

7.1.6 Steam sterilization cycles

Cycle no. 1 and 2 is recommended for full load fabric packs, instrument trays, utensils.

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
1	PreVac 270 long Dry	270°F	4 minutes	30 minutes	Double-Wrapped Instrument Trays or Fabric Packs, For amount refer to <i>Table 7-2</i>
2	PreVac 270 short Dry	270°F	4 minutes	6 minutes	Fabric Packs, For amount refer to <i>Table 7-2</i>

Table 7-1: Recommended solid cycles

Do not sterilize a mix load of instrument trays and fabric packs.

Model	Sterilizer Chamber Size	Wrapped Instrument Trays 20"x10" max. 17lb each	Wrapped Instrument Trays 17"x12" max. 17lb each	Fabric Packs 11"x11"x9" max. 6.6lb each	Fabric Packs 9"x9"x6" max. 3.3lb each
5-5-9 VS1 5-5-9 VS2	21"x 21" x 38" (535 x 535 x 965) mm	2	4	6	12

Table 7-2: Recommended solid load

Ensure that the load is suitable for the intended sterilization temperature (temperature resistant) before selecting a cycle.

7.1.7 Material qualification

Fabrics:

Use only fabrics designed for steam sterilization with 134°C.

Life time of fabrics depend mainly on washing and ironing the fabrics.

Please note the manufacturer's declaration.

Instruments:

Use only instruments designed for steam sterilization with 134°C.

Instrument life time depends mainly on use and cleaning with chemicals.

100% ASTM 316 Stainless steel instruments can generally be used for steam sterilization.

For other materials note the manufacturer's declaration.

Other materials:

Before sterilizing other materials refer to the manufacturer's recommendations and restrictions (steam, temperature, max. sterilization time).

7.2 Guidelines for preparing and sterilizing liquids

WARNING - EXPLOSION HAZARD:

This sterilizer is not designed to process flammable compounds

WARNING - BURN HAZARD:

When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedures:

Use Liquid cycle only; no other cycle is safe for processing liquids.

Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.

Use only Type I borosilicate glass bottles, do not use ordinary glass bottles or any container not designed for sterilization.

Do not allow hot bottles to be jolted; this can cause hot-bottle explosions. Do not move bottles if any boiling or bubbling is present.

Avoid sudden opening of door at end of cycle. Wait at least 10 minutes before door opening and unloading the sterilizer.

Sterilizer and rack/shelves will be HOT after cycle is run. Always wear protective gloves, apron and face shield when removing a processed load. Protective gloves and apron should also be worn when reloading sterilizer following previous operation.

WARNING - Restrictions

It is inappropriate for a health care facility to sterilize liquids for direct patient contact.

7.2.1 Preparing load

- Use only Type I borosilicate glass bottles, do not use ordinary glass bottles or any container not designed for steam sterilization
- The maximum allowed bottle size is 1000ml
- Do not use flammable liquids
- Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.
- Place the bottles into baskets to prevent falling down during transportation
- Place chemical indicator in a hermetically sealed vial in the center of a bottle.




Figure 7-3: Loading shelves with liquid load

7.2.2 Liquid cycle

Only use cycle no. 6 liquid cycle

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
10	Liquid 250	250°F	45 minutes	--	Max 1000ml of Liquid per bottle. Use only vented Type I borosilicate glass bottles. For amount refer to <i>Table 7-4</i>

Table 7-3: 250°F Liquid Cycle

Model	Chamber Size (mm)	Volume of Liquid in one bottle	Number of bottles
5-5-9 VS1	21"x 21" x 38"	1000 ml	36
5-5-9 VS2	(535 x 535 x 965) mm		

Table 7-4: Recommended liquid loads

This sterilizer is not designed to process flammable liquids nor liquids in containers that are not designed for sterilization.

7.3 Appropriated use of Flash cycle

The express cycle is used to sterilize a single unwrapped instrument tray **for immediate** use (e.g. a dropped instrument). There is **no storage or shelf life** of Flash sterilized items.

This cycle may only be uses if the sterilizers non operating end is within the surgery theatre.

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
3	Flash 270 single Ins	270°F	3 minutes	1 minute	Unwrapped Instrument Tray with a single instrument
4	Flash 270 multi Ins	270°F	10 minutes	1 minute	Unwrapped Instrument Tray with non porous multiple instruments (maximum weight of 17lbs)

Table 7-5: 270°F Flash Cycles

Note: The Belimed Steam Sterilizer TOP 5000 is usually installed in the Central Sterilization Supply Department CSSD.

7.4 Appropriated use of Express cycle

The express cycle is used to sterilize a **single wrapped** instrument tray with a non porous single instrument for **immediate** use. The Express cycle is useful in providing quick turnaround of an instrument using the wrapped technique for transport from sterilizer to the point of use. A single wrapped item sterilized with the express cycle does not have a shelf life.

No.	Cycle	Sterilize temp	Sterilize time	Dry time	Recommended load
5	Express 270	270°F	4 minutes	3 minutes	Single Wrapped Instrument Tray with non porous single instrument

Table 7-6: 270°F Express Cycle

7.5 Guideline for determining wet packs

7.5.1 Introduction

Aseptic techniques for surgical procedures require that all supplies coming in contact with the surgical field have to be sterile.

An inherent, inseparable quality for sterility of supplies is a 'state of dryness'. Wet materials transmit bacteria; therefore, a 'state of wetness' would compromise the sterility of process packs and instruments presented to sterile field.

7.5.2 Evaluation of wet packs

A load must be examined for wet packs for three conditions:

- Water droplets on the exterior of a pack
- Water droplets within a pack
- Absorbed moisture in a pack

7.5.3 Moisture retention acceptance criteria

The AAMI ST8: 2001 sets the following acceptance criteria:

7.5.3.1 Fabric packs

Moisture retained by the fabric test pack must cause no more than 3% increase in pre-sterilization test pack weight, and the pack must not exhibit wet spots.

Wrapped instruments

Upon completion of the recommended cycle, the wrapped instrument pack must have no wet spots on the outer wrappers. Moisture retained in the load must cause no more than a 20% increase in pre-sterilization weight of the towel.

7.5.3.2 Summary

- External droplets or visible moisture on the exterior pack, or on the tape, are unacceptable unless the wrap is completely impermeable to water.

Note: This should be inspected while unloading the sterilizer

- Water droplets on the interior of a wrap, or on the items within the pack, are unacceptable
- A pack is unacceptable if the pack is damp or wet when opened for use.
- A general guideline is that the pack be completely dry after cooling to room temperature (i.e. 21°C and 50% relative humidity) for a minimum of one hour following unloading from the sterilizer.

7.5.4 Solving wet pack problems

If moisture is detected on the interior of the wrap, in the huck towel or on the item the causes can be:

- Insufficient drying time
- Chamber overloaded; size, density or weight of specific packs
- Excessively large hard goods items being sterilized
- No absorbent huck towel used
- Materials not being permitted to cool and equilibrate to room temperature
- Wet steam supply
- Defective vacuum drying system
- Air leak

7.6 Sterility maintenance

7.6.1 Cooling

All items removed from the sterilizer should remain on transfer carriage until adequately cooled.

- They should not be touched within cooling process
- A cooling time of 30-90 minutes is recommended
- During cooling the load should be placed in a low-traffic area where there are no air-conditioning or other cold-air vents in close proximity

7.6.2 Inspection

As items are removed from the sterilizer, they should be visually inspected. Do not use:

- Any items with torn packaging
- Packaging that appears to be wet

7.6.3 Sterile storage

Sterile materials should be stored:

- at least 8-10 Inches above floor
- at least 18 inches below ceiling
- at least 2 inches from outside walls

Items should be positioned so that packaging:

- is not crushed
- is not bent
- is not compressed
- can't become wet

Closed or covered cabinets are recommended for the storage of seldom-used supplies.

7.6.4 Expiration dating

Each item should be labeled with a control date for stock rotation and the following statement:

'Contents sterile unless package is open or damaged. Please check before using'.

If the product contains material that degrades over time, the product package should be labeled with a clearly identifiable expiration date that takes this degradation into account and/or that is based on the device manufacturer's instructions.

8 Operating device TOP5000

8.1 User administration

Each user is defined by a:

- user name (alphanumeric, max. 20 characters)
- user access level (numeric 1..4, for user rights see chapter 9.2)
- Password (alphanumeric, max.10 characters)

8.1.1 User's rights

Users with four different access levels can be defined on the TOP5000:

Level 1: Operator

Level 2: Group Leader

Level 3: Service

Level 4: User Admin


The Administrator only can define new users or delete users.

Additionally one user (level 5) is reserved for Belimed Service and one user (level 9) is reserved for the Belimed Administrator.

Important

In some cases, only buttons for which the operator has read or write access are displayed.

Table 8-1 shows the TOP5000 menus with the corresponding user rights.






Table 8-1: User menus and rights

8.1.2 Direct login

In the case of visible buttons for which the user has no authorization, the user is prompted directly to enter the user name and the Password. The user can log in with his or her Password and then receives access authorization.




Figure 8-1: Password entry

8.1.3 Administrate user

The Administrator (Level 4) has permission to create new users, delete users or change Passwords. Press the 'Menu' button. The Mask menu is then displayed (refer to Figure 8-2).




Figure 8-2: Menu display

Press the 'User Administration' button (Figure 8-3). A list of all users with Password, Group and Level is then displayed. Only users with Level (Group) equal to or lower than the user logged in are displayed. Additionally a automatic logoff time (in minutes) can be defined, where 0 means no automatic logoff.




Figure 8-3: User administration

8.1.3.1 Create New User

A new user can be defined by making entries in fields "User", "Password" and "Group". The user name may contain maximum 20 characters and the Password may contain maximum 10 characters. User names and Passwords may never be the same. Enter the automatic logoff time in minutes, if required.

8.1.3.2 Delete User

Fields "User" and "Password" must be deleted in order to delete a user from the list. Button "BSP" can be used for this.

8.1.3.3 Change Password

The Password can be changed easily in field "Password".

8.2 Info

Each user can open the information display.

Press 'menu' button. The menu functions are displayed (refer to Figure 8-4):

P1 chamber	T1 chamber	status
886 mbar	85.3 °C	standby
menu		
log off	info	readings
set-up	print last cycle	alarms
cycle intervention	batch data	user administration
		maintenance

Figure 8-4: Menu display

Press 'info' button

The display shows the following information (refer to Figure 8-5):

- actual time and date
- User name
- Time since cycle start
- Remaining time in the actual phase
- Date of last Leak Test and last Bowie Dick Test
- Cycle counter
- Status of doors and steam generator (if configured)

P1 chamber	T1 chamber	status	door locked
983 mbar	121.9 °C	standby	
date / time	dd/mm/yy h:m:s	07/02/2006 13:49:30	
user name	Belimed_Service		
time since cycle start	0.0 min		
end of phase in	0.0 min		
last Leak Test	01/01/1990		
last Bowie-Dick Test	07/02/2006		
cycle counter	62	door side 1	closed
steam gen.	off	door side 2	closed

Figure 8-5: Info display

8.3 Readings

Each user can open the Readings display.

Press 'menu' button. The menu functions are displayed (refer to Figure 8-6):




Figure 8-6: Menu display

Press 'readings' button.

The display shows the actual readings of all temperature- and pressure sensors, used in the sterilizer (refer to Figure 8-7):




Figure 8-7: Readings display

- | | | | |
|----|--|----|-------------------------------------|
| P1 | Chamber pressure - control | P2 | Chamber pressure - documentation |
| P3 | Jacket pressure | P4 | Steam supply pressure |
| T1 | Chamber temperature - control | T2 | Chamber temperature - documentation |
| T4 | Condenser temperature - chamber outlet cooling | | |
| T5 | Circulation tank temperature - vacuum pump cooling | | |

8.4 Set-up

Press 'menu' button. The menu functions are displayed (refer to Figure 8-8):




Figure 8-8: Menu display

Press 'set-up' button.

The display shows the 'set-up' functions (refer to Figure 8-9):

- Early start and night-time shutdown (access level 1)
- Set clock (access level 2)
- Select Language (access level 1)
- Steam generator on/off (access level 1 / if steam generator connected)
- Drain steam generator (access level 2 / option: internal steam generator)
- Configure cycle selection (access level 2)
- Set screen contrast (access level 1)
- Activate on-line documentation




Figure 8-9: Set-up display

8.4.1 Early start / night-time shutdown

Press 'early start / night shutd.' button. The display shows the input-mask (refer to Figure 8-10):



Figure 8-10: Early start / night-time shutdown display

8.4.1.1 Night-time shutdown

If you press '**night-time shutdown**' button, the sterilizer goes into a sleep mode at the end of the next cycle. All valves are closed, an optional connected steam generator is turned off.

To restart the sterilizer turn key switch off and on again.

8.4.1.2 Early start

Conditions for early start:

The Sterilizer must be in the initial position and the chamber doors must be closed.

If you want to start a cycle early in the morning (usually the warm-up and leak test cycle), set '**early start**' date and time, enter '**prog. no.**' and activate '**early start**' button (left side). With the keys + and - the starting date can be selected gradually.

The sterilizer goes into a sleep mode (the display is turned off). At the pre-set early start time the sterilizer goes into an active mode, the operating panel display is turned on and the pre-selected cycle starts. The system suggests the date of the next day and the last pre-set time.

If the sterilizer is in sleep mode (early start or night-time shutdown) and you want to operate, turn the key-switch off and on again. The early start and night-time shutdown are turned off.

8.4.2 Set clock

Usually twice a year (change of winter- to summertime and summer- to wintertime) the time has to be adjusted.

Press '**Set clock**' button (refer to Figure 8-9).

Press the date/time button and enter the time and date in the format: MM/DD/YY hh:mm:ss

Legend:

MM: Month (1..12) DD: Day of the month (1..31) YY: Year (00..99)
hh: hour 1..12 mm: minute (0..59) ss: second (0..59) xx: AM or PM

Press '**save date/time**' button (refer to Figure 8-11) to set the new date/time




Figure 8-11: Set clock display

8.4.3 Language and units

Each operator (access level 1) can select the language he wants to have displayed.

Press 'language' button in the set-up menu and then press the operation language select button '▼' (refer to Figure 8-12). All predefined languages are displayed. Press the language you need in the display field. The display changes to the selected language.

The administrator (access level 4) has the right to change the units:

Temperature:

- Degree Celsius (°C)
- Degree Fahrenheit (°F)

Press the temp. unit select button '▼' and choose the temperature unit.

Pressure:

- Millibar absolute (mbar)
- Bar gauge (bar)
- Kilopascal absolute (kPa)
- Pounds/square inch absolute (psia)
- Pounds/square inch gauge for pressure (psig) and inches of Mercury for vacuum (in.Hg)

Press the pressure unit select button '▼' and choose the pressure unit.



Figure 8-12: Language and units display

8.4.4 Steam generator

If the steam generator can be activated from the sterilizer (option) you can turn on and off the steam generator by pressing the 'steam generator' button (refer to Figure 8-14).

8.4.5 Drain steam generator

If an internal steam generator (option) is installed, it can be emptied by pressing the 'drain steam generator' button (refer to Figure 8-14).

8.4.6 Configure cycle selection

The administrator (level 4) can define:

- The selectable cycles by a user can be activated or locked.
- Cycle start only after entering a password by the operator (set: yes)
- Automatic logon with user name 'common' using level 1 (set: yes).

To activate a cycle or a function or press corresponding button (refer to Figure 8-13).




Figure 8-13: Configure cycle selection

8.4.7 Display contrast

All logged on users can adjust the display contrast (Refer to Figure 8-14).

To increase the contrast of the display, press 'contrast +' button.

To decrease the contrast of the display, press 'contrast -' button.




Figure 8-14: Set-up display

8.4.8 Online docu

If there is an external PC- based documentation system 8535 installed, the built-in printer can be turned off (refer to Figure 8-14). In case of a breakdown of the system 8535, the built in printer can be activated by the group leader.

8.5 Print last cycle

The cycle documentation is printed at the end of the cycle.

It is possible to print a copy of the cycle documentation until the next cycle is started.

Press 'menu' button. The menu functions are displayed (refer to Figure 8-15):

Press 'print last cycle' button and a cycle documentation copy is printed.




Figure 8-15: Menu print last cycle display

8.6 Alarm messages

All important process parameters are continuously monitored. If there is a deviation, an alarm button is displayed in the main mask.(refer to Figure 8-17).




Figure 8-16: Out-of-cycle display with active alarm




Figure 8-17: In-cycle display with active alarm

Press 'alarm' button. All active alarms are displayed (refer to Figure 8-18).

The alarm mask has 3 buttons:

- Help (□)
- ACK
- history

The alarm is displayed with number/time/date and the alarm message.

Each alarm must be acknowledged by an operator (access level 1)

- To acknowledge an alarm press the 'ACK' button.
- To display the history of the last 250 alarms press 'history' button
- To get more information about the active alarm press Help button 'H'



Figure 8-18: Alarm display

The operator should follow the Help instructions on the alarm Display. If these instructions fail to clear the alarm, consult your department supervisor or a trained service technician before using the sterilizer further.

8.6.1 Critical alarms

Only the critical alarms are printed out in the batch report. If a critical alarm has occurred during a cycle, the door can only be released on the operating end. A message 'cycle failed' is displayed on the operating panel and printed on the batch documentation. The start of a new cycle is not possible if a critical alarm is pending.

8.6.2 Alarm list

Table 8-2 shows the list of alarms and the possible causes.

No.	1st line: Alarm message text 2nd line: Cause	Critical alarm
1	Emergency-Stop Emergency-Stop button on the sterilizer has been pressed.	X
2	Vacuum pump Motor protection switch vacuum pump has tripped.	X
3	Door motor Motor protection switch of the door motors has tripped.	
4	Cycle aborted Function 'cycle abort' has been activated by an operator or by the system itself.	X
5	Power failure Loss of electric power	X
6	Manual advance The 'manual advance' button has been pressed by on operator	X
7	Fault of PLC battery PLC battery empty (do not switch-off electric power switch); change battery	

No.	1st line: Alarm message text 2nd line: Cause	Critical alarm
8	On-Off Key switch off The key switch was turned off during cycle running	X
9	Too long in step The monitoring time of a cycle phase has been exceeded.	X
11	Sensor T1 failed Temperature sensor T1 signal is out of range (0..150°C)	X
12	Sensor T2 failed Temperature sensor T2 signal is out of range (0..150°C)	X
14	Sensor T4 failed Temperature sensor T4 signal is out of range (0..150°C)	X
15	Sensor T5 failed Temperature sensor T5 signal is out of range (0..150°C)	X
18	Sensor P1 failed Pressure sensor P1 signal is out of range (4..20mA)	X
19	Sensor P2 failed Pressure sensor P2 signal is out of range (4..20mA)	X
20	Sensor P3 failed Pressure sensor P3 signal is out of range (4..20mA)	X
21	Sensor P4 failed Pressure sensor P4 signal is out of range (4..20mA)	
22	Undertemperature T1 Chamber temperature T1 is below alarm set-point during sterilization phase	X
25	Overtemperature T1 Chamber temperature T1 is above max. alarm set-point during sterilization phase	X
27	Overtemperature T5 Temperature in circulation tank is >50°C	X
28	Divergence T1-T2 During sterilization phase the difference T1-T2 is >±1°C	X
30	Max. chamber pressure Maximum chamber pressure exceeded, steam to chamber valve not closed / leaky	X
31	Max. jacket pressure Maximum jacket pressure exceeded, steam to jacket valve not closed / leaky	X
32	Low jacket pressure Jacket set pressure not reached within pre-set time, steam to jacket valve not open or low steam supply pressure	X

No.	1st line: Alarm message text 2nd line: Cause	Critical alarm
33	Steam pressure low Steam supply pressure is below set-point	
34	Correlation low temp. Temp. T1 is lower than the saturated steam temperature calculated from P1	X
35	Correlation high temp. Temp. T1 is higher than the saturated steam temperature calculated from P1	X
36	Leak rate test failed The chamber is leak or wet	X
37	Divergence P1-P2 During sterilization phase the difference P1-P2 is $> \pm 100\text{mbar}$	X
38	Door 1 time overrun Door on the operating end hasn't reached the end position within 20 seconds	
39	Door 2 time overrun Door on non operating end hasn't reached the end position within 20 seconds	
42	Door 1 not closed Limit switch door closed on operating end is not actuated during a cycle	X
43	Door 2 not closed Limit switch door closed on non operating end is not actuated during a cycle	X
48	Door not unlocked Door seal didn't unlock during door- unlock step	
49	Door not locked Minimum door seal pressure not reached	X
50	Ext. steam generator Alarm signal from external steam generator	
58	Atmospheric pressure monitoring Atmospheric pressure switch difference to sensor P1	
59	Deviation from time base The real time clock deviates from the PLC time-base	X
61	Overtemperature T4 Temperature in condenser is $>80^{\circ}\text{C}$	X
62	Lack of compressed air Compressed air pressure is $<4\text{ bar g}$	X
76	Connection to printer The PLC doesn't receive telegrams form printer interface	
78	No connection to system 8535 The PLC doesn't receive telegrams form documentation system 8535	

Table 8-2: Alarm list

8.7 Cycle intervention

The cycle intervention touch screen button is used to

- Abort a cycle before it is finished normally (by operator, access level 1)
- Advance the cycle into the next phase (for maintenance and qualification, access level 2) or
- Continue an interrupted cycle (by emergency stop or motor starter failure, access level 1)

Press 'menu' button. The menu functions are displayed (refer to Figure 8-19):




Figure 8-19: Menu display

Press 'cycle intervention' button (refer to Figure 8-20).




Figure 8-20: Cycle intervention display

8.7.1 Cycle advance

Press 'cycle advance' button to advance the cycle into the next phase. An alarm message 'manual advance' is generated.

8.7.2 Cycle abort

Press 'cycle abort' button. This causes the sterilizer-chamber to depressurise (if pressurised), evacuate the chamber (to dry) and to air break. The door opening deactivates, the control prompts the operator to acknowledge the 'cycle aborted' message and to Log on.

If a cycle has been aborted, the door on the non operating is locked. The door on the operating end can be opened. The material is not sterile. Unload the sterilizer. The load must be wrapped again with dry material before sterilizing.

8.7.3 Cycle continue

Press 'cycle continue' to continue an interrupted cycle. This function only works if the alarm has already gone.

8.8 Batch data

8.8.1 Manual data entry

Batch data entry must be activated in the set-up menu.

Press 'menu' button and then 'batch data' button (refer to Figure 8-19):

The following batch data can be entered (refer to Figure 8-21):

Load details (alphanumeric 12 characters).

batch No. (12 digits).

Batch data is printed on the cycle documentation printer.




Figure 8-21: Manual batch data entry

8.8.2 Batch data entry with Documentation system 8535-BC

An optional PC- based cycle documentation system 8535-BC can be supplied for electronic batch data acquisition.

Before loading the sterilizer, the load ID- numbers are to be read with a barcode scanner.

A maximum of 18 items can be entered. The ID's are displayed on the display (refer to Figure 8-22).

The bar code type EAN13 is used, whereby 12 characters are displayed. The 13th character is the check sum and is not displayed. If necessary, the last or all entries can be deleted. At cycle start the data is transferred to the system 8535 and the table on the sterilizer is deleted afterwards.




Figure 8-22: Scanner based batch data entry

At cycle start the data is transferred to the system 8535 and the table on the sterilizer is deleted afterwards. The data can also be transferred several times to system 8535 with the "transfer batch data" button even before program start

8.8.3 Batch data input via scanner with ICS 8665 Documentation System

The data is **not** transferred automatically to system 8565 on program start. The data must be transferred to system 8565 before program start with the "transfer batch data" button. The program is then selected automatically by system 8565. This is indicated in a dialog box. After this, no other program may be selected.




Figure 8-23: Automatic program selection by ICS 8565

8.9 Maintenance messages

After a predefined number of cycles a maintenance message is displayed in the Out-of-cycle display.

P1 chamber	T1 chamber	status	door locked
14.75 psia	258.1 °F	standby	
more ...		maintenance	alarm
1 PreVac 270 long Dry		2 PreVac 270 short Dry	
		4 Flash 270 single Ins	
5 Flash 270 multi Ins		cycle start	
PreVac 270 long Dry		open door	menu
cycle counter	1		

Figure 8-24: Out-of-cycle display with maintenance message.

The maintenance message can be acknowledged by maintenance personnel (access level 3).

Press 'maintenance' button in the out-of cycle mask. If the maintenance has been done press 'ACK' button in the maintenance mask (refer to Figure 8-25).

P1 chamber	T1 chamber	status
1029 mbar	30.4 °C	standby
<p>maintenance level: 1</p> <p>Maintenance level required</p> <p>Refer to routine maintenance manual. Acknowledgement only by maintenance personal possible.</p>		
←		confirm

Figure 8-25: Maintenance message display.

8.10 Maintenance functions

Press 'menu' button. The menu functions are displayed (refer to Figure 8-26):

P1 chamber	T1 chamber	status												
886 mbar	85.3 °C	standby												
<p>menu</p> <table border="1"> <tr> <td>log off</td> <td>info</td> <td>readings</td> </tr> <tr> <td>set-up</td> <td>print last cycle</td> <td>alarms</td> </tr> <tr> <td>cycle intervention</td> <td>batch data</td> <td>user administration</td> </tr> <tr> <td colspan="2">←</td> <td>maintenance</td> </tr> </table>			log off	info	readings	set-up	print last cycle	alarms	cycle intervention	batch data	user administration	←		maintenance
log off	info	readings												
set-up	print last cycle	alarms												
cycle intervention	batch data	user administration												
←		maintenance												

Figure 8-26: Menu display

Press 'maintenance' button. The following functions can be selected (refer to Figure 8-27):

- Replace door seal at operating end (access level 3)
- Calibrate touch screen (access level 3)
- Calibrate sensors (access level 5)
- Machine set-up (access level 5)
- Maintenance intervals (access level 5)
- Printer self test (access level 5)
- Cycle parameters (access level 5)
- PLC functions (access level 5)
- I/O Test (Level 5)
- Air Detector (option, access level 5)




Figure 8-27: Maintenance menu display

8.10.1 Replace door seal

The maintenance personnel can replace the door seal. on the operating end by pressing 'replace door seal' button while the chamber door is open (refer to Figure 8-28).

Press 'blow out door seal' button to remove door seal on the operating end.

Press 'pull back door seal' button and the vacuum pump sucks a vacuum behind the door seal.

Press 'standby' button to set the sterilizer to initial (load) position

Press 'complete' button to set the sterilizer to 'complete' (unload) position to replace the door seal on the non operation end.

For replacing the door seal refer to 'maintenance manual' chapter 1.7




Figure 8-28: Replace door seal display

8.10.2 Calibrate touch screen

If the touch screen adjustment is inaccurate, adjust the touch screen.

Press 'calibrate touch screen' button and follow the instructions (refer to Figure 8-27).

8.10.3 Sensors adjustment

See Maintenance guidelines

8.10.4 Machine set-up

The Belimed maintenance personnel must set the machine set-up before running the first cycle.

The first display (refer to Figure 8-29) must be adjusted:

- Enter hospital name into the 'hospital name' field .
- Enter department name into the 'department' field.
- Enter the machine number (serial number) into the 'machine No.' field.
- Enter the ambient pressure (in mbar absolute) into the 'norm ambient pressure' field.

The ambient (atmospheric) pressure depends on the altitude of the sterilizer. Enter the following values:

Altitude (m)	ambient pressure (mbar)	Altitude (m)	ambient pressure (mbar)
0	1013	1200	877
200	989	1400	856
400	966	1600	835
600	943	1800	815
800	921	2000	795
1000	898	2200	775

Table 9-3: Ambient pressure for altitudes up to 2200m.




Figure 8-29: Machine set-up display 1

The second display (refer to Figure 8-30) functions are:

- The low steam pressure alarm delay time can be set (default: 300s)
- Pressure difference between set steam pressure and alarm pressure (default: 500mbar)
- Steam set pressure in out-of cycle position (default: 3900mbar)
- Steam set pressure in- cycle (default: 700mbar above chamber pressure during sterilize phase)




Figure 8-30: Machine set-up display 2

The third display (refer to Figure 8-31) functions are:

- The compensation time before the sterilization is started can be set (default: 4 seconds).
 - The control temperature of the circulation tank (default: 40/40°C*)
 - the control temperature of the circulation tank during vacuum steps (prevac, dry, default:21/30°C*)
 - The control temperature of the condenser if the chamber pressure is >500mbar (default: 50/30°C*)
 - The control temperature of the condenser if the chamber pressure is <500mbar (default: 35/15°C*)
- *) without / with cooling water loop.




Figure 8-31: Machine set-up display 3

The fourth display (refer to Figure 8-32) functions are:

- The system can request batch data before a cycle is started
- A buzzer can be activated at an incoming alarm (until alarm is acknowledged)
- A buzzer can be activated at cycle end (until a chamber door is opened)
- Display mode: The remain process time or the chamber temperature on side 2 will be displayed.
- The door opening / closing function can be set to an inching mode (door moves only as long as the door open or door close button is pressed)
 - 1o: open door on operating end
 - 1c: close door on operating end
 - 2o: open door on non operating end
 - 2c: close door on non operating end




Figure 8-32: Machine set-up display 4

The fifth display (refer to Figure 8-33) functions are (only with internal steam generator):

- Blow-down interval (default 1800s)
- Blow down pulse (default 2s)
- degassing time (default 2 sec); depending on water quality.
- to reduce peak current, the second heater can be delayed.
 - during heat-up phase (default =2 sec)
 - during standby phase (default=20 sec)




Figure 8-33: Machine set-up display 5

The sixth display (refer to Figure 8-34) functions are:

- Start point of the reduced steam pressure increase rate below sterilize pressure (default: 200mbar)
- Pressure increase rate before sterilization (reduced ramp, default: 200mbar/min) to prevent an over-temperature on the top of the chamber
- Open time steam valve gravitation: the steam inlet and, thus, also the steam outlet can be reduced by reducing the open time in the range between 10 and 100 %.
- Maximum chamber pressure during gravity: the chamber steam valve is closed when the maximum chamber pressure is exceeded.
- Print interval and data store during sterilize phase on built-in printer (default: 60 s)
- The software version of the panel is displayed




Figure 8-34: Machine set-up display 6

The seventh display (refer to Figure 8-35) functions are:

- Define the external documentation system (none / 8535 / 8535 BC/8565)
(BC: Barcode reader for recording the load)

- If an external documentation system type 8535 is connected via serial interface (RS485), enter the sterilizer system address communication port address, range 3..20)




Figure 8-35: Machine set-up ICS 8535

- IP address, documentation system 8535, with Ethernet connection: if the sterilizer is connected to an 8535 documentation system, each unit must have a TCP/IP address. Click on the "Set IP address" button to accept this address after it has been set.




Figure 8-36: Machine set-up ICS 8535 TCP/IP

8.10.5 Maintenance interval

The factory-set maintenance intervals can be changed by Belimed maintenance personnel when special circumstances require shorter intervals (e.g. quality of utilities do not fulfill our specifications).

The standard settings are defined in the 'Maintenance Manual' (Refer to Figure 8-37). The current values are displayed at the right. You can reset the maintenance counters at any time by pressing the "reset" buttons.




Figure 8-37: Maintenance interval display

8.10.6 Printer self-test

The printer system-test forces a printout (refer to Figure 8-38). Thereby the following are tested:

- Serial communication between PLC and Printer interface
- Serial communication between Printer interface and printer
- Printer interface
- Accuracy of Sensor inputs on printer interface (T2, P2), if a reference RTD (T2=134°C) and a reference input current (P2= 12mA) are connected to the sensor inputs.
- Printer function

A record is printed out with machine type and no., date, time, set values , actual values, allowed tolerance and result.

Additionally to the system- test, a printout of the alarm texts and the cycle phase texts can be forced.




Figure 8-38: Printer self test

8.10.7 Modify Cycle Parameters and copy program

The Belimed maintenance personnel can modify the following cycle parameters (refer to Figure 8-39):

- Cycle name
- Sterilize time
- Dry time

IMPORTANT:

The cycles listed in Table 6-1 have been validated using the techniques documented in AAMI ST-8 and ST-37. If different cycle parameters (sterilize time or dry time) are required, it is the responsibility of the health care facility to validate the cycle. Reference appropriate AAMI guidelines for validating sterilization cycles to assure the proper Sterility Assurance Level (SAL) as well as moisture retention acceptance criteria.

After modifying parameters press 'save parameters' button. A new version is generated on the basis of the date if the parameters are changed.




Figure 8-39: Modify cycle parameters display 1

All parameters can be displayed and modified in a further mask.




Figure 8-40: Modify cycle parameters display 2

Procedure for displaying or changing program parameters:

- Enter the parameter number in the "**parameter number**" field. Belimed Service or the technician will be able to provide you with information on the corresponding number.
- The current parameter is displayed without unit in the "**current value**" field
- In the case of modification: enter the new value in the "**new value**" field.
- In the case of modification: press the '**save parameters**' button. If the parameters are changed, a new version is generated on the basis of the date.
- The new value is displayed immediately in the "**current value**" field.

None of the parameter values have a unit. This means:

- Temperatures are displayed in tenths of 1 °C
- Pressure is displayed in mbar absolute
- Times are displayed in tenths of a minute.

In a third screen a complete program can be copied under a new number




Figure 8-41: Copy program display

Guidance to copy a program:

- Select the program to be copied. (See chapter.5.4.2 Select a cycle)
- Enter the number of the new program in field 'new program number'

Warning: Any existing program under that number will be overwritten !

- Press '**copy program**' button to copy the program.
- To modify the program name and parameters, the new program needs to be selected.

8.10.8 PLC functions

The PLC functions allow the Belimed maintenance personnel to read any data from the PLC.

Warning:

Write entries are allowed only to Belimed maintenance personnel, who are authorized by the sterilizer manufacturer. The functions, and addresses must be announced by the manufacturer before.

The Belimed maintenance personnel has access to the PLC data memory (refer to Figure 8-42).




Figure 8-42: PLC functions

Reading a value from PLC:

- Select the data type
- Enter the PLC memory address and DB- No.
- Push button 'read value'

Change a value in the PLC:

- Read the value, you want to change as described above
- enter the new value in the field 'value'
- Press button 'write value'

The new value is displayed on the left side.

8.10.9 I/O- Test

The Belimed maintenance personnel can read the digital PLC inputs for diagnostics and wiring tests.

He also can temporary set an digital output (refer to Figure 8-43).

Leaving the Mask resets the (manual set) outputs

Reading a digital input from PLC:

- Enter the Input byte address in the field 'input byte'

Values with 24V are displayed 'High'

Writing a digital output:

- Enter the output byte address in the field 'output byte'
- Select the desired bit in the lower row. The button changes to 'high'
- Push the 'write outputs' button. The corresponding output is set.

Warning:

Digital outputs are set without any conditions (not limit switches are checked). The outputs stay set until a reset command is set or the mask is closed. Only the emergency stop button stays in action.




Figure 8-43: I-O Test

9 Routine Monitoring

9.1 Biological monitoring

As part of the operator's verification of the sterilization process, biological indicators may be used to demonstrate that sterilization conditions have been met.

A live spore test utilizing *Geobazillus stearothermophilus* is the most reliable form of biological monitoring.

To verify the process insert the biological indicator in a test pack (refer to AAMI ST-8 guidelines) and place pack on the lower shelf. Run test pack through a typical Pre-vacuum cycle (Cycle no. 1 or 2).

On completion, forward test pack and monitor to appropriate personnel for evaluation.

9.2 Testing for pre-vacuum efficiency

9.2.1 Leak Test

The leak-rate-test measures the integrity of the sealed pressure vessel and associated piping to assure air is not being admitted to the sterilizer during vacuum steps.

After running a Leak Test, the leak rate is printed on the cycle documentation.

If the leak rate is out of range, additionally an alarm message is printed.

The malfunction must be reported to the supervisor, who will take appropriate action to determine the cause of the problem. Sterilizer should not be used during this time.

The leak rate value will help define a trend over a period of time if the integrity of the system begins to deteriorate.

9.2.2 Bowie-Dick Test

Daily run a Bowie-Dick Test cycle before processing any loads. This cycle should be used to test the adequacy of air removal from chamber and load, so that steam can penetrate the load. It is not a test for adequate exposure to heat in terms of time-at-temperature.

Tests such as Bowie-Dick are designed to document the removal of residual air from a sample challenge load.

In case of these tests, following exposure in a Bowie-Dick Test cycle, the pack is opened, the indicator examined and conclusions are drawn as to the pattern of residual air, if any, that remained in the pack during the sterilization cycle. Any indication of a malfunction must be reported to the supervisor, who will take appropriate action to determine the cause of the problem. Sterilizer should not be used during this time.

10 Routine Maintenance

CAUTION

Use only original BELIMED replacement parts for maintenance purposes: in all other cases the manufacturer can no more guaranty the specified performances of the equipment nor apply the warranty conditions.

10.1 Replacement of ribbon or paper roll

The printer is built in a strong plastic housing with removable cover. For replacing the ribbon or the paper roll, the cover must be removed from the housing.

Step 1: Remove the cover

Shift the cover while placing the two thumbs on the grooved rail and pressing into the opposite direction (Fig. 10-1).




Figure 10-1: Remove Cover

Step 2: replace the ribbon cassette

For replacing the printing ribbon press onto the edge of the ribbon cassette, where the words „PUSH“ and „EJECT“ appear (figure 11-2). After that the ribbon cassette is unlocked on the right hand side and it can be removed. Tighten the ribbon of the new cassette by turning the small wheel on the right hand side in the arrow-indicated direction (figure 11-3). Place the ribbon cassette over the paper. The paper must be between the textile ribbon and the plastic rail. Place the cassette in its definitive location until it engages clearly (figure 11-4)




Figure 10-2: Replace ribbon 1




Figure 10-3: Replace ribbon 2




Figure 10-4: Replace ribbon 3

Step 3: replace the paper roll

Paper rolls with a diameter of 38mm (1.5") are compatible with the printer housing.

Remove the ribbon cassette according to step 2. Remove the spindle with the core of the empty paper roll and place it into the core of the new roll. Place the roll in this way that the paper will be carried down and forwards

(figure 11-5). If necessary cut off the paper end straightly. Introduce the end of the paper from underneath into the corresponding slot at the printer until there is a resistance. In order to get the paper in front of the textile ribbon press the paper advance button until about 5 cm come out of the printer.

Introduce the ribbon cassette in its location as explained in step 2 (figure 11-6 to 11-8).




Figure 10-5: Replace paper 1




Figure 10-6: Replace paper 2




Figure 10-7: Replace paper 3




Figure 10-8: Replace paper 4

Step 4: Place the cover back in its position

Introduce the paper from the back side through the cover. Place the cover on the housing with a distance of about 1 cm from the edge and shift it as far as it goes (Fig. 11-9).




Figure 10-9: Place back printer cover

10.2 Cleaning of Chamber and Front panels

Clean chamber and front panels if dirty but at least once per week.

WARNING

Before entering the chamber (e.g. for cleaning) turn the key switch off (position '0') and keep the key in your pocket as long as you work in the chamber.

CAUTION

Clean the chamber in cold condition only.

Chamber


Clean with cotton towel or sponge. Use water or washing-up liquid. See to it that no residues are left.

Front panels

Clean with soft towel. Clean dry or use washing-up liquid or stainless steel cleaner EAB TOP (Belimed part no. 3-175).

Belimed Steam Sterilizer MST-V TOP 5000

TECHNICAL MANUAL



Model 5-5-9
535x535x965 mm

11 General Safety Precautions

It is essential that the following safety regulations be followed:

- This manual contains important information on proper maintenance of this sterilizer.
- Maintenance and servicing may be carried out only by BELIMED-trained staff authorized to do so. Follow the instructions in this the manual when carrying out maintenance and servicing work.
- The system may be placed into operation only in perfect working order and in compliance with the Operating Manual.
- Do not remove or bypass safety facilities.
- The pressure equipment must be inspected regularly pursuant to the applicable legal regulations.
- Please inform your superior in the event of faults or damage. Proceed on the basis of the troubleshooting chart for fault clearance.
- Since the system is operated with saturated or superheated steam, the chamber body, material being sterilized and loading equipment must be considered as hot. Use appropriate protective equipment.
- Take great care, particularly when opening the sterilizer door. The edges of the door and inside walls of the chamber are hot. The chamber interior and any batch trolleys used are also hot (>50°C).
- In the service compartment, all valves, pipes and chamber components must be considered hot until the operator has established otherwise.
- Keep the door's movement zone unobstructed when opening and closing the door.
- The owner must fit all pipes with shut-off valves.
- The owner must ensure that the media have the required operating pressures.
- Always switch off the master switch before carrying out maintenance or servicing work. Depressurize all tanks, containers and pipes before opening screw unions.
- Use only original spare parts. Make no modifications to the system. This would otherwise render warranty null and void.
- Only BELIMED's own staff or others under the supervision of BELIMED staff may move (i.e. load, unload or relocate) the system.
- The system must be disposed of professionally and in line with national regulations.
- Block the chamber door on the closed position by turning the door safety device 90 degrees, so that the door cannot open before you make any changes on the door drive (chain, clutch, motor).

11.1 Summary of warnings and cautions

The following is a summary of safety precautions to be observed when operating or servicing this unit.

WARNINGS indicate the potential for danger to personnel and CAUTIONS indicate the potential for damage to equipment. The precautions may be repeated where applicable throughout the manual. The following is a list of all safety precautions to be taken. Carefully read them before proceeding to use or service the unit.

WARNINGS

- **BURN AND SHOCK HAZARD:** Repairs and adjustments should be attempted only by authorized persons fully acquainted with this equipment. Use of inexperienced, unqualified persons to work on the equipment or the installation of unauthorized parts could cause personal injury or result in costly damage!
- Observe the label "Dangerous voltage": turn off main switch before opening
- **BURN HAZARD:** Allow sterilizer, generator (if applicable) and accessories to cool to room temperature before performing any cleaning or maintenance procedures.
- This sterilizer is not designed to process any liquids.
- **BURN HAZARD:** Sterilizer and rack/shelves will be HOT after cycle is run. Always wear protective gloves and apron (also face shield if processing liquids) when removing a processed load. Protective gloves and apron should also be worn when reloading sterilizer following previous operation.
- **FALL HAZARD:** To prevent falls keep floors dry by immediately wiping up any liquids or condensation in sterilizer loading or unloading area.
- **EXPLOSION HAZARD:** This sterilizer is not designed to process flammable liquids.
- **BURN HAZARD:** A steam supply malfunction may cause the sterilizer chamber to fill with scalding water. Do not open chamber door if the unit fails to complete an automatic cycle or if water leaks past the door gasket upon unlocking the door.
- **STERILITY ASSURANCE HAZARD:** According to EN285 standards, a measured leak rate greater 1.3 mbar/minute indicates a problem with the sterilizer. Refer the situation to a qualified maintenance technician before using the sterilizer further.
- Before entering the chamber (e.g. for cleaning), turn the key switch off (position '0') and keep the key in your pocket as long as you work in the chamber.
- Block the chamber door on the closed position by turning the door safety device 90 degrees, so that the door cannot open before you make any changes on the door drive (chain, clutch, motor).

CAUTIONS

- **CAUTION:** Avoid letting moisture get into chamber insulation, as it will cause rusting of the outer jacket.
- **CAUTION:** Never use sharp tools to push gasket into groove.
- **CAUTION:** Never use wire brush or steel wool on door and chamber assembly.
- **CAUTION:** Observe the Electrostatic Precautions outlined in the maintenance manual. Always wear a grounding wrist strap when removing or replacing PC boards or ICs.
- **CAUTION:** Solenoid valves are equipped with a special material which can be attacked by oils and grease. When replacing entire valve, wipe threads clean of cuttings oils and use Teflon tape to seal pipe joints.
- **CAUTION:** Handle siphon and bellows assembly gently to avoid damage.
- **CAUTION:** Keep the movement area of the door unobstructed during door opening and closing.
- **CAUTION:** Any alteration of the sterilizer which affects its operation will void the warranty and could violate state and local regulations and jeopardize insurance coverage.
- **CAUTION:** Use only original Belimed replacement parts for maintenance purposes: in all other cases the manufacturer can no more guaranty the specified performances of the equipment nor apply the warranty conditions.
- **CAUTION:** In order to ensure the proper function of the sterilizer, run a leak test cycle daily.
- **CAUTION:** The Bowie-Dick Test Cycle has to be performed daily.
- **CAUTION:** Sterilize only goods which are declared to be steam-sterilizable by their manufacturer and compatible with the chosen cycle parameters (maximum temperature)
- **CAUTION:** Respect the maximum number of sterilization's given by the manufacturer of the goods to be sterilized. Check their proper function according to the manufacturer's information of the goods (e.g. instruments)
- **CAUTION: STERILITY ASSURANCE HAZARD:** Load sterility may be compromised if the chemical or biological indicator or the Bowie-Dick Test indicate a potential problem. If these indicators show a potential problem, refer the situation to a qualified maintenance technician before using the sterilizer further.

12 Maintenance / periodic tests

12.1 Safety information

You must always do the following before starting maintenance:

- Switch off the key-operated switch on the operating panel and remove the key.
- Switch off the master switch on the switch cabinet and lock it with a padlock.
- Shut off all energy supplies (water, steam and compressed air) at the shut-off valve.

Before entering the chamber (e.g. for cleaning), turn the key switch off (position '0') and keep the key in your pocket as long as you work in the chamber. Enter the chamber from load side (operating end) only.

12.2 Maintenance intervals

The maintenance intervals are defined in the TOP5000 maintenance menu. 3 maintenance levels can be set. A maintenance interval can be defined on the basis of the operating hours. If a number of operating hours (operating hours of the PLC) are elapsed, the maintenance prompt is displayed on the screen when the first event has occurred. If a maintenance level is reset, the operating hours of all lower levels are also reset. The higher levels must always be defined as a multiple of the next layer down.

Maintenance messages not displayed on the TOP 5000 screen:

- Daily
- Weekly 50h (10 hours/day, 5 days a week)

Maintenance messages displayed on the TOP 5000 screen:

- Level 1: Quarterly 650h
- Level 2: Yearly 2600h
- Level 3: 2-Years 5200h

12.3 Maintenance personnel

The test person must be authorized and trained by the sterilizer manufacturer.

12.4 Maintenance plan

12.4.1 Daily Maintenance (before start sterilizing)

No maintenance message is displayed for daily tests.

Item to be serviced	Work to be performed	Performed by
Cleaning	Clean: <ul style="list-style-type: none">• Sterilizer surface, if dirty	user
Printer	Printer paper: If red paper marking is visible on printout, replace the paper roll. If the printout is poor, replace printer ribbon	user
Leak Test	Run a Warm- Up & Leak Test cycle. Any indication of a malfunction, (leak rate > 1.3mbar/Min) must be reported to the supervisor, who will take appropriate action to determine the cause of the problem. Sterilizer should not be used during this time.	user
Steam penetration test	Run a Bowie- Dick Test. If the sterilizer was turned off before (cold chamber), run a Warm-up cycle first. Any indication of a malfunction must be reported to the supervisor, who will take appropriate action to determine the cause of the problem. Sterilizer should not be used during this time.	user

12.4.2 Weekly Maintenance

No maintenance message is displayed for weekly tests.

Item to be serviced	Work to be performed	Performed by
Cleaning	Clean: Chamber Chamber drain strainer Sterilizer surface	user
Door seals	Examine the door seal: <ul style="list-style-type: none"> • Door seal must be completely sucked back in the frame • Door seal must not be mechanically damaged. Any indication of a malfunction must be reported to the supervisor, who will take appropriate action to determine the cause of the problem. Sterilizer should not be used during this time. Door seal must be completely sucked back in the frame. If not, clean dirty door seal and grease.	user test person
Option motorized door: Door closing force	While closing the door you must be able to stop the movement by pressing down the door Alternatively: Place a test weight of 15 kg (at Belimed available) on the open door. The door must not close, when the button 'close door' is pressed.	user
Piping system	Open service doors and check for leaks In case of a leak: <ul style="list-style-type: none"> • Reseal leaking points and replace seals or exchange parts if necessary 	user test person

12.4.3 Maintenance level 1

Maintenance interval 650 hours (Quarterly)

Item to be serviced	Work to be performed	Performed by
Cleaning	Clean: <ul style="list-style-type: none"> • Chamber • Chamber drain strainer • Sterilizer surface 	user
Door seals	<ul style="list-style-type: none"> • Examine the door seal. • Clean door seal and grease. • Replace door seal, if mechanically damaged. 	test person
Option motorized door: Door closing force	<ul style="list-style-type: none"> • Place a test weight of 15 kg (at Belimed available) on the open door. • The door must not close, when the button 'close door' is pressed. 	test person
EMERGENCY-STOP	<ul style="list-style-type: none"> • Check for proper operation. • 2- door sterilizers: check function on operating and non-operating end 	test person
Piping system	<ul style="list-style-type: none"> • Check for leaks. • Reseal leaking points and replace seals or exchange parts if necessary 	test person
Compressed air filter	<ul style="list-style-type: none"> • Drain off the condensate in the filter casing using the drain plug 	test person
Sterile filter	<ul style="list-style-type: none"> • Unscrew sterile filter and sterilize filter with 134°C Fabric Packs cycle. 	test person
Control panel fan	<ul style="list-style-type: none"> • Check the air inlet fan for proper operation 	test person
Leak Test	<ul style="list-style-type: none"> • Run a Warm- Up & Leak Test. <p>Any indication of a malfunction, (leak rate > 1.3mbar/Min) must be reported to the supervisor, who will take appropriate action to determine the cause of the problem. Sterilizer should not be used during this time.</p>	user
Steam penetration test	<ul style="list-style-type: none"> • Run a Bowie- Dick Test. <p>If the sterilizer was turned off before (cold chamber), run a Warm-up cycle first. Any indication of a malfunction must be reported to the supervisor, who will take appropriate action to determine the cause of the problem. Sterilizer should not be used during this time.</p>	user

12.4.4 Maintenance level 2

Maintenance interval 2600 hours (Yearly)

Item to be serviced	Work to be performed	Performed by
Cleaning	Clean: <ul style="list-style-type: none"> • Chamber • Chamber drain strainer • Sterilizer surface 	user
Door seal	<ul style="list-style-type: none"> • Replace door seal 	test person
Door drive	<ul style="list-style-type: none"> • Grease chain, as required 	test person
Option motorized door: Door drive	<ul style="list-style-type: none"> • As required, clean chain and grease • Check chain tension, as required 	test person
EMERGENCY-STOP	<ul style="list-style-type: none"> • Check for proper operation • 2- door sterilizer: check on operating end and non operating end 	test person
Condensate steam piping jacket	<ul style="list-style-type: none"> • Replace steam trap capsule elements 	test person
Piping system	<ul style="list-style-type: none"> • Check for leaks: • Reseal leaking points, replace seals or exchange parts if necessary 	test person
Strainer in the pipes, dirt trap	<ul style="list-style-type: none"> • Clean strainer inserts: • Remove and clean the strainer on the marked pipe screw unions or on the dirt trap 	test person
Orifices, nozzles	<ul style="list-style-type: none"> • Clean orifice plates and nozzles • Remove and clean orifices on the pipe 	test person
Sterile filter	<ul style="list-style-type: none"> • Replace sterile filter 	test person
Compressed air filter	<ul style="list-style-type: none"> • Drain filter • Drain off the condensate in the filter casing using the drain plug 	test person
Control panel fan	<ul style="list-style-type: none"> • Check the air inlet fan for proper operation • Exchange the fan inlet- and outlet filters 	test person
Solenoid valves	<ul style="list-style-type: none"> • Check for proper operation 	test person
Angle-seat valves	<ul style="list-style-type: none"> • Check for leak • Replace piston actuator or valve in case of a leak or malfunction 	test person
Vacuum pump	<ul style="list-style-type: none"> • Check suction capacity with dry and empty chamber: Evacuation time from atmosphere to 70mbar must not exceed: Model 6-6-6: 3 min.; model 6-6-9: 4 min.; model 6-6-12: 5min. 	test person
Gauges	<ul style="list-style-type: none"> • Check for proper operation 	test person
One-way restrictor on V623 and V624	<ul style="list-style-type: none"> • Check adjustment: opening time: ~1-2 seconds 	test person
pressure vessel steam generator (option)	<ul style="list-style-type: none"> • Take and analyze Steam condensate (external laboratory). • Examine if necessary feed water quality and clean pressure vessel 	test person

12.4.5 Maintenance level 3

Maintenance interval 5200 hours (2-year interval)

Item to be serviced	Work to be performed	Performed by
All level 2 items	All level 2 work as described above	
PLC Battery	<p>The batteries must be replaced on the PLC</p> <ul style="list-style-type: none"> use Lithium battery, Type CR2032 <p>Warning: Do not switch off electrical power while changing the battery</p>	test person

12.4.6 Yearly safety check

The Safety check is to log in a maintenance record.

No.	Item to be tested	Work to be performed	refer to chapter	Performed by
1	Door seal	<ul style="list-style-type: none"> Examine the door seal 	16.2	test person
2	Atmospheric pressure switch	<ul style="list-style-type: none"> Check switch point 	17.3.1	test person
3	Door lock pressure switch	<ul style="list-style-type: none"> Check switch point g 	17.3.2	test person
4	Door Limit switches	<ul style="list-style-type: none"> Check switch points Check for firm attachment (they may have come loose) 	17.3.4	test person
5	Temperature sensors in Chamber drain	<ul style="list-style-type: none"> Calibration verification Calibrate or replace if required 	14.1, 15.1	test person
6	Verification of calibration Pressure sensors	<ul style="list-style-type: none"> Calibration verification Calibrate or replace pressure sensors, if required 	14.2, 14.3, 15.2	test person
7	Pressure regulator door seal	<ul style="list-style-type: none"> Check setting 	18.1	test person
8	Safety valve chamber	<ul style="list-style-type: none"> Check safety valve 	17.2	test person
9	Safety valve jacket	<ul style="list-style-type: none"> Check safety valve 	17.2	test person
11	EMERGENCY-STOP	<ul style="list-style-type: none"> Check for proper operation 	17.1	test person
12	Option motorized door: Door closing force	<ul style="list-style-type: none"> Place a test weight of 15 kg (at Belimed available) on the open door. The door must not close, when the button 'close door' is pressed. 	16.3	test person
16	Chamber steam valve	<ul style="list-style-type: none"> Function test 	17.3.3	test person
17	Compressed air pressure switch	<ul style="list-style-type: none"> Function test (pressure set-point) 	18.1	test person

12.4.7 Periodic pressure vessel check

The pressure vessel must be checked according to the national or local regulations.

13 Cleaning

Warning

Before entering the chamber (e.g. for cleaning), turn the key switch on the operating panel off (position '0') and keep the key with you as long as you are in the chamber.

Caution

Clean the chamber in cold condition only.

Caution

Do not use any cleaning agent containing salt or chloride or any other cleaning agent which is not suitable for stainless steel.

Chamber

Clean with cotton towel or sponge. Use water or washing-up liquid only. Make sure that no residues remain on the surface.

Front panels

Clean with soft towel. Clean dry or use washing-up liquid or suitable stainless steel cleaner.

We recommend EAB-TOP, Belimed part no.: 350 07473.

Touch panel

Before cleaning the touch screen make sure the sterilizer is turned off with the key switch. Clean the touch screen with a damp towel. Apply only water or washing-up liquid or a suitable screen cleaning foam to moisturize cleaning towel. Do not spray any cleaning agents directly on the screen. Never use aggressive or abrasive detergents.

14 Calibration verification

Calibration involves recording and documenting the difference between indicated value and correct value of the process-related parameters on the measuring facility .

The temperature and pressure measuring chain must be calibrated at least once per year.

The user is responsible for calibration. The user may amend the adjustment points suggested below by own instructions. It is important that calibration be performed at the sterilization temperature(s) used.

For calibration use official calibrated and certified temperature and pressure references. Calibration certificate must not be expired.

14.1 Temperature calibration verification

P1 chamber	T1 chamber	status
886 mbar	85.3 °C	standby

menu		
log off	info	readings
set-up	print last cycle	alarms
cycle intervention	batch data	user administration
←		maintenance

Figure 14-1: Menu

P1 chamber	T1 chamber	status
959 mbar	25.0 °C	standby

readings		
P1 chamber	P2 cham. doc.	P3 jacket
959 mbar	959 mbar	2127 mbar
P4 steam	T1 chamber	T2 cham. doc.
4000 mbar	25.0 °C	25.0 °C
T4 condenser	T5 circ. tank	
25.0 °C	25.0 °C	

←

Figure 14-2: Readings

Calibrate offset

- Select on the operating panel the menu 'readings' (refer to Figure 14-1 and Figure 14-2).
- Insert all sensors in a temperature calibration unit (or bucket). Sensors must be immersed at least 95% (minimum 3 ½") into the temperature reference.
- Set the reference temperature to approx. 20°C (room temperature).
- Wait until the temperature reading is stable for at least 5 minutes.
- If the temperature T1 or T2 shown on the display differs by more than ± 0.2 degrees from the actual value (certified temperature reference), the offset point is to be calibrated.
- If the temperature T4,T5 shown on the display differs by more than ± 2 degrees from the actual value (certified temperature reference), the sensor is to be replaced.

Calibrate gain

- Set the reference temperature to approx. 134°C.
- Insert all sensors in a temperature calibration unit (e.g. oil bath). Sensors must be immersed at least 95% (minimum 3 ½") into the temperature reference.
- Wait until the temperature reading is stable for at least 5 minutes.
- If the temperature T1 or T2 indicated on the display differs by more than ± 0.2 degrees from the actual value, the gain is to be calibrated.
- If the temperature T4,T5 shown on the display differs by more than ± 2 degrees from the actual value (certified temperature reference), the sensor is to be replaced.

14.2 Chamber pressure transducer calibration verification

- Connect a calibrated absolute-value vacuum-gauge to the chamber.
- Select the Leak test program and start it. Wait until the program enters stabilize phase.
- Compare the pressure readings on the display with the pressure reading of the pressure gauge.
- If the chamber vacuum of P1 or P2 differs by more than ± 15 mbar from the effective value, the pressure transducer must be calibrated.
- Abort the cycle

14.3 Jacket and steam pressure transducer calibration verification

Warning

Before disconnecting the pressure transducers ensure that the steam supply valve is turned off and that there is no overpressure in the jacket and the steam supply pipe.

Offset calibration

- Disconnect jacket pressure transducer P3 and steam pressure transducer P4.
- Use a calibrated pressure transducer (Range 0..4000mbara, accuracy $\pm 0.15\%$ or better) as a pressure reference.
- Compare the pressure readings of P3 and P4 on the display with the pressure reference (ambient pressure).
- If the pressure P3 or P4 differs by more than ± 100 mbar from the reference value, the faulty pressure transducer should be replaced.

Remark

A vacuum calibration is not required because there is always an overpressure in steam supply and jacket during normal operation.

Gain calibration:

- Install the transducers P3 and P4 and the pressure reference to a pressure test tank.
- Set the pressure in the test box (at least 1 gallon) to 3000mbar absolute by compressed air. Close the compressed air inlet valve.
- Wait one minute for pressure stabilization.
- If the pressure P3 or P4 differs by more than ± 100 mbar from the reference value, the faulty pressure transducer should be replaced.
- Release pressure in test tank to ambient pressure.
- Reinstall the pressure transducers P3 and P4.

15 Sensors adjustment

Calibration involves setting or adjusting the measuring facility with the aim of maintaining the difference between indicated and correct value of the process-related parameters as small as possible (DIN 58950-1).

The chamber temperature or pressure measuring chain has to be re-calibrated if deviations of more than $\pm 0.2^{\circ}\text{C}$ in the case of temperature or more than ± 20 mbar in the case of pressure are determined during calibration.

15.1 Temperature sensors adjustment

The temperature calibration must be done in degree Celsius ($^{\circ}\text{C}$) (refer to Figure)

15.1.1 Offset calibration:

- Put the temperature sensor T1, T2 into a temperature reference (accuracy $\pm 0.1^{\circ}\text{C}$).
- Set the reference temperature to 0°C or to ambient temperature (20 - 25°C).
- Wait until the temperature T1 and T2 are stable (at least 5 minutes).
- Enter the actual value in 'reference temperature' field.
- Choose the corresponding temperature sensor T1, T2 or, optionally, T3 (the selected temperature sensor lights green).
- The uncorrected value of the selected temperature sensor is displayed
- Press 'offset' button
- Press 'confirm values' button

15.1.2 Gain calibration:

- The gain calibration must be done after the offset calibration.
- Set the reference temperature to 134°C .
- Wait until the temperature T1 and T2 are stable.
- Enter the actual (reference) value in 'reference temperature' field.
- Choose the corresponding temperature sensor T1, T2 or, optionally, T3 (the selected temperature sensor lights green).
- The uncorrected value of the selected temperature sensor is displayed
- Press 'gain' button
- Press 'confirm values' button




Figure15-1: Calibrate temperature sensors display

15.2 Pressure sensors adjustment

The pressure calibration must be done in mbar absolute (mbar) (refer to Figure)

15.2.1 Offset calibration:

- Install a calibrated vacuum transducer (range 0..1000mbara, accuracy 0.2%) on the chamber validation flange..
- Start a leak test cycle
- Wait until the leak test starts .
- Enter the actual pressure value (from reference transducer) in '**reference pressure**' field. The pressure should be between 50mbar and 80mbar.
- Choose the corresponding pressure sensor P1 or P2 (the selected pressure sensor lights green). The uncorrected value of the selected pressure sensor is displayed
- Press '**offset**' button
- Press '**confirm values**' button
- Abort the cycle

15.2.2 Gain calibration:

- Install a calibrated pressure transducer (Range 0..4000mbara, accuracy 0.2%) and a shut-off valve V1 on the chamber validation flange. Install compressed air (> 30psi) to the other side of the shut-off valve V1)
- Start a leak test cycle
- Disconnect the connector X43 on the print 2005-07 (all valves are closed, door seals are locked)
- Open shut-off valve V1. When the chamber pressure reaches 3200mbar, close V1.
- Wait until the pressure is stable
- Enter the actual pressure value (from reference transducer) in '**reference pressure**' field.
- Choose the corresponding pressure sensor P1 or P2 (the selected pressure sensor lights green). The uncorrected value of the selected pressure sensor is displayed
- Press '**gain**' button
- Press '**confirm value**' button
- Connect the connector X43
- Abort the cycle




Figure 15-2: Calibrate pressure sensors display

16 Servicing instructions

16.1 Opening control panel

The pullout control panel is placed below the chamber. Open the lower maintenance door. Then push the handle on the left side of the panel (Figure 16-1) and pull out the control panel. Lock the panel in the open position with the 2 handles on the lower left and lower right side of the panel (Figure 16-2).




Figure 16-1 : unlock control panel




Figure 16-2: Lock panel in open position

16.2 Cleaning and replacement of door seal

16.2.1 Door seal on operating end

To replace the door seal, press button 'replace door seal' in the 'maintenance' menu:




Figure 16-3 : Menu functions




Figure 16-4: Maintenance functions




Figure 16-5: replace door seal

- Fully open the door, the sterilizer must be in standby position.
- The door seal can be removed with the 'replace door seal' function in the operating panel maintenance menu (refer to Figure 16-5):
- Press 'menu' button.
- Press 'maintenance' button.
- Press 'replace door seal' button.
- Press 'blow out door seal' button.
- The door seal is then pushed out a little.
- Remove the seal.
- Clean the seal and grease it with grease Belimed part no. 350 07118.
- Insert a seal as follows (refer to Figure16-6)
- Insert the seal in all 4 corners
- Insert the seal in the center at the top and bottom
- Then insert it half way up at the left and right side
- Then insert the seal between two sections
- Then press the rest of the seal into the groove. (It is always best to do this at the center, between two sections already in the groove).
- When the seal has been fitted completely, press '**pull back door seal**' button.
- The seal is sucked into the groove by vacuum for approx. 15 seconds.
- If the seal is not yet neatly fitted in the groove, you can press '**pull back door seal**' button again.
- Close the door




Figure16-6: replace door seal

16.2.2 Door seal on non operating end

- Fully open the door on the non operating end (the sterilizer must be in 'cycle complete' position)
- If the cycle is in the standby position, you can set the sterilizer to 'complete' position by pressing '**complete**' button.

- The door seal can be removed by pressing the 'open door' and the 'close door' buttons together on the non operating end.
- The door seal is then pushed out a little.
- Remove the seal.
- Clean the seal and grease it with grease Belimed part no. 350 07118
- Insert a seal as described in chapter 16.2.1.
- When the seal has been fitted completely, press 'close door' button.
- The seal is sucked into the groove by vacuum for approx. 15 seconds.
- Close the door, press button 'close door'.
- Press 'complete' button to go back to standby.

16.3 Adjust door sliding clutch (option motorized doors)

The clutch must be adjusted, so that the door closing-force is between 75N and 100N.

The force must not exceed 150N.




Figure 16-7: sliding clutch on door motor

For door sliding clutch adjustment proceed as follows:

16.3.1 Check the closing force

- With the door open, place a total test weight of **15 Kgs** (2 x 7.5 Kgs, available at Belimed) on the top of the door frame. Press the button „close door“. The door must not close.
- Place a test weight of **7.5 Kgs** over the door. Press the button 'close door'. The door must close.

If one of this two conditions is not fulfilled, you have to adjust the sliding clutch.

16.3.2 Door sliding clutch adjustment

16.3.2.1 Preparation

- Close the door.
- Switch off the sterilizer and pull out the key of the switch. Carry the switch with you in your pocket during the whole maintenance work.
- Open the maintenance service door below the chamber.
- Block the chamber door at the closed position by turning the door safety device 90 degrees, so that the door cannot open when the sliding clutch is released.
- Open the Hexagon screw 1 (Figure 16-7).

16.3.2.2 Increase of the closing pressure:

Rotate the tension screw 2 (Figure 16-7) on the sliding clutch by approx. 1/8 turn clockwise.

Tighten the Hexagon screw 1 (Figure 16-7).

16.3.2.3 Decrease of the closing pressure:

Rotate the tension screw 2 (Figure 16-7) on the by approx. 1/8 turn counterclockwise.

Tighten the Hexagon screw 1 (Figure 16-7) .

16.3.2.4 Door release

- Turn back the door safety device parallel to the chamber door.
- Close the maintenance service door below the chamber.
- Switch the sterilizer on again.
- Open the chamber door.
- Check the closing pressure according chapter 16.3.1

16.4 Replacing the batteries (accumulators)

- The PLC controller must be switched on to prevent data loss.
- Remove MPI- plug (to operating panel)
- Remove the controller cover (PLC).
- You can now exchange the batteries.
Use only the lithium battery 3V. type CR2032.
- Clip the controller cover back on.
- Put the MPI- plug in and fasten with 2 connector- screws

16.5 Arrangements to be taken before and after an operation interruption

Arrangements before and after a longer interruption of the sterilizer.

All utilities shut-off valves have to be closed.

The electric main switch must be switched off.

Prior to a unit interruption of more than 1 month, the vacuum pump has to be preserved with a water soluble anti-corrosion preservation agent. Depending on the size of the vacuum pump, a quantity of 1-2 lt. preservation fluid has to be filled in to the pump.

- After an interruption of more than 1 month, the rotation capability of the vacuum pump must be checked before starting the operation. (In case of the pump sticking, the grip shall be loosen by rotating the pump shaft).

In case of frost (ice) danger :

- the pressure sensors and the pressure switches are to be removed.
- all the pipes are to be drained.
- all components (pumps, heat exchangers, valves) are to be drained.

17 Tests and inspections

Safety-related monitoring functions are also implemented independently of the PLC in order to enhance operating and system safety.

17.1 Emergency-Stop (double door sterilizers only)

The Emergency-Stop switch A925 on non operating end disconnects the supply voltage of the digital outputs.

The following components are disconnected from the electrical power supply:

- Door motors (Option motorized doors)
- Vacuum pump
- All angle seat valves

The following components remain connected to the electrical power supply:

- PLC control
- Operating panels
- Solenoid valves for cooling down the condensate (to prevent hot water can get to drain)
- Pressure and temperature sensors
- All limit switches, pressure switches and pushbuttons
- On the operating end the main disconnect switch disconnects the whole power (208V/60Hz and 24VDC).

17.2 Safety valves

Caution!

Safety valves may be inspected only by trained servicing personnel. Such inspections involve manual interventions which require a precise knowledge of the system!

17.2.1 General:

Safety valves are set and lead-sealed by approved and authorized agencies resp. companies. The official acceptance testing authority conducts an initial and periodic inspection of the safety valves. An interim inspection should be conducted only when absolutely necessary. After an inspection of the safety valve, there is a risk that the valve seat may leak as the result of fouling, deformation or damage, thus no longer guaranteeing perfect sterilization. Consequently, after each inspection, check the safety valve for leaks by removing the blow-off line and inspecting the safety valve at the outlet connector. Leaking safety valves must be replaced! The permitted operating gauge pressures can be read off on the rating plate.

17.2.2 Inspection in fitted condition

17.2.2.1 Safety valve, jacket Y 702

BURN HAZARD

Waste steam streams out through the safety valve. Safety valve drain must be connected to a pipework leading outside maintenance area. Use safety goggles.

- Select the warm-up cycle so that the jacket is heated.
- Monitor the jacket pressure on the pressure gauge PI 804. Under no circumstances may the pressure in the jacket increase by more than 15% above the permitted operating gauge pressure.

- The safety valve must open at the specified response pressure +/- 10%. Activate the jacket steam valve manually via the related pilot valve V622.1 (refer to Figure 20-3) until the safety valve blows off or the pressure is max. 15% above the permitted operating gauge pressure.
- The steam pressure in the steam supply line may need to be increased.
- If this is not possible, inspect the safety valve after removing it.
- If the safety valve does not blow off within the specified limits, it must be replaced or readjusted and lead-sealed by an approved agency.

17.2.2.2 Safety valve, chamber Y 700

- Start the heating program and wait until the program switches to the sterilization phase.
- Monitor the chamber pressure on the pressure gauge PI 800. Under no circumstances may the pressure in the chamber increase by more than 15% above the permitted operating gauge pressure.
- The safety valve must open at the specified response pressure +/- 10%.
- Activate the chamber steam valve manually via the related pilot valve V623.1 (refer to Figure 20-3) until the safety valve blows off or the pressure is max. 15 % above the permitted operating gauge pressure.
- The steam pressure in the steam supply line may need to be increased.
- If this is not possible, the safety valve must be inspected after removing it.
- If the safety valve does not blow off within the specified limits, it must be replaced or readjusted and lead-sealed by an approved agency.

17.2.3 Inspection in removed condition

If the required test pressure is not available, the safety valve must be removed and inspected by an external test agency. Ensure that the vessels (chamber, jacket) are depressurized before removing the safety valves.

17.3 Safety equipment, door lock

17.3.1 Chamber atmospheric pressure, safety switch

17.3.1.1 Checking the switch point, chamber atmospheric pressure

- Start a Warm-up & Leak Test cycle (P4).
- Wait until the program is in the heat-up phase and the chamber pressure is at least +0.5 bar g (1500mbar a).
- LED D6 on the I/O- print (in control panel) must not lit, as long as the chamber pressure is above 0.12 bar g.
- Abort the cycle.
- The LED D6 must lit green, as soon as the chamber pressure drops below 0.1 bar g (refer to Figure).




Figure 17-1: Ambient pressure display

17.3.1.2 Atmospheric pressure switch function test

- Start a Warm-up & Leak Test cycle.
- Wait until the program is in the steam phase and the chamber pressure is at least +0.5 bar g (1500mbara).
- The LED D6 on the I/O- print must not lit (refer to Figure).
- Disconnect the plug (pin 1) from pressure switch S335.
- The alarm message 'Atmospheric pressure monitoring' must be displayed on operating panel.
- Abort the cycle.

17.3.2 Door seal pressure

- Start a Leak Test cycle.
- Safety pilot valve Y334 must open and release compressed air to pilot valve Y623 (V623.1).

Remark

The position indicator juts out approx. 2mm, if a pilot valve is actuated (refer to Figure 20-3).

- Reduce door seal air pressure on reducing valve V300 to 2 bar g (refer to Figure 20-3).
- The door seal pressure slowly drops.
- When the door seal pressure drops below 2.2 bar (pressure switch S331), an alarm message 'Door not locked' must be displayed on the operating panel and the cycle aborts.
- The pilot valve Y334 must close and interrupt compressed air to steam to chamber pilot valve Y623 (V623.1).

Remark

The position indicator juts out approx. 2mm, if a pilot valve is actuated

- Acknowledge the alarms.
- Start a Leak Test cycle again and wait 10 seconds.
- Increase door seal pressure on reducing valve V300 to 3 bar g.
- Acknowledge the alarms and press 'open door' button.




Figure 17-2: compressed air pressure reducing valve V300




Figure 17-3 pilot valve battery

17.3.3 Door limit switches and steam supply to chamber monitoring

Steam to chamber must be locked until the chamber doors are closed and locked.

- Start a Leak Test cycle
- Safety pilot valve Y334 must open and release compressed air to pilot valve Y623 (V623.1).

Remark

The position indicator juts out approx. 2mm, if a pilot valve is actuated (refer to Figure 20-3).

- Actuate door limit switch S332 on operating end (on the right side below chamber door)

- An alarm message 'Door 1 not closed 'must be displayed on the operating panel and the cycle aborts.
- The pilot valve Y334 must close and interrupt compressed air to steam to chamber pilot valve Y623 (V623.1).

17.3.3.1 Two door sterilizers:

- Start a Leak Test cycle
- Actuate door limit switch S333 on non operating (on the left side below chamber door)
- An alarm message 'Door 1 not closed 'must be displayed on the operating panel and the cycle aborts.
- The pilot valve Y334 must close and interrupt compressed air to steam to chamber pilot valve Y623 (V623.1).
- Open the chamber door and ensure, that no steam enters into the chamber (V623 must not be leaky).

17.3.4 Door limit switch adjustment

Test feature	Procedure/Measuring equipment	Set-point
Limit switch S332 „Door 1 closed“ on operating end	<ul style="list-style-type: none"> • Open door 10 cm on operating end. • Close door manually. Limit switch S332 shall switch off (metering point: terminal X33.12), when distance A between door and limit stop is between 1-3mm. • Slowly open door. Limit switch S332 shall switch on (metering point: terminal X33.12), when distance A between door and limit stop is ≤ 10mm. 	
Limit switch S342 „Door 1 open“ on operating end (Option motorized door)	<ul style="list-style-type: none"> • Open door. • Limit switch S342 shall switch on (metering point: terminal X33.16 =24V), when distance B between lower chamber border and upper door border is 20-30 mm. 	
Limit switch S333 „Door 2 closed“ on non operating end (option 2-door sterilizer)	<ul style="list-style-type: none"> • Manually open door 10 cm. • Close door manually. Limit switch S333 shall switch off (metering point: terminal X34.5), when distance A between door and limit stop is between 1-3mm. • Slowly open door. Limit switch S333 shall switch on (metering point: terminal X34.5), when the distance A between door and limit stop is ≤ 10mm. 	
Limit switch S343 „Door 2 open“ on non operating end (option 2-door sterilizer with motorized doors)	<ul style="list-style-type: none"> • Open door. • Limit switch S343 shall switch on (metering point: terminal X33.9 =24V), when distance B between lower chamber border and upper door border is between 20-30 mm. 	

18 Settings

18.1 Pressure switches, pressure regulators and safety valves

The settings must be checked in accordance with the table below. In the case of deviations, the corresponding elements must be readjusted, adjusted or replaced.

Designation:	Description:	Setting:
S 331	Pressure switch, door seal	with increasing pressure 2.6 ± 0.2 bar gauge pressure with decreasing pressure 2.2 ± 0.2 bar gauge
S 335	Pressure switch chamber atmospheric pressure	with increasing pressure 0.12 ± 0.01 bar gauge with decreasing pressure 0.10 ± 0.01 bar gauge
S 346	Pressure switch Compressed air	with increasing pressure 4.5 ± 0.3 bar gauge with decreasing pressure 4.0 ± 0.3 bar gauge
V 300	Pressure regulator, compressed air door seal	with increasing pressure on PI301 3.0 ± 0.15 bar gauge
Y 700	Safety valve, chamber	40 ± 0.4 psi gauge
Y 702	Safety valve, jacket	40 ± 0.4 psi gauge

19 Tests after electrical repair

An insulation test and a PE wire test must be conducted after exchanging low-voltage components (motors, contactors and heaters).

The following minimum requirements must be complied with:

- Insulation resistance: >1Mohm
- Continuity of the PE wire system connection:

• Continuity of the PE wire system connection: Minimum cross-section of the PE wire [mm²]	Max. measured voltage drop [V] (at 10 A test current)
1.0	3.3
1.5	2.6
2.5	1.9
4	1.4
≥6	1.0

Figure 6-4: Cycle Graph - Flash cycles	48
Figure 6-5: Cycle Graph - DART / Bowie-Dick Test cycle	49
Figure 6-6: Cycle Graph - Warm-up & Leak Test cycle	50
Figure 7-1: Double-wrapping: Square fold	52
Figure 7-2: Double-wrapping: Envelope fold	54
Figure 7-3: Loading shelves with liquid load	56
Figure 8-1: Password entry	62
Figure 8-2: Menu display	63
Figure 8-3: User administration	63
Figure 8-4: Menu display	64
Figure 8-5: Info display	64
Figure 8-6: Menu display	65
Figure 8-7: Readings display	65
Figure 8-8: Menu display	66
Figure 8-9: Set-up display	66
Figure 8-10: Early start / night-time shutdown display	67
Figure 8-11: Set clock display	68
Figure 8-12: Language and units display	68
Figure 8-13: Configure cycle selection	69
Figure 8-14: Set-up display	69
Figure 8-15: Menu print last cycle display	70
Figure 8-16: Out-of-cycle display with active alarm	70
Figure 8-17: In-cycle display with active alarm	70
Figure 8-18: Alarm display	71
Figure 8-19: Menu display	74
Figure 8-20: Cycle intervention display	74
Figure 8-21: Manual batch data entry	75
Figure 8-22: Scanner based batch data entry	76
Figure 8-23: Automatic program selection by ICS 8565	76
Figure 8-24: Out-of-cycle display with maintenance message	77
Figure 8-25: Maintenance message display	77
Figure 8-26: Menu display	77
Figure 8-27: Maintenance menu display	78
Figure 8-28: Replace door seal display	79
Figure 8-29: Machine set-up display 1	80
Figure 8-30: Machine set-up display 2	80
Figure 8-31: Machine set-up display 3	81
Figure 8-32: Machine set-up display 4	81
Figure 8-33: Machine set-up display 5	82
Figure 8-34: Machine set-up display 6	82
Figure 8-35: Machine set-up ICS 8535	83


Figure 8-36: Machine set-up ICS 8535 TCP/IP.....	83
Figure 8-37: Maintenance interval display	84
Figure 8-38: Printer self test.....	84
Figure 8-39: Modify cycle parameters display 1.....	85
Figure 8-40: Modify cycle parameters display 2.....	85
Figure 8-41: Copy program display.....	86
Figure 8-42: PLC functions	87
Figure 8-43: I-O Test.....	88
Figure 10-1: Remove Cover.....	90
Figure 10-2: Replace ribbon 1	90
Figure 10-3: Replace ribbon 2	90
Figure 10-4: Replace ribbon 3	90
Figure 10-5: Replace paper 1	91
Figure 10-6: Replace paper 2	91
Figure 10-7: Replace paper 3	91
Figure 10-8: Replace paper 4	91
Figure 10-9: Place back printer cover	91
Figure 14-1: Menu.....	104
Figure 14-2: Readings	104
Figure 15-1: Calibrate temperature sensors display	106
Figure 15-2: Calibrate pressure sensors display.....	107
Figure 16-1 : unlock control panel.....	108
Figure 16-2: Lock panel in open position	108
Figure 16-3 : Menu functions	108
Figure 16-4: Maintenance functions.....	108
Figure 16-5: replace door seal	108
Figure 16-6: replace door seal	109
Figure 16-7: sliding clutch on door motor.....	110
Figure 17-1: Ambient pressure display	113
Figure 17-2: compressed air pressure reducing valve V300.....	114
Figure 17-3 pilot valve battery.....	114

22 List of Tables

Table 2-1: Sterilizer sizes.....	17
Table 2-2: Utilities for chamber size 5-5-9	18
Table 2-3: Recommended feed water q quality for sterilizers	22
Table 3-1: Size and Designation.....	24
Table 4-1: Pressure indicators in service area.....	28
Table 6-1 Factory set cycles	43
Table 6-2: 270°F Pre-vacuum Cycles.....	44
Table 6-3: 250°F Liquid Cycle	46

Table 6-4: 270°F Express Cycle	47
Table 6-5: 270°F Flash Cycles.....	48
Table 6-6: DART / Bowie-Dick Test cycle	49
Table 6-7: Warm-up & Leak Test cycle	50
Table 7-1: Recommended solid cycles	55
Table 7-2: Recommended solid load.....	55
Table 7-3: 250°F Liquid Cycle.....	57
Table 7-4: Recommended liquid loads.....	57
Table 7-5: 270°F Flash Cycles.....	57
Table 7-6: 270°F Express Cycle	57
Table 8-1: User menus and rights.....	62
Table 8-2: Alarm list	73

CIRCUIT DIAGRAM



Belimed
Infection Control

Standard: UL 61010A-1
 Project manager: P. Herde
 Connection data:
 Rated voltage: 3x208V, 60Hz
 Rated current/achievement: 6.0A / 1.6kW
 Extra low voltage: 24V/DC
 Back-up fuse: 20A
 Rated insulation voltage: 1000V
 Short circuit capacity: 10kA
 Mains system: TT

Series circuit diagram
MST-V 5-5-9

Belimed
Infection Control

MST-V 5-5-9_V1.2

Title sheet

No.	33
Sheet	1
Quantity	1

Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

Date	11.12.2009
User	AFA
Project	FEB
Secret	SECRET
Drawn	BECDL/EN/0004
Checked	
Released	

For more documents and data, please contact your Belimed representative.
 Belimed Sauter AG, CH-8583 Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure without the express written permission of Belimed Sauter AG is strictly prohibited.

Nous réservons tous droits dans ce document et les informations qu'il contient. Toute réimpression, utilisation ou divulgation sans autorisation expresse écrite de Belimed Sauter AG est strictement interdite.

Alle Rechte vorbehalten. Nachdruck, Vervielfältigung und Verbreitung, auch auszugsweise, ist ohne schriftliche Genehmigung der Belimed Sauter AG.


<p><u>Terminal strips</u></p> <p>Control panel</p> <ul style="list-style-type: none"> -XT1 Supply lead -XT2 Three-phase motor -XT3 Consumer 208V -XT9 Measuring terminals -XT11 Power supply 24V DC, emergency-off 	<p>On the plant</p> <ul style="list-style-type: none"> -XT23 Valves -XT24 Sensors -XT29 Measuring terminals 	<p><u>Structure of the destination device</u></p> <p>/ 21. A7</p> <ul style="list-style-type: none"> Separator Sheet number Plan section <p><u>Legend of operating material identification</u></p> <p>+SS -K119</p> <ul style="list-style-type: none"> Assembly place Class of operating material Counting number of a sheet Sheet number 	<p><u>Assembly place</u></p> <ul style="list-style-type: none"> +B1 Operating Panel side 1 +KA Chamber +KR Installation rack chamber +SS Control panel +EXT External 	<p><u>Line colors by IEC 204-1</u></p> <p>Main wiring system 3x208V black</p> <p>Protection- potential equalization wire green-yellow</p> <p>Control circuit 208 V without transformer black</p> <p>Control circuit with transformer < 48V green</p> <p>Control circuit with transformer > 48V brown</p> <p>Continuous current + red</p> <p>- pink</p> <p>blue</p> <p>blue-white</p>	<p><u>Line cross section</u></p> <p>Cross section of all no-designated lines: AWG19</p> <p><u>Modification administration</u></p> <p>Description</p> <p>Date</p> <p>Signature</p>	<p>Belimed Infection Control</p> <p>Legend</p> <p>Series circuit diagram MST-V 5-5-9</p> <p>MST-V 5-5-9_V1.2</p> <table border="1"> <tr> <td>Date</td> <td>11.12.2009</td> </tr> <tr> <td>User</td> <td>AFA</td> </tr> <tr> <td>Posid</td> <td>REB</td> </tr> <tr> <td>Datedel</td> <td>IEC204-EN6204</td> </tr> <tr> <td>Name</td> <td></td> </tr> <tr> <td>Change</td> <td></td> </tr> <tr> <td>Sheet</td> <td>3</td> </tr> <tr> <td>Blatt</td> <td>4</td> </tr> </table>	Date	11.12.2009	User	AFA	Posid	REB	Datedel	IEC204-EN6204	Name		Change		Sheet	3	Blatt	4
Date	11.12.2009																					
User	AFA																					
Posid	REB																					
Datedel	IEC204-EN6204																					
Name																						
Change																						
Sheet	3																					
Blatt	4																					

Für dieses Dokument und den darin dargestellten Gegenstand behält Belimed Sauter AG alle Rechte vor. Vervielfältigung, Verbreitung, auch auszugsweise, ist ohne schriftliche Genehmigung von Belimed Sauter AG, CH-Sulgen.

We reserve all rights in this document and in the information it contains. Reproduction, use or disclosure of this document without express authority is strictly forbidden.


Nous réservons tous droits pour ce document et les informations que il contient. Réimpression, utilisation, diffusion sans le consentement écrit de Belimed Sauter AG, CH-Sulgen est strictement interdite.

Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00




Belimed Sauter AG
+41 71 644 85 00
Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Verbreitung, auch auszugsweise, ist ohne schriftliche Genehmigung Belimed Sauter AG, CH-Sulgen.
We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden.
© Belimed Sauter AG, CH-Sulgen
est strictement interdite
© Belimed Sauter AG, CH-Sulgen.
Nous réservons tous droits pour ce document et les informations qu'il contient. Réimpression, utilisation, divulgation ou autre sans autorisation expresse Belimed Sauter AG, CH-Sulgen.

Date	11.12.2009
User	AFA
Project	FIEB
Standard	IEC60204/EN60204
Change	
Date	
Name	
Series circuit diagram	
MST-V 5-5-9	
Control cabinet layout	
MST-V 5-5-9 V1.2	
Sheet	4
Drawing	10



Status		Change	Date	Name	Standard	1
Belimed Saiter AG CH-8583 Sulgen +41 71 644 85 00						2
Für diese Dokumente und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Verantwortlich für die Inhalte sind die Autoren. Die Rechte an den Inhalten dieser Dokumente sind Eigentum der Belimed Saiter AG. Die Rechte an den Inhalten dieser Dokumente sind Eigentum der Belimed Saiter AG. Die Rechte an den Inhalten dieser Dokumente sind Eigentum der Belimed Saiter AG.						3
Belimed Saiter AG Infektion Control		Series circuit diagram MST-V 5-5-9		MST-V 5-5-9_V1.2		4
Belimed		Supply		MST-V 5-5-9_V1.2		5
-SS		Sheet 10		Following 15		6
						7
						8




Belimed Sauter AG CH-8583 Sulgen +41 71 644 85 00		Date: 11.12.2009 User: AFA		Series circuit diagram MST-V 5-5-9		Door motor		MST-V 5-5-9_V1.2	
Change	Date	Name	Standard	IEC60470/EN60204	REB	+SS			
Sheet				15			8		
Blanking				16			16		


We reserve all rights in this document and the information contained therein, express authority is strictly limited to third parties without express authorisation. No part of this document may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or by any information storage and retrieval system, without the prior written permission of Belimed Sauter AG, CH-Sulgen.

Nous réservons tous droits pour ce document et les informations qu'il contient. Toute réimpression, utilisation, reproduction ou utilisation dans un système de stockage et de récupération d'informations sans autorisation expresse de Belimed Sauter AG, CH-Sulgen, est strictement interdite.

© Belimed Sauter AG, CH-Sulgen



Status	Change	Date	Name	Standard	IEC204/EN6204
CH-8583 Sulger AG +41 71 644 85 00		Date	11.12.2009	Umr	AFA
Belimed Sauter AG		Preced	REB		
Series circuit diagram MST-V 5-5-9					
Belimed Infection Control		Vacuum pump		MST-V 5-5-9_V1.2	
					+SS
		Sheet	16		
		Numbering	30		




On/Off Key-switch

emergency-stop S2

Watch dog

Single door JZ on

Belimed Sauter Ag CH-8583 Sulgen +41 71 644 85 00	Date: 11.12.2009 User: AFA Prepared: REB Standard: IEC204-EN60204	Series circuit diagram		MST-V 5-5-9		Power supply		MST-V 5-5-9_V1.2		+SS	
		Revision	Date	Name	Date	Change	Date	Name	Date	Sheet	Total
Belimed Sauter Ag, CH-8583 Sulgen Für dieses Dokument und den darin dargestellten Gegenstand behält Belimed Sauter Ag alle Rechte vor. Verweigerung, Bekanntheit an Dritte oder Verwertung seines Inhaltes sind ohne unsere ausdrückliche Zustimmung verb.											
We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden.											
Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, diffusion à une tierce partie sans autorisation expresse est strictement interdite.											
© Belimed Sauter Ag, CH-8583 Sulgen											



Belimed Sauter AG
+41 71 644 85 00
CH-8583 Sulgen

Status	Change	Date	11.12.2009
Name	REB	User	AFB
Drawn	REB	Drawn	AFB
Checked	REB	Checked	AFB
Approved	REB	Approved	AFB
Released	REB	Released	AFB

Series circuit diagram
MST-V 5-5-9

Belimed
Infection Control

PLC, overview

MST-V 5-5-9_V1.2

-SS-

Sheet 37
Drawing 38

8

7

6

5




4

3

2

1

Option ICS 8535



0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
---	---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----

A

B

C

D

E

F

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48

49

50

51

52

53

54

55

56

57

58

59

60

61

62

63

64

65

66

67

68

69

70

71

72

73

74

75

76

77

78

79

80

81

82

83

84

85

86

87

88

89

90

91

92

93

94

95

96

97

98

99

100

101

102

103

104

105

106

107

108

109

110

111

112

113

114

115

116

117

118

119

120

121

122

123

124

125

126

127

128

129

130

131

132

133

134

135

136

137

138

139

140

141

142

143

144

145

146

147

148

149

150

151

152

153

154

155

156

157

158

159

160

161

162

163

164

165

166

167

168

169

170

171

172

173

174

175

176

177

178

179

180

181

182

183

184

185

186

187

188

189

190

191

192

193

194

195

196

197

198

199

200

201

202

203

204

205

206

207

208

209

210

211

212

213

214

215

216

217

218

219

220

221

222

223

224

225

226

227

228

229

230

231

232

233

234

235

236

237

238

239

240

241

242

243

244

245

246

247

248

249

250

251


252


253

254

255

256





	Belimed	PLC, digital outputs	MST-V 5-5-9 V1.2
Series circuit diagram	Infection Control		
MST-V 5-5-9			
Date: 11.12.2009	Layer: AFA	Project: FEB	Sheet: 61
Drawn: [Signature]	Name: [Signature]	Checked: [Signature]	Following: 62


Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. **Verweigerung** Bekanngabe an Dritte oder Verwertung seines Inhalts ist ohne unsere ausdrückliche Zustimmung verb. We reserve all rights in this document and in the relevant content thereof. **Refusal** to third parties without express authorisation is strictly forbidden.

© Belimed Sauter AG, CH-Sulgen
© Belimed Sauter AG, CH-Sulgen
© Belimed Sauter AG, CH-Sulgen

Neus reverts tous droits pour ce document et les informations que il contient. **Reproduction, utilisation, divulgation** à une tierce partie sans autorisation expresse est strictement interdite.

© Belimed Sauter AG, CH-Sulgen



Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

Für dieses Dokument und den darin dargestellten Gegenstand behälten wir uns alle Rechte vor. Vervielfältigung, Reproduktion an Dritte oder Verwertung seines Inhaltes ist ohne unsere ausdrückliche Zustimmung verb. © Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly prohibited. © Belimed Sauter AG, CH-Sulgen

Nous réservons tous droits pour ce document et les informations qu'il contient. Reproduction, utilisation et diffusion à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-Sulgen

Change	Date	Name	Standard
	11.12.2009	AFA	FREQ
			IEC60470:EN60204

Series circuit diagram
MST-V 5-5-9

Belimed
Inflection Control


PLC, digital outputs

MST-V 5-5-9 V1.2

Sheet	of
1	1

refer to schematic steam generator page 63

Option steam generator V30



Date		1.12.2009	Series circuit diagram		MST-V 5-5-9_V1.2		Sheet 63	
User		AFA	MST-V 5-5-9		-SS		Following 80	
Proved		FEB	PLC, digital outputs					
Standard		IEC204EN60204	Belimed		Infection Control			
Change	Date	Name						


Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Reproduktion, Verbreitung oder die öffentliche Zugänglichkeit ist ohne schriftliche Zustimmung von Belimed Sauter AG, CH Sulgen, Schweiz, ist ausdrücklich untersagt.

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden.

Nous réservons tous droits pour ce document et les informations que il contient. Réimpression, utilisation, diffusion ou l'accès par tiers sans autorisation expresse est formellement interdite.

© Belimed Sauter AG, CH Sulgen



Belimed Sauter AG
CH-8583 Sülgen
+41 71 644 85 00


Change	Date	Name	Status
	11.12.2009	AFA	Final

Series circuit diagram
MST-V 5-5-9

Belimed
Infection Control
PLC, analogue inputs


MST-V 5-5-9_V1.2

Sheet	80
Belimed	81




Belimed Sauter AG +41 71 644 85 00	Date: 11.12.2009	User: AFA	Project: REB	Name: IEC2004EN60204	Series circuit diagram		MST-V 5-5-9		PLC, analogue inputs		MST-V 5-5-9 V1.2		+SS	Sheet number: 81	8
					MST-V 5-5-9		PLC, analogue inputs		MST-V 5-5-9 V1.2						

Belimed
Infection Control




Status	Change	Date	Name	Prepared	REB	IEC204/EN60204
		11.12.2009	APA	REB		
Series circuit diagram MST-V 5-5-9						
Belimed Infection Control						
Operating Panel						
MST-V 5-5-9_V1.2						
						+SS
						Sheet 95
						Library 95



Status	Change	Date	Name	Drawn	Standard	IEC2004/EN60204
				AFB	REB	
Belimed Sauter AG CH-8583 Sulgen +41 71 644 85 00		Date		11.12.2009		
		User		AFA		
		Date				
		Name				
		Series circuit diagram				
		MST-V 5-5-9				
		Potential-free contacts				
		MST-V 5-5-9_V1.2				
		+SS				
		Sheet		96		
		Drawing		1103		

Me reserve all rights in this document and in the drawings. No change therein, reproduction, use or disclosure to third parties without express authority is strictly forbidden.
 © Belimed Sauter AG, CH-Sulgen
 Nous reservons tous droits pour ce document et les illustrations que il contient. Reproduction, utilisation et diffusion à tierce partie sans autorisation expresse est strictement interdite.
 © Belimed Sauter AG, CH-Sulgen

Für dieses Dokument und den darin dargestellten Gegenstand behält wir uns alle Rechte vor. Nachdruck, Verbreitung oder Weitergabe an Dritte ohne ausdrückliche Zustimmung verb. ist ohne unsere ausdrückliche Zustimmung verb.
 © Belimed Sauter AG, CH-Sulgen
 Belimed behält alle Rechte vor. Nachdruck, Verbreitung oder Weitergabe an Dritte ohne ausdrückliche Zustimmung verb. ist ohne unsere ausdrückliche Zustimmung verb.
 © Belimed Sauter AG, CH-Sulgen



Belimed Saurer AG CH-8583 Sülgen +41 71 644 85 00	11.12.2009 User AFA	Series circuit diagram MST-V 5-5-9	Belimed Infection Control	MST-V 5-5-9_V1.2	+SS	Sheet 103	8
Status	Change	Date	Name	Startat.	IECC04EN0004	108	

terminal strip

-XT2

		target external			potential		terminal number		jumper		target internal		Placement in Schematics		cable descriptions		
PE	PE	+KA -X154	4	PE	PE								/15.D4				
	PE	+KA -X156	4	#1533	PE	*							/15.D4				
	PE	+KA -M162	PE	PE	PE	PE							/15.D6				
					PE	*							/16.D2				
					#1614	PE				+SS -A381	-X30:3			/30.D2			
					#1454	*				+SS -G302	PE			/30.D2			

Terminal type

terminal ZPE2,5/4AN | PE - PE | _____

terminal _____

terminal _____

MST-V 5-5-9_V1.2



Series circuit diagram

MST-V 5-5-9

Date	11.12.2009
Umr	AFA
Prod	REB
Standard	IEC304-EN6204

Change	Date	Name
Status		

Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

terminal strip

-XT8

target external	potential	terminal number	jumper	target internal	Placement in Schematics	cable designation	cable type
	#1264	371			/37.C2		
	#1262	*		+SS -A370	/36		NS
	#1262			+SS -R370			
		371			/37.C2		
	#1265	372			/37.C2		
	#1263	*		+SS -A370	/37		JS
	#1263			+SS -R370			
		372			/37.C2		

Terminal type

terminal ZDU2.5/4AN | 371 - 372 |

terminal _____

terminal _____

MST-V 5-5-9_V1.2

Belimed Contents

Series circuit diagram

MST-V 5-5-9

Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

Für dieses Dokument sind den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. **Vervielfältigung, Bekanntgabe an Dritte oder Verwertung seines Inhaltes** sind ohne unsere ausdrückliche Zustimmung verb. © Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden. © Belimed Sauter AG, CH-Sulgen

Tous réservoirs tous droits pour ce document et les informations que il contient. Reproduction, utilisation, divulgation à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-Sulgen

terminal strip

-XT26

	target external	potential	terminal number	jumper	target internal	Placement in Schematics	cable designation
	+SS-X35	:8	#1593	35		+SS/43.B7	
	+SS-X35	:1	#1596	431		+SS/43.B2	
	+SS-X35	:2	#1597	432		+SS/43.B2	
	+SS-X35	:3	#1598	433		+SS/43.B3	
	+SS-X35	:4	#1599	434		+SS/43.B4	
	+SS-X35	:5	#1595	435		+SS/43.B5	
	+SS-X35	:6	#1594	437		+SS/43.B6	
	+SS-X35	:10	#1605	632		+SS/63.D2	
	+SS-X35	:11	#1606	633		+SS/63.D3	
	+SS-X35	:9	#1607	634		+SS/63.D4	
	+SS-X35	:12	#1604	635		+SS/63.D5	
	+SS-X35	:15	#1603	15		+SS/63.D6	
	+SS-X35	:13	#1602	636		+SS/63.D6	
	+SS-X35	:14	#1600	637		+SS/63.D7	
	+SS-X35	:7	#1608	638		+SS/63.D8	
PE	+SS-X35	:16	#1601	PE		+SS/63.D7	

Terminal type terminal 1608510000 [15 - PE]
 terminal _____
 terminal _____

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Veröffentlichung, Bekanntgabe an Dritte oder Verwertung seines Inhaltes sind ohne unsere ausdrückliche Zustimmung verb. © Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden. © Belimed Sauter AG, CH-Sulgen

Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, divulgation à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-Sulgen

Belimed Sauter AG
 CH-8583 Sulgen
 +41 71 644 85 00

Belimed Contents

Infection Control

Series circuit diagram
 MST-V 5-5-9

MST-V 5-5-9 V1.2

+ELD

Sheet 109
 Drawing 110

terminal strip

-XT29

		target external	potential	terminal number	jumper	target internal	Placement in Schematics	cable designation	cable type
		+KA -R811	.1	#1580	811		+SS/81.E1		
		+SS -XT9	.811	#922	*		+SS/81.E1		
		+KA -R811	.2	#1581	812		+SS/81.E2		
		+SS -XT9	.812	COM	*		+SS/81.E2		
		+KA -R812	.1	#1579	813		+SS/81.E2		
		+SS -XT9	.813	DI_9	*		+SS/81.E2		
				COM	812		+SS/81.E4		
		+KA -R812	.2	#1582	814		+SS/81.E3		
		+SS -XT9	.814	DI_10	*		+SS/81.E3		
				COM	812		+SS/81.E4		
		+KA -R815	.1	#1583	816		+SS/81.E5		
		+SS -XT9	.816	#924	*		+SS/81.E5		
		+KA -R815	.2	#1584	812		+SS/81.E6		
				COM	*		+SS/81.E6		
		+KA -R817	.1	#1585	817		+SS/81.E7		
		+SS -XT9	.817	#925	*		+SS/81.E7		
		+KA -R817	.2	#1586	812		+SS/81.E8		
				COM	*		+SS/81.E8		
		+EXT -S338	.bn	#1617	20		+SS/33.D8		
		+SS -A381	-X34:12	#1491	*		+SS/33.D8		
		+EXT -S338	.bu	#1537	338		+SS/33.D8		
		+SS -A381	-X34:13	#1618	*		+SS/33.D8		

Terminal type terminal ZDK-1.5 [20 - 817]
 terminal _____
 terminal _____

MST-V 5-5-9_V1.2

Belimed Contents
 Infection Control

Series circuit diagram
 MST-V 5-5-9

Date 11.12.2009
 User AFA
 Project REB
 Standort IEC304/EN61204

Sheet 110
 of 109

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Bekanntgabe an Dritte oder Verwertung seines Inhaltes sind ohne unsere ausdrückliche Zustimmung verb. © Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden. © Belimed Sauter AG, CH-Sulgen

Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, divulgation à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-Sulgen

Belimed Sauter AG
 CH-8583 Sulgen
 +41 71 644 85 00

Cable List: MST-V 5-5-9_V1.2

		Cable	Type			Length (mm)	
		+SS-B1	-W031	Multinorm-AWM2587-6LPE0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A921	+		1	-A381	-X30.1		
-A921	-		2	-A381	-X30.2		
-S303	14		3	-A381	-X31.3		
-S303	13		4	-A381	-X31.4		
			5				
			6				
-A921	GND		PE	-A381	-X31.7		
		Cable	Type			Length (mm)	
		+SS-B2	-W032	Multinorm-AWM2587-11LPE0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-X32	1		1	-X60	1		
-X32	2		2	-X60	2		
-X32	3		3	-X60	3		
-X32	4		4	-X60	4		
-X32	5		5	-X60	5		
-X32	6		6	-X60	6		
-X32	7		7	-X60	7		
-X32	8		8	-X60	8		
-X32	9		9	-X60	9		
-X32	10		10	-X60	12		
-X32	11		11	-X60	13		
-X32	12		PE	-X60	14		
		Cable	Type			Length (mm)	
		+SS-ELD	-W035	Multinorm-AWM2587-15LPE0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-X35	1		1	-XT26	431		
-X35	2		2	-XT26	432		
-X35	3		3	-XT26	433		
-X35	4		4	-XT26	434		
-X35	5		5	-XT26	435		
-X35	6		6	-XT26	437		
-X35	7		7	-XT26	638		
-X35	8		8	-XT26	35		
-X35	9		9	-XT26	634		
-X35	10		10	-XT26	632		
-X35	11		11	-XT26	633		
-X35	12		12	-XT26	635		
-X35	13		13	-XT26	636		
-X35	14		14	-XT26	637		
-X35	15		15	-XT26	15		
-X35	16		PE	-XT26	PE		

Cable List: MST-V 5-5-9_V1.2

		Cable	Type			Length (mm)	
		+SS-KA	-W154	Multinorm-AWM2587-3LPE2.5-AWG16			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-K611	2		1	-X154	1		
-K611	4		2	-X154	2		
-K611	6		3	-X154	3		
-XT2	PE		PE	-X154	4		
		Cable	Type			Length (mm)	
		+SS-KA	-W156	Multinorm-AWM2587-3LPE2.5-AWG16			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-K613	2		1	-X156	1		
-K613	4		2	-X156	2		
-K613	6		3	-X156	3		
-XT2	PE		PE	-X156	4		
		Cable	Type			Length (mm)	
		+SS-KA	-W162	Multinorm-AWM2587-3LPE2.5-AWG16			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-Q162	2		1	-M162	U		
-Q162	4		2	-M162	V		
-Q162	6		3	-M162	W		
-XT2	PE		PE	-M162	PE		
		Cable	Type			Length (mm)	
		+SS-KA	-W331	Multinorm-AWM2587-3L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X33.8		1	-S331	1		
-A381	-X33.9		2	-S331	2		
-A381	-X33.10		3	-S331	3		
		Cable	Type			Length (mm)	
		+SS-KA	-W332	Multinorm-AWM2587-4L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X33.11		1	-S332	13		
-A381	-X33.12		2	-S332	14		
-A381	-X33.13		3	-S332	21		
-A381	-X33.14		4	-S332	22		
		Cable	Type			Length (mm)	
		+SS-KA	-W333	Multinorm-AWM2587-4L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X34.4		1	-S333	13		
-A381	-X34.5		2	-S333	14		
-A381	-X34.6		3	-S333	21		
-A381	-X34.7		4	-S333	22		
		Cable	Type			Length (mm)	
		+SS-KA	-W335	Multinorm-AWM2587-3L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X33.1		1	-S335	1		
-A381	-X33.2		2	-S335	2		
-A381	-X33.3		3	-S335	3		

Cable List: MST-V 5-5-9_V1.2

		Cable	Type			Length (mm)	
		+SS-EXT	-W337	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X34.10		1	-S337	bn		
-A381	-X34.11		2	-S337	bu		
		Cable	Type			Length (mm)	
		+SS-KR	-W338	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X34.12		1	-XT29	20		
-A381	-X34.13		2	-XT29	338		
		Cable	Type			Length (mm)	
		+SS-KA	-W342	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X33.15		1	-S342	21		
-A381	-X33.16		2	-S342	22		
		Cable	Type			Length (mm)	
		+SS-KA	-W343	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X34.8		1	-S343	21		
-A381	-X34.9		2	-S343	22		
		Cable	Type			Length (mm)	
		+SS-KA	-W344	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X33.4		1	-S344	1		
-A381	-X33.5		2	-S344	3		
		Cable	Type			Length (mm)	
		+SS-KA	-W345	Multinorm-AWM2587-3L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X34.1		1	-S345	1		
-A381	-X34.2		2	-S345	2		
-A381	-X34.3		3	-S345	3		
		Cable	Type			Length (mm)	
		+SS-KA	-W346	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X33.6		1	-S346	1		
-A381	-X33.7		2	-S346	3		
		Cable	Type			Length (mm)	
		+SS-SS	-W351	Volltron-Twist-CY-2*2*0.34			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A370	11		ws	-X351	2		
-A370	12		bn	-X351	3		
-A370	15		gn	-X351	5		
			gb				

Cable List: MST-V 5-5-9_V1.2

		Cable	Type			Length (mm)	
		+SS-B1	-W355	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X39.5		1	-P355	+		
-A381	-X39.6		2	-P355	-		
		Cable	Type			Length (mm)	
		+SS-SS	-W371	Profibuskabel			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT8	371		ws	-A370	36		
-XT8	372		bn	-A370	37		
		Cable	Type			Length (mm)	
		+SS-KA	-W607	Multinorm-AWM2587-3L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X37.1		1	-Y607	B		
-A381	-X37.2		2	-Y607	C		
-A381	-x30.3		3	-Y607	A		
		Cable	Type			Length (mm)	
		+SS-KA	-W616	Multinorm-AWM2587-3L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X37.10		1	-Y616	B		
-A381	-X37.11		2	-Y616	C		
-A381	-X37.12		3	-Y616	A		
		Cable	Type			Length (mm)	
		+SS-KA	-W618	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X37.4		1	-Y618	1		
-A381	-X37.5		2	-Y618	2		
		Cable	Type			Length (mm)	
		+SS-KA	-W621	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X37.6		1	-Y621	1		
-A381	-X37.7		2	-Y621	2		
		Cable	Type			Length (mm)	
		+SS-KA	-W627	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-A381	-X37.8		1	-Y627	1		
-A381	-X37.9		2	-Y627	2		

Cable List: MST-V 5-5-9_V1.2

		Cable	Type			Length (mm)	
		+SS-KR	-W800	Multinorm-AWM2464-7L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT9	813		1	-XT29	813		
-XT9	814		2	-XT29	814		
-XT9	812		3	-XT29	812		
-XT9	811		4	-XT29	811		
-XT9	816		5	-XT29	816		
-XT9	817		6	-XT29	817		
			7				
-W800	SH		SH1				
		Cable	Type			Length (mm)	
		+SS-KA	-W801	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT9	20		1	-B801	1		
-XT9	801		2	-B801	2		
		Cable	Type			Length (mm)	
		+SS-KA	-W802	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT9	802		1	-B802	1		
-XT9	803		2	-B802	2		
		Cable	Type			Length (mm)	
		+SS-KA	-W803	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT9	20		1	-B803	1		
-XT9	804		2	-B803	2		
		Cable	Type			Length (mm)	
		+SS-KA	-W805	Multinorm-AWM2587-2L0.5-AWG20			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT9	20		1	-B805	1		
-XT9	805		2	-B805	2		
		Cable	Type			Length (mm)	
		+SS-B1	-W922	Volltron-Twist-CY-2*2*0.34			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-X922	3		ws	-X923	3		
-X922	8		bn	-X923	8		
			gn				
			gb				
			SH1				
		Cable	Type			Length (mm)	
		+SS-EXT	-W999	Multinorm-AWM2587-3LPE4.0-AWG12			
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-F154	1		1	-Q101	1		
-F154	3		2	-Q101	3		
-F154	5		3	-Q101	5		
-XT1	PE		PE	-Q101	PE		

MST-V 5-5-9_V1.2_SS

Bill of Material

Quant.	Ident.No.	Description	Supplier	Device Names	Article number
1	4-502-7224-0	Flash-EPROM 512 kBytes /4MBit	SAIA	-A370	5601SOD-134
1	PCD2.F120	Functionmodul; PCD2.F120	SAIA	-A370	5601SOD-108
1	PCD2.F520	Functionmodul; PCD2.F520	SAIA	-A370	5601SOD-93
1	PCD2.M157	CPU enclosure; PCD2.M157; CPU	SAIA	-A370	5600SOD-12
1	5-5-x-2005-07	I/OPrint-Input	Sauter AG	-A381	850 22135
2	PCD2.E160	Digital Input; PCD2.E160; 16 Inputs	SAIA	-D371,-D372	5601SOD-138
2	PCD2.A460	Digital Output; PCD2.A460; 16 Outputs	SAIA	-D373,-D374	5601SOD-139
1	PCD2.W340	Analog Input; PCD2.W340; 8 Inputs	SAIA	-D378	5601SOD-150
1	FAZ-C10/3	Circuitbreaker 3/10A; Characteristic C	Möller AG	-F154	5230KLM-82
1	TPP100-24	Power supply 24V; 100W	Lambda	-G302	5320-26
4	DIL-EM-01-G-C	CONTACTOR; AC-3 4 KW/400 V, 1 NC, DC 24 V, 3-POLE, SIZE S00, SCREW CONNECTION	Möller AG	-K611,-K612,-K613,-K614	5250KLM-117
1	DIL-EM-10-G	CONTACTOR; AC-3 4 KW/400 V, 1 NO, DC 24 V, 3-POLE, SIZE S00, SCREW CONNECTION	Möller AG	-K628	5250KLM-81
1	SFC10	Ferrite core;SFC10	Sibalco	-L151	5900-258
1	20906	Printerinterface; 4 Languages	Sauter AG	-P352	5731-274
1	56.02.30	Buzzer Piezo 24V	Piezo	-P615	5214-10
1	ZE-6	Thermic overload relay; adjustable range of current is 4A-6A, auxiliary contactors 1NO +1NC	Moeller	-Q162	5240KLM-65
1	150 Ohm/1%/0.6W	Precision resistor 150 Ohm; 1%; 0.6 W	Dätwyler	-R370	5810-165
1	FM4/TS35	Mounting plate for shield terminal	Weidmüller	-W800	5971WDM-253
1	KLBÜ4-13.5SC	Shield terminal	Weidmüller	-W800	5971WDM253
1	BLZF3.5/12	Connector 12 contacts; female contacts 3.5	Weidmüller	-X32	5971WDM-328
1	BLZF3.5/16	Connector 16 contacts; female contacts 3.5	Weidmüller	-X35	5971WDM-305
1	ZPE6	2-CONDUCTOR EARTH TERMINAL BLOCK; CENTRE MARKING (ASYMMETRIC LEFT); FRONT-ENTRY; green-yellow	Weidmüller	-XT1	5970WDM-79
2	ZPE2.5/4AN	4-CONDUCTOR THROUGH TERMINAL BLOCK; CENTRE MARKING (ASYMMETRIC LEFT); FRONT-ENTRY; green-yellow	Weidmüller	-XT2(2)	5970WDM-109
2	ZDU2.5/4AN	4-CONDUCTOR THROUGH TERMINAL BLOCK; CENTRE MARKING (ASYMMETRIC LEFT); FRONT-ENTRY; grey	Weidmüller	-XT8(2)	5970WDM-72

Bill of Material

MST-V 5-5-9_V1.2_SS

Quant.	Ident-No.	Description	Supplier	Device Names	Article number
14	ZDU2.5	2-CONDUCTOR THROUGH TERMINAL BLOCK; CENTRE MARKING (ASYMMETRIC LEFT); FRONT- ENTRY; grey	Weidmüller	-XT9(14)	5970WDM-67
7	PS1E21162BQ	Pilot valve	Parker	-Y334,-Y622,-Y623,-Y624,-Y625, -Y626,-Y631	121000-04

Bill of Material MST-V 5-5-9_V1.2_KR

Quant.	Ident-No.	Description	Supplier	Device Names	Article number
6	ZDK-1.5	DOUBLE DECK TERMINAL BLOCK; THROUGH/THROUGH; FRONT-ENTRY	Weidmüller	-XT29	5970WDM-85

Bill of Material

MST-V 5-5-9_V1.2_EXT

Quant.	Ident-No.	Description	Supplier	Device Names	Article number
1	0172001	Main switch; 25 A/3p; 3NO; Compact construction; Rotary drive	Baco	-Q101	5212-03
1	0172185	Logic element	Baco	-Q101	5211-143
2	0172167	Clamp cover	Baco	-Q101(2)	5211-144
2	PML 6A7230VCA	Foot switch;1NO	EHS	-S337,-S338	5218-53

Bill of Material

MST-V 5-5-9_V1.2_ELD

Quant.	Ident-No.	Description	Supplier	Device Names	Article number
16	1608510000	Terminal Z-series _ ZDU 2.5; Feed-through terminal; 2.50 mm ² ; dark beige; tension spring connection	Weidmüller	-XT26	5970WDM-67

Bill of Material

MST-V 5-5-9_V1.1_B1

Quant.	Ident-No.	Description	Supplier	Device Names	Article number
1	TP177B	Panel;TOP5000	Sauter AG	-A921	5601SIE-77
1	FB-58-S	Ribbon to matrixprinter	Gfm	-P355	5731-260
1	NP-58	Paper to matrixprinter	Gfm	-P355	5731-261
1	NP183-200-S2R	Matrixprinter;NP183-200-S2 R	Gfm	-P355	5730-101
1	130.245.201/0700	Key switch 0-1; 1NO	RAFI	-S303	5218-59
1	500.100.145	Logic element, 1NO	RAFI	-S303	5211-137
1	505.510.935	Coupler	RAFI	-S303	5211-136

Bill of Material

MST-V 5-5-9_V1.1_B2

Quant.	Ident-No.	Description	Supplier	Device Names	Article number
1	BT2/04	Panel S2	Sauter AG	-A925	350 06157
1	BLZF3.5/14FSN	Connector 14 contacts; female contacts 3.5	Weidmüller	-X60	5971WDM-331

CIRCUIT DIAGRAM

Standard: UL508A

Project manager:

Connection data:

Rated voltage: 3~480V, 60Hz

Rated current/power: 37A / 31kW

Extra low voltage: 120V/AC+24V/DC

Largest motor FLA: 0.35A / 0.16kW

Branch circuit protection: 70A

Rated insulation voltage: 600V

Short circuit capacity: 5kA

Belimed Title page

Infection Control

Steam generator V30

Date 12.07.2011

AGL

REB

IEC204:EN42204

User

Provid

Name

Standard

Change

Date

Name

Standard

CH-8583 Sulgen

+41 71 644 85 00

Belimed Sauter AG

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Bekanntheit an Dritte oder Verwertung seines Inhalts ist ohne unsere ausdrückliche Zustimmung verboten.

© Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden.

© Belimed Sauter AG, CH-Sulgen

Belimed Sauter AG, CH-Sulgen

est strictement interdite.

© Belimed Sauter AG, CH-Sulgen

Nous réservons tous droits pour ce document et les informations qu'il contient. Reproduction, utilisation, divulgation à une tierce partie sans autorisation expresse est strictement interdite.










© Belimed Sauter AG, CH-Sulgen


Contents

Sheet	Description	date	User	Document type
1	Title page	12.07.2011	AGL	01-Title page
3	Legend	12.07.2011	AGL	03-Circuit diagram
4	Control cabinet layout	12.07.2011	AGL	03-Circuit diagram
5	Control cabinet layout	12.07.2011	AGL	03-Circuit diagram
6	Elements heating	12.07.2011	AGL	03-Circuit diagram
10	Supply	12.07.2011	AGL	03-Circuit diagram
11	Control 120V/AC steam generator	12.07.2011	AGL	03-Circuit diagram
12	Control 120V/AC steam generator	12.07.2011	AGL	03-Circuit diagram
43	Power supply 24V	12.07.2011	AGL	03-Circuit diagram
63	Power supply 24V	12.07.2011	AGL	03-Circuit diagram
105	Terminal diagram (-X100)	12.07.2011	AGL	04-Terminal diagram
106	Terminal diagram (-X126)	12.07.2011	AGL	04-Terminal diagram
107	Terminal diagram (-X191)	12.07.2011	AGL	04-Terminal diagram

Belimed Sauter AG CH-8583 Sulgen +41 71 644 85 00	Date	12.07.2011	Steam generator V30		Contents	ELD30JUL_V1.4
	User	AGL				
	Printed	FREE				
	Standard	IEC2004/EN62024				
Change	Date	Name				
Status			Infection Control			

<p><u>Terminal strips</u></p> <p>In the steam generator cabinet</p> <ul style="list-style-type: none"> -XT00 Supply -XT26 Power supply 24V -XT91 Power supply 120V 	<p>On the plant</p> <ul style="list-style-type: none"> -XT113 Temperature limiter 	<p><u>Structure of the destination device</u></p> <p>/ 21. A7</p> <p>Separator</p> <p>Sheet number</p> <p>Plan section</p>	<p><u>Legend of operating material identification</u></p> <p>+SS -K119</p> <p>Assembly place</p> <p>Class of operating material</p> <p>Counting number of a sheet</p> <p>Sheet number</p>	<p><u>Line cross section</u></p> <p>Cross section of all no-designated lines: 0.75mm², AWG19</p>	<p><u>Change management</u></p> <p>Description</p> <p>Date</p> <p>Signature</p>			<p>Belimed Sauter AG +41 71 644 85 00</p> <p>CH-8583 Sulgen</p> <p>Belimed Sauter AG, CH-8583 Sulgen Bekanntgabe an Dritte oder Verwendungszweck sind ohne unsere ausdrückliche Zustimmung verb to third parties without express authori is strictly forbidden. © Belimed Sauter AG, CH-8583 Sulgen</p> <p>Belimed Sauter AG, CH-8583 Sulgen We reserve all rights in this document and in the inform ation contained therein. Reproduction, use or disclosure to third parties without express authorisation is strictly forbidden. © Belimed Sauter AG, CH-8583 Sulgen</p> <p>Belimed Sauter AG, CH-8583 Sulgen Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, diffusion à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-8583 Sulgen</p>
<p><u>Assembly place</u></p> <ul style="list-style-type: none"> +ED Steam generator +ELD In the steam generator cabinet +KD Terminal box 	<p><u>Line colors by EN 60204-1</u></p> <ul style="list-style-type: none"> Main wiring system 3x208/480V black Protection- potential equalization wire green-yellow Control circuit 208 V without transformer black Control circuit with transformer < 48V green Control circuit with transformer > 48V brown Direct current + white - blue blue-white 	<p><u>Legend</u></p> <p>Belimed Infection Control</p>	<p>Steam generator V30</p>	<p>ELDV30UL_V1.4</p>	<p>12.07.2011</p> <p>User: AGL</p> <p>Prepared: REB</p> <p>Standard: IEC/EN 60204-1</p>	<p>Change</p> <p>Date</p> <p>Name</p>	<p>Sheet numbering</p> <p>3</p> <p>4</p>	<p>8</p>

1	2	3	4	5	6	7	8				
A	B	C	D	E	F	G	H				
				 <p style="text-align: center;">Turn off main switch before opening</p>	 <p style="text-align: center;">Danger overpressure! Switch off the machine and depressurize and empty the system before maintenance</p>	 <p style="text-align: center;">Danger of scalding! Switch off the machine and empty the system before maintenance</p>					
								<p style="text-align: center;">Steam generator V30</p>	<p style="text-align: center;">Control cabinet layout</p>	<p style="text-align: center;">ELD V30JUL_V1.4</p>	
Status	Change	Date	Name	Signature	Date	User	AGL	Date	AGL	Date	AGL
<p>Belimed Sauter AG +41 71 644 85 00 CH-8583 Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>	
<p>Belimed Sauter AG +41 71 644 85 00 CH-8583 Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>		<p>© Belimed Sauter AG, CH-Sulgen</p>	



ELDVS30UL_V1.4

Control cabinet layout

Belimed
Infection Control


Steam generator V30

Date	User	AGL	REB	Standard	IEC204/EN60204
12.07.2011					


Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

Steam generator

30KW



Date		12.07.2011		AGL		ELDV30UL_V1.4					
User		AGL		AGL		AGL					
Printed		FEB		FEB		FEB					
Standard		IEC204/EN60204		IEC204/EN60204		IEC204/EN60204					
Name		Name		Name		Name					
Date		Date		Date		Date					
Change		Change		Change		Change					
Belimed Sauter AG CH-8583 Sulgen +41 71 644 85 00				Steam generator V30				Belimed Infection Control			
Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Bekanngabe an Dritte oder Vernetzung seines Inhaltes sind ohne unsere ausdrückliche Zustimmung verb.				© Belimed Sauter AG, CH-Sulgen				© Belimed Sauter AG, CH-Sulgen			
We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden.				© Belimed Sauter AG, CH-Sulgen				© Belimed Sauter AG, CH-Sulgen			
Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, diffusion à une tierce partie sans autorisation expresse est strictement interdite.				© Belimed Sauter AG, CH-Sulgen				© Belimed Sauter AG, CH-Sulgen			



Belimed Supply
 Infection Control

Steam generator V30

ELDVS00UL_V1.4

Status	Change	Date	Name	Standard	Prepared	User	AGL	AGL	AGL	AGL	AGL
				IEC204/EN60204							


Date	12.07.2011
User	AGL
Prepared	REB
Standard	IEC204/EN60204

FIELD WIRING:
 USE COPPER WIRE ONLY,
 MIN. 75°C 600V AWG4,
 TORQUE VALUE: Marked
 on Device 1,7 lbs (2,3Nm)


Branch circuit protection max. 70A
 to be provided by customer
 ICU max. 5KA

Step 1: 30kW

Sheet	8
Labeling	10
	11



Belimed Sauter AG		CH-8583 Sulgen		+41 71 644 85 00	
Date	12.07.2011	Date	12.07.2011	Date	12.07.2011
User	AGL	User	AGL	User	AGL
Drawn	PREB	Drawn	PREB	Drawn	PREB
Standard	IEC60364-EN60204	Standard	IEC60364-EN60204	Standard	IEC60364-EN60204
Sheet	11	Sheet	11	Sheet	11
Change	2	Change	2	Change	2
Steam generator V30		Control 120V/AC steam generator		ELD30UL_V1.4	
Infestation Control		Infestation Control		Infestation Control	




Belimed Sauter AG
 CH-8583 Sülgen
 +41 71 644 85 00

Für diesen Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Wiederverbreitung, Kopie, Reproduktion, Verbreitung oder sonstiger Gebrauch ist ohne schriftliche Genehmigung von Belimed Sauter AG, CH-Sülgen.


We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden.

Belimed Sauter AG, CH-Sülgen
 est établissement industriel
 est établissement industriel
 est établissement industriel



Belimed Sauter AG CH-8583 Sulgen +41 71 644 85 00		Date: 12.07.2011 User: AGL		Steam generator V30		Belimed Power supply 24V		ELDY30UL_V1.4		Sheet 43 of 63	
Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Verbreitung, Reproduktion, Nutzung oder disclosure ohne unsere ausdrückliche Zustimmung verb.		We reserve all rights in this document and the information contained therein. Reproduction, use or disclosure without express authority is strictly forbidden.		© Belimed Sauter AG, CH-Sulgen		© Belimed Sauter AG, CH-Sulgen		© Belimed Sauter AG, CH-Sulgen		© Belimed Sauter AG, CH-Sulgen	
Change	Date	Name	Standard	IEC204/EN60204							

refer to circuit diagram steam sterilizer page 63



Status	Change	Date	Name	Signature	IEC2004EN02004
Steam generator V30			Belimed		
Power supply 24V			ELD30UL_V1.4		
+41 71 644 85 00			Belimed Sauter AG, CH-Sülgen		
12.07.2011			12.07.2011		
User: AQL			User: AQL		
Prep: REB			Prep: REB		
Date			Date		
Name			Name		
Signature			Signature		
Sheet 03			Sheet 03		
Drawing 105			Drawing 105		

terminal strip

-XT00

*ELD-ED-W107
cable designation

PE +ED-M107 :PE #1529 PE

target
external

potential

terminal number

jumps

target
internal

Placement in
Schematics

cable designation

+Customer-Q101	:1	#1975	1			
+Customer-Q101	:3	#1976	2			
+Customer-Q101	:5	#1977	3			
		#1847	\$		/I0.D1	
		#1847	PE		/I0.D1	
+Customer-Q101	:PE	#1847	\$		/I0.D1	
		#1847	PE		/I0.D1	
+ED-R101	:PE	#1524	PE		/I0.D2	
+ED-R102	:PE	#1524	PE		/I0.D2	
			\$			
PE +ED-M107	:PE	#1529	PE		/I0.D2	

+ELD-K113	:2				/I0.D1	
+ELD-K113	:4				/I0.D1	
+ELD-K113	:6				/I0.D1	

Terminal type

terminal 1020500000 [1-3]

terminal 0661460000 [\$,\$,PE,PE]

terminal 7907400000 [\$,PE]

Belimed
Infection Control

Steam generator V30

12.07.2011
AGL
REB
IEC20M/EN60204

Date Name
Printed Standard

Charge Status

ELD
+ELD

ELD V30UL_V1.4
Sheet 105
Following 106

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Bekanntgabe an Dritte oder Verwertung seines Inhaltes, sind ohne unsere ausdrückliche Zustimmung verb. © Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden. © Belimed Sauter AG, CH-Sulgen

Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, divulgation à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-Sulgen

Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

terminal strip

-XT26

		target external				jumper		target internal				
			potenti	terminal number				Placement in Schematics				
1	1	+SS -X46	:2	#2021	431	•	+ELD -K119	:14	/43.D2	•		
		Verbindung ist nicht eindeutig.					#2021		/43.D2			
		+SS -X35	:2	#2026	432		+ELD -K114	:14	/43.D2			
		+SS -X35	:3	#2022	433		+ELD -K112	:14	/43.D3			
		+SS -X35	:4	#2023	434		+ELD -K122	:6	/43.D4			
		+SS -X35	:5	#2019	435		+ELD -K116	:14	/43.D5			
		Verbindung ist nicht eindeutig.					#2019		/43.D5			
		+SS -X35	:6	#2024	437		+ELD -Q107	:14	/43.D7			
		+SS -X35	:8	24V_ELD	35		+ELD -Q107	:13	/43.D8			
		24V_ELD							/43.D8			
		+SS -X35	:10	#2039	632				/63.B2			
		+ED -Y632	:1	#2039					/63.B2			
		24V_ELD					35		/63.D1			
		24V_ELD							/63.D1			
		Verbindung ist nicht eindeutig.					#2040		/63.B3			
		+SS -X35	:11	#2040	633				/63.B3			
+ED -Y633	:1	#2040				/63.B3						
+SS -X35	:9	#2029	634	+ELD -K634	:A1	/63.B4						
+SS -X35	:13	#2031	636	+ELD -K636	:A1	/63.B6						
+SS -X35	:12	#2030	635	+ELD -K635	:A1	/63.B5						
+SS -X35	:7	#2033	638	+ELD -K638	:A1	/63.B8						
+SS -X35	:14	#2032	637	+ELD -K637	:A1	/63.B7						
+SS -X35	:16	#2038	PE			/63.B6						
+SS -X35	:15	0V_ELD	15	Verbindung ist nicht eindeutig.			/63.B5					
+ED -Y632	:2	0V_ELD	15	Verbindung ist nicht eindeutig.			/63.D2					
+ED -Y633	:2	0V_ELD	15	Verbindung ist nicht eindeutig.			/63.D3					
		0V_ELD	15	+ELD -M631			/63.D1					
		0V_ELD	15	Verbindung ist nicht eindeutig.			/63.D1					

Terminal type terminal 1791100000 [15 - PE]
 terminal _____
 terminal _____

Belimed Terminal diagram

Steam generator V30

ELDV30UL_V1.4

+ELD

Infection Control

Date	12.07.2011
User	AGL
Project	REB
Specialist	IEC2004/EN02004
Change	
Date	
Name	
Status	

Sheet 106
Drawing 107

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Bekanntgabe an Dritte oder Verwertung seines Inhaltes sind ohne unsere ausdrückliche Zustimmung verb. © Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden. © Belimed Sauter AG, CH-Sulgen

Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, divulgation à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-Sulgen

Belimed Sauter AG
 CH-8583 Sulgen
 +41 71 644 85 00

terminal strip

-XT91

MultiForm-AMM2587-3L PE-A
MultiForm-AMM2587-3L PE-D
MultiForm-AMM2587-3L PE-O
cable designation

+ELD-ED-W115
+ELD-ED-W113
+ELD-ED-W112
cable designation

PE
PE
PE

target external

potential

terminal number

jumpers

target internal

Placement in Schematics

cable designation

+ED-S112	I.1	#2004	112
+ED-S112	II.4	#4293	113
+ED-S112	II.2	#2099	114
+ED-S113	I.1	#2099	115
+ED-S113	I.2	#4292	116
+ED-S113	I.4	#4293	117
+ED-S115	I.1	#2013	118
+ED-S115	I.2	#2100	119
		#2100	PE
		N_ELD	PE
		N_ELD	
		N_ELD	
+ED-S112	I. PE	#1999	PE
+ED-S113	I. PE	#2098	PE
+ED-S115	I. PE	#2011	PE

+ELD-K635	I.11	/I1.C2
+ELD-K112	I.A1	/I1.D2
		/I1.D3
		/I1.D3
+ELD-K125	I.1	/I1.B3
		/I1.B3
+ELD-K638	I.22	/I1.B5
+ELD-K636	I.14	/I1.C5
+ELD-K119	I.A1	/I1.C5
+ELD-T111	I. PE	/I1.A1
+ELD-T111	I.0	/I1.A1
Verbindung ist nicht eindeutig.		/I1.A1
		/I1.C2
		/I1.D3
		/I1.B5

Terminal type

terminal 1608570000 [I12 - I19]
terminal 1608660000 [PE - PE]
terminal

Für dieses Dokument und den darin dargestellten Gegenstand behalten wir uns alle Rechte vor. Vervielfältigung, Bekanntgabe an Dritte oder Verwertung seines Inhaltes sind ohne unsere ausdrückliche Zustimmung verb. © Belimed Sauter AG, CH-Sulgen

We reserve all rights in this document and in the information contained therein. Reproduction, use or disclosure to third parties without express authority is strictly forbidden. © Belimed Sauter AG, CH-Sulgen

Nous réservons tous droits pour ce document et les informations que il contient. Reproduction, utilisation, divulgation à une tierce partie sans autorisation expresse est strictement interdite. © Belimed Sauter AG, CH-Sulgen

Belimed Sauter AG
CH-8583 Sulgen
+41 71 644 85 00

ELDVS0UL_V1.4

Terminal diagram

Belimed
Infection Control

Steam generator V30

Date	12.07.2011
User	AGL
Project	REB
Standard	IEC29041ENB0204
Name	
Date	
Change	
Status	

+ELD

Sheet 107
Drawing

Cable List: ELDV30UL_V1.4

		Cable	Type			Length (mm)	
		+SS-ELD	-W035	Multinorm-AWM2587-17LPE0.75-AWG19		3800	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT26	431		1	-X46	2		
-XT26	432		2	-X35	2		
-XT26	433		3	-X35	3		
-XT26	434		4	-X35	4		
-XT26	435		5	-X35	5		
-XT26	437		6	-X35	6		
-XT26	638		7	-X35	7		
-XT26	35		8	-X35	8		
-XT26	634		9	-X35	9		
-XT26	632		10	-X35	10		
-XT26	633		11	-X35	11		
-XT26	635		12	-X35	12		
-XT26	636		13	-X35	13		
-XT26	637		14	-X35	14		
-XT26	15		15	-X35	15		
			16				
			17				
-XT26	PE		PE	-X35	16		
		Cable	Type			Length (mm)	
		+ELD-ED	-W107	Multinorm-AWM2587-3LPE1.0-AWG18		4500	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-K634	2		1	-M107	U		
-K634	4		2	-M107	V		
-K634	6		3	-M107	W		
-XT00	PE		PE	-M107	PE		
		Cable	Type			Length (mm)	
		+ELD-ED	-W112	Multinorm-AWM2587-3LPE1.0-AWG18		2250	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT91	112		1	-S112	I.1		
-XT91	114		2	-S112	II.2		
-XT91	113		3	-S112	II.4		
-XT91	PE		PE	-S112	PE		
		Cable	Type			Length (mm)	
		+ELD-ED	-W113	Multinorm-AWM2587-3LPE1.0-AWG18		2450	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT91	115		1	-S113	1		
-XT91	116		2	-S113	2		
-XT91	117		3	-S113	4		
-XT91	PE		PE	-S113	PE		
		Cable	Type			Length (mm)	
		+ELD-ED	-W115	Multinorm-AWM2587-2LPE1-AWG18		2050	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT91	118		1	-S115	1		
-XT91	119		2	-S115	2		
-XT91	PE		PE	-S115	PE		

Cable List: ELDV30UL_V1.4

		Cable	Type			Length (mm)	
		+ELD-ED	-W122	Multinorm-AWM2587-4L0.5-AWG20		2300	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-K122	5		1	-B122			
-K122	3		2	-B122			
-K125	11		3	-B125			
-K125	1		4	-B125			
		Cable	Type			Length (mm)	
		+ELD-ED	-W632	Multinorm-AWM2587-2L0.5-AWG20		2850	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT26	632		1	-Y632	1		
-XT26	15		2	-Y632	2		
		Cable	Type			Length (mm)	
		+ELD-ED	-W633	Multinorm-AWM2587-2L0.5-AWG20		1750	
Device Name	Pin	Fitting Part	Core	Device Name	Pin	Fitting Part	Length (mm)
-XT26	633		1	-Y633	1		
-XT26	15		2	-Y633	2		

Bill of Material **ELDV30UL_V1.4_ELD**


Quant.	Ident.-No.	Description	Supplier	Device Names	Article number
2	CC3/30A	Fuse holder 3/max. 30A	AMBUS	-F102,-F103	5221-04
2	HCLR25	Fuse class CC 25A	Edison	-F102,-F103	5220-24
1	CC2/1A	Circuitbreaker 2/1A	AMBUS	-F108	5221-05
1	HCTR1	Fuse class CC 1A	Edison	-F108	5220-25
3	55.32.8.120.0040	Relay; 120VAC; 2NC, 2NO	Finder	-K112,-K114,-K119	5264-15
3	94.82	Relay socket 8 pole	Finder	-K112,-K114,-K119	5261-10
1	DILM65(120V)	CONTACTOR; lth 80A; AC120V60HZ	Möller AG	-K113	5250KLM-137
2	DILM1000-XHIC11-SI	Auxiliary contact, 1NC, 1NO	Möller AG	-K113(2)	5251KLM-27
2	DILMC25-10(120V)	CONTACTOR; lth 35A; 1S; AC120V60HZ	Möller AG	-K115,-K116	5250KLM-133
1	16HMW1A00101	Level relay	Warrick	-K122	5263-03
1	EC-8	Relay socket 8 pole; EC-8	Comat	-K122	5261WIS-01
1	EC-11	Relay socket 11 pole; EC-11	Comat	-K125	5261-11
1	LLC6	Level relay	ABB	-K125	5217-30
1	DIL-EM-01-G-C	CONTACTOR, AC-3 4 KW/400 V, 1 NC, DC 24 V, 3-POLE, SIZE S00, SCREW CONNECTION	Möller AG	-K634	5250KLM-117
4	Haltebügel	Retainer to relays RT424024; 3-1393161-3	Schrack	-K635,-K636,-K637,-K638	5261-05
4	RT424024	Relay; 24VDC; 2NC, 2NO	Schrack	-K635,-K636,-K637,-K638	5264-03
4	Socket2Pol	Relay socket 2 pole; 2-1419108-2	Schrack	-K635,-K636,-K637,-K638	5261-04
1	FPF08KRD24BR	Fan; 24V; 105x105mm	Simpex	-M631	5160-13
1	FPF08KRG	Fan kit 105x105mm	Simpex	-M631	5161-08
1	BK25-E	Terminalblock	Möller AG	-Q107	5241KLM-14
1	PKZMO-0.4-C+NHI-11	CIRCUIT-BREAKER FOR MOTOR PROTECTION; 0.25-0.4A; 0.06kW400V	Möller AG	-Q107	5240KLM-52
1	EEJ175/EJ96/C59.75-380-3	Transformer 460V/480V -120V;175VA	Wagner&Grimm	-T111	5310-28
1	0661460000	Terminal EK _EK 35/35; PE terminal; 35.00 mm²; green/yellow; screw connection	Weidmüller	-XT00	5970WDM-24
2	7907400000	Terminal Z-series _ZPE 6/3AN; PE terminal; 6.00 mm²; green/yellow; tension spring connection	Weidmüller	-XT00(2)	5970WDM-99
3	1020500000	Terminal WDU _WDU 35; Feed-through terminal; 35.00 mm²; dark beige; screw connection	Weidmüller	-XT00(3)	5970WDM-86
5	1608660000	Terminal Z-series _ZPE 2.5/4AN; PE terminal; 2.50 mm²; green/yellow; tension spring connection	Weidmüller	-XT00,-XT91(4)	5970WDM-109

Bill of Material

ELD30UL_V1.4_ELD

Quant.	Ident-No.	Description	Supplier	Device Names	Article number
10	1791100000	Terminal Z-series _ZDK 1.5; Double terminal; 1.50 mm ² ; dark beige; tension spring connection	Weidmüller	-XT26(10)	5970WDM-85
8	1608570000	Terminal Z-series _ZDU 2.5/4AN; Feed-through terminal; 2.50 mm ² ; dark beige; tension spring connection	Weidmüller	-XT91(8)	5970WDM-72
1	26183	RC element; A0	Murr	-Z113	5810-75
1	DILM32-XSPR240	RC element	Möller AG	-Z115	5251KLM-25

27 P & I Diagram Sterilizer without cooling circuit



The copyright in this document shall inhere in Belimed SA. This documentation is personal and confidential to the customer. It may not be reproduced, stored in a retrieval system, or used in any way without the prior written approval of Belimed SA. Also it is not allowed to make it available to any third party, even in third parties.

Replaces by:	Replacement Part:	Height (opt.):
Project:	MST-V 5-5-9 VSI/WS2 USA	Scale:
Drawn by:	23.08.11	SSR
Checked by:	23.08.11	RRR
Released by:	01.10.12	TRU
Approved by:	01.10.12	PAL
Priority Water:	short	
Part Number:	60_300_17579	

Belimed Belimed Sauter AG
 CH-6883 Sulgen
 Tel.: +41 71 664 80 00
 Fax: +41 71 664 80 00

28 Bill of Material to P & I Diagram Sterilizer without cooling circuit

page: 1 / 3

LIST OF SPARE PARTS / PROCESS TECHNOLOGY


DESIGN.	EACH	TERM	DIMENSION	MATERIAL	TYPE	PART NO.	MANUFACTURER/SUPPLIER
A 502	1	Steam trap	1/2"	MS	BPT 13S SUB (F=24K	102000-48	Spirax Sarco
A 504	1	Steam trap	1/2"	MS	BPT 13S SUB (F=24K	102000-48	Spirax Sarco
B 801	1	Pressure transducer	1/2"	316L	MIDAS S05/0-4bar/4	5510JUM-17	Jumo
B 802	1	Pressure transducer	1/2"	316L	MIDAS S05/0-4bar/4	5510JUM-17	Jumo
B 803	1	Pressure transducer	1/4"	1.4571	MIDAS S05/0-6bar/4	5510JUM-19	Jumo
B 805	1	Pressure transducer	1/4"	1.4571	MIDAS S05/0-6bar/4	5510JUM-19	Jumo
F 200	1	Strainer	3/4" zu ø 1/2"	1.4404	03-000.5691	03-000.5691	Sauter Sul. /Künzli AG
F 300	1	Filter	1/4"	AL/ZN	Multi-Fix F.11	136001MF-07	Bachofen / Knocks
F 305	1	Filter	1/2"	KS	ZGPB.01V-E8	103001-82	Crustag/d.hunter
F 400	1	Strainer	1" zu ø 3/4"	1.4404	03-000.5690	03-000.5690	Sauter Sul. /Künzli AG
F 502	1	Strainer	3/4" zu ø 1/2"	1.4404	03-000.5691	03-000.5691	Sauter Sul. /Künzli AG
P4	1	Vacuum pump	-	GG-25	VZ 30	5130SPK-78	Speck /E.W. Müller
PI 100	1	Mano-vacuometer	1/4"	1.4301	MGS18.1/A-63 -1..+	111000D063-57	Moesch K. AG
PI 200	1	Manometer	1/4"	MS 58	MGS10.3/A-63/0..+1	111000D063-32	Nuova Fima
PI 300	1	Manometer	1/4"	MS	Multi-Fix G.50-10	136001MF-06	Bachofen / Knocks
PI 301	1	Mano-vacuometer	1/4"	MS	MGS10.3/A-63SEC/-1	111000D063-31	Nuova Fima
PI 302	1	Mano-vacuometer	1/4"	MS	MGS10.3/A-63SEC/-1	111000D063-31	Nuova Fima
PI 800	1	Mano-vacuometer	1/4"	1.4301	MGS18.A-63,-1,+5ba	111000D063-58	Nuova Fima
PI 804	1	Mano-vacuometer	1/4"	MS	MGS10.3/A63SEC/-1.	111000D063-33	Nuova Fima
R 300	1	Spring-loaded check valve	1/4"	AL/ZN	Multi-Fix V.11	136001MF-10	Bachofen / Knocks
R 305	1	Spring-loaded check valve	1/2"	1.4301	Small 0,02-0,03 ba	142A20015-02	Hamberger / Cromax
R 402	1	Spring-loaded check valve	1/4"	1.4301	Small 0,02-0,03 ba	142A20008	Hamberger / Cromax
R 623	1	Butterfly flap valve	1/8"x4.0mm	MS verni.	A327927	121001-26	Bachofen AG
R 624	1	Butterfly flap valve	1/8"x4.0mm	MS verni.	A327927	121001-26	Bachofen AG
R 811/812	1	Resistance thermometers	5.0mm	1.4571	2xPT1000 5.0/120/3	5550JUM-108	Jumo
R 815	1	Resistance thermometers	4.5mm	1.4571	1xPT1000/1L-4.5-20	5550JUM-19	Jumo
R 817	1	Resistance thermometers	4.5mm	1.4571	1xPT1000/1L-4.5-20	5550JUM-19	Jumo
S 331	1	Manometric switch	1/4"	MS	901.51 111L1 +2.6b	5215-25	Pewatron
S 335	1	Manometric switch	1/4"	MS verni.	901.11.122X1 +0.12	5215-30	Pewatron
S 345	1	Manometric switch	1/4"	MS	901.51 111L1 +2.6b	5215-25	Pewatron

DESIGN.	EACH	TERM	DIMENSION	MATERIAL	TYPE	PART NO.	MANUFACTURER/SUPPLIER
S 346	1	Manometric switch	1/4"	MS	901.51 111L1 +4.5b	5215-28	Pewatron / Beck GmbH
V 300	1	Pressure controller	1/4"	AL	R.11 / Serie 1 / M	136001MF-09	Bachofen / Knocks
V 334	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 607	1	Solenoid valve 4-way	1/4"	AL	135MF0/24VDC/13W	121000LU002-01	Bachofen / Lucifer
V 616	1	Solenoid valve 4-way	1/4"	AL	135MF0/24VDC/13W	121000LU002-01	Bachofen / Lucifer
V 618	1	Solenoid valve 2-way	1/4"	MS	EV210A 3.5B/24V DC	121MS0DA008-01	Danfoss
V 621	1	Solenoid valve 2-way	1/2"	MS	EV220B 15B/24V DC/	121MS0DA015-02	Danfoss
V 622	1	Slanted seat valve	1/2"	1.4408	2000-A-013,0-E-VAG	122A40BU015-01	Bürkert
V 622.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 623	1	Slanted seat valve	1/2"	1.4581	2000-A-015,0-E-VAS	122A40BU015	Bürkert
V 623.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 624	1	Slanted seat valve	3/4"	RG	2000-A-020,0-E-RGG	121000-04	Bürkert
V 624.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 625	1	Slanted seat valve	1/2"	1.4408	2000-A-013,0-E-VAG	122A40BU015-01	Bürkert
V 625.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 626	1	Slanted seat valve	1/2"	RG	2000-A-013,0-E-RGG	122RG0BU015	Bürkert
V 626.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 627	1	Solenoid valve 2-way	1/4"	1.4401	21L2KIT30-2/2Weg-1	121A0SA008	ODE
W4	1	Heat exchanger	1"	1.4401	SL 34-TL-BB-14 / A	04-002.7451	ETS
X 100	1	Orifice	G 3/4" zu ø 1/2"	1.4435	zum Aufbohren	750 05362	Marton / Zweifel
X 110	1	Orifice	zu ø 12.0mm	1.4435	D=14.5 d=5.0 s=1.5	750 05133	Marton / Zweifel
X 205	1	Tube	8.0mmx13.5mm	1.4435	8x14 mm PVC-transp	350 00727	Maagtechnic
X 221	1	Orifice	zu ø 15.0mm	1.4435	D=18.5 d=4.0 s=1.5	750 00052	Marton / Zweifel
X 300	1	Pipe	ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 309	1	Pipe	ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 309.1	1	End plate	6.0mm	KS	PS1-E102	121001-04	Girmatic / Parker
X 309.2	1	Sound absorber	6.0mm	KS	0671 06 00CH	121001-05	Parker
X 309.3	1	Ball	3.0mm	ST chem.verni.	HY219 488	121001-07	Girmatic / Parker
X 309.4	1	Sound absorber	6.0mm	KS	0671 06 00CH	121001-05	Parker
X 400	1	Compensator	1"	1.4404	Assiwell DN 25	106A40-61	Angst-Pfister

* DRAW NO. : 60_300.17579
 * PROJECT : MST-V 5-5-9 VSI CSD1 o. KK USA
 * DESIGNATION 1 : P & ID 5-5-9 VSI/VS2 CSD1
 * DESIGNATION 2 : MST-V
 * ORIGINATOR : GGR
 * DATE : 23.08.11
 * CHECKED : ABR
 * DATE : 23.08.11
 * REVISED : IRU
 * DATE : 01.10.12
 * CHECKED : PAL
 * DATE : 01.10.12
 * INDEX : B
 * FILE : MSTV_009_A
 * LIST TYPE: EFWTXA

DESIGN	EACH	TERM	DIMENSION	MATERIAL	TYPE	PART NO.	MANUFACTURER/SUPPLIER
X 400.1	1	Compensator	1"	1.4404	Assiwell DN 25	106A40-61	Angst+Pfister
X 400.2	1	Expansion joint	3/4"	1.4404	FMH840.020.25	300 02805	Torgen GmbH
X 402	1	Pipe	Ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 402.1	1	Pipe	Ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 405	1	Orifice	zu Ø 8.0mm	1.4435	D=10.0 d=0.6 s=1.5	350 07094	Marton / Zweifel
X 405.1	1	Pipe	Ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 500	1	Orifice	zu Ø 15.0mm	1.4435	D=18.5 d=2.0 s=1.5	750 00049	Marton / Zweifel
X 600	1	Compensator	1"	1.4404	Assiwell DN 25	106A40-61	Angst+Pfister
X 705	1	Sound absorber	1"	Sinterbronze	06700034-1"	121001-25	Legris
Y 700	1	Safety valve	3/4"-1"	RG	0041EDA-40 / Stamp	133RG000020-32	Spence/USA
Y 702	1	Safety valve	3/4"-1"	RG	0041EDA-40 / Stamp	133RG000020-32	Spence/USA

29 P & I Diagram Sterilizer with cooling circuit



Revised by:	Replacment for:	Weight:	kg/l:
Project:	MST-V 5-5-9 VS1/WS2 USA	Sheet:	IRU
	P & ID 5-5-9 VS1/WS2 CSD1	Issued:	23.09.11
	MST-V	Revised:	-
		Sheet:	60_300_17500

Belimed Sauter AG
Tel: +41(0)71 644 85 00
Fax: +41(0)71 644 86 00

30 Bill of Material to P & I Diagram Sterilizer with cooling circuit

LIST OF SPARE PARTS / PROCESS TECHNOLOGY

page: 1 / 3

DESIGN	EACH	TERM	DIMENSION	MATERIAL	TYPE	PART NO.	MANUFACTURER/SUPPLIER
A 502	1	Steam trap	1/2"	MS	BPT 13S SUB (F=24K	102000-48	Spirax Sarco
A 504	1	Steam trap	1/2"	MS	BPT 13S SUB (F=24K	102000-48	Spirax Sarco
B 801	1	Pressure transducer	1/2"	316L	MIDAS S05/0-4bar/4	5510JUM-17	Jumo
B 802	1	Pressure transducer	1/2"	316L	MIDAS S05/0-4bar/4	5510JUM-17	Jumo
B 803	1	Pressure transducer	1/4"	1.4571	MIDAS S05/0-6bar/4	5510JUM-19	Jumo
B 805	1	Pressure transducer	1/4"	1.4571	MIDAS S05/0-6bar/4	5510JUM-19	Jumo
F 200	1	Strainer	3/4" zu ø 1/2"	1.4404	03-000.5691	03-000.5691	Sauter Sul. /Künzli AG
F 220	1	Strainer	3/4"	Bronze	SF12.1	101RG0-04	Spirax Sarco
F 300	1	Filter	1/4"	AL/ZN	Multi-Fix F.11	136001MF-07	Bachofen / Knocks
F 305	1	Filter	1/2"	KS	ZGPB.01V-E8	103001-82	Crustag/d.hunter
F 400	1	Strainer	1" zu ø 3/4"	1.4404	03-000.5690	03-000.5690	Sauter Sul. /Künzli AG
F 502	1	Strainer	3/4" zu ø 1/2"	1.4404	03-000.5691	03-000.5691	Sauter Sul. /Künzli AG
P4	1	Vacuum pump		GG-25	VZ 30	5130SPK-78	Speck /E.W. Müller
PI 100	1	Mano-vacuummeter	1/4"	1.4301	MGS18.1/A-63 -1..+	111000D063-57	Moesch K. AG
PI 200	1	Manometer	1/4"	MS 58	MGS10.3/A-63/0..+1	111000D063-32	Nuova Fima
PI 220	1	Manometer	1/4"	MS 58	MGS10.3/A-63/0..+1	111000D063-32	Nuova Fima
PI 300	1	Manometer	1/4"	MS	Multi-Fix G.50-10	136001MF-06	Bachofen / Knocks
PI 301	1	Mano-vacuummeter	1/4"	MS	MGS10.3/A-63SEC/-1	111000D063-31	Nuova Fima
PI 302	1	Mano-vacuummeter	1/4"	MS	MGS10.3/A-63SEC/-1	111000D063-31	Nuova Fima
PI 800	1	Mano-vacuummeter	1/4"	1.4301	MGS18.A-63,-1.+5ba	111000D063-58	Nuova Fima
PI 804	1	Mano-vacuummeter	1/4"	MS	MGS10.3/A63SEC/-1.	111000D063-33	Nuova Fima
R 300	1	Spring-loaded check valve	1/4"	AL/ZN	Multi-Fix V.11	136001MF-10	Bachofen / Knocks
R 305	1	Spring-loaded check valve	1/2"	1.4301	Small 0,02-0,03 ba	142A20015-02	Hamberger / Cromax
R 402	1	Spring-loaded check valve	1/4"	1.4301	Small 0,02-0,03 ba	142A20008	Hamberger / Cromax
R 623	1	Butterfly flap valve	1/8"x4.0mm	MS verni.	A327927	121001-26	Bachofen AG
R 624	1	Butterfly flap valve	1/8"x4.0mm	MS verni.	A327927	121001-26	Bachofen AG
R 811/812	1	Resistance thermometers	5.0mm	1.4571	2xPT1000 5.0/120/3	5550JUM-108	Jumo
R 815	1	Resistance thermometers	4.5mm	1.4571	1xPT1000/1L-4.5-20	5550JUM-19	Jumo
R 817	1	Resistance thermometers	4.5mm	1.4571	1xPT1000/1L-4.5-20	5550JUM-19	Jumo
S 331	1	Manometric switch	1/4"	MS	901.51 111L1 +2.6b	5215-25	Pewatron

DRAW NO. : 60_300.17580
 PROJECT : Standard MST-V 5-5-9 VS1 CSD1 m. KK USA DATE : 08.01.13
 DESIGNATION 1 : P & ID 5-5-9 VS1/VS2 CSD1
 DESIGNATION 2 : MST-V DATE : 23.08.11

DESIGN.	EACH	TERM	DIMENSION	MATERIAL	TYPE	PART NO.	MANUFACTURER/SUPPLIER
S 335	1	Manometric switch	1/4"	MS verni.	901.11.122X1 +0.12	5215-30	Pewatron
S 345	1	Manometric switch	1/4"	MS	901.51 111L1 +2.6b	5215-25	Pewatron
S 346	1	Manometric switch	1/4"	MS	901.51 111L1 +4.5b	5215-28	Pewatron / Beck GmbH
V 300	1	Pressure controller	1/4"	AL	R.11 / Serie 1 / M	136001MF-09	Bachofen / Knocks
V 334	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 607	1	Solenoid valve 4-way	1/4"	AL	135MF0/24VDC/13W	121000LU002-01	Bachofen / Lucifer
V 616	1	Solenoid valve 4-way	1/4"	AL	135MF0/24VDC/13W	121000LU002-01	Bachofen / Lucifer
V 618	1	Solenoid valve 2-way	1/4"	MS	EV210A 3.5B/24V DC	121MS0DA008-01	Danfoss
V 621	1	Slanted seat valve	1/2"	RG	2000-A-013, 0-E-RGG	122RG0BU015	Bürkert
V 622	1	Slanted seat valve	1/2"	1.4408	2000-A-013, 0-E-VAG	122A40BU015-01	Bürkert
V 622.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 623	1	Slanted seat valve	1/2"	1.4581	2000-A-015, 0-E-VAS	122A40BU015	Bürkert
V 623.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 624	1	Slanted seat valve	3/4"	RG	2000-A-020, 0-E-RGG	122RG0BU020-01	Bürkert
V 624.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 625	1	Slanted seat valve	1/2"	1.4408	2000-A-013, 0-E-VAG	122A40BU015-01	Bürkert
V 625.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 626	1	Slanted seat valve	1/2"	RG	2000-A-013, 0-E-RGG	122RG0BU015	Bürkert
V 626.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 627	1	Solenoid valve 2-way	1/4"	1.4401	21L2K1T30-2/2Weg-1	121A0SA008	ODE
W4	1	Heat exchanger	1"	1.4401	SL 34-TL-BB-14 / A	04-002.7451	ETS
W5	1	Heat exchanger	3/4"	1.4401	WP-2-20 m. Bolzen	05-300.13842	Javorka Thermochnik
X 100	1	Orifice	G 3/4" zu ø 1/2"	1.4435	zum Aufbohren	750 05362	Marton / Zweifel
X 110	1	Orifice	zu ø 12.0mm	1.4435	D=14.5 d=5.0 s=1.5	750 05133	Marton / Zweifel
X 205	1	Tube	8.0mmx13.5mm	PVC-Ny.gefl.	8x14 mm PVC-transp	350 00727	Maagtechnic
X 205.1	1	Tube	8.0mmx13.5mm	PVC-Ny.gefl.	8x14 mm PVC-transp	350 00727	Maagtechnic
X 300	1	Pipe	ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 309	1	Pipe	ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 309.1	1	End plate	6.0mm	KS	PS1-E102	121001-04	Girmatic / Parker
X 309.2	1	Sound absorber	6.0mm	KS	0671 06 00CH	121001-05	Parker

DRAW NO. : 60_300.17580
 PROJECT : Standard MST-V 5-5-9 VS1 CSD1 m. KK USA
 DESIGNATION 1 : P & ID 5-5-9 VS1/VS2 CSD1
 DESIGNATION 2 : MST-V
 ORIGINATOR : IRU
 DATE : 08.01.13
 CHECKED : PBA
 DATE : 23.08.11
 REVISED : -
 DATE : -
 CHECKED : -
 DATE : -
 INDEX : -
 FILE : MSTV_010_A
 LIST TYPE: EFMTXA

DESIGN.	EACH	TERM	DIMENSION	MATERIAL	TYPE	PART NO.	MANUFACTURER/SUPPLIER
X 309.3	1	Ball	3.0mm	ST chem.verni.	HY219 488	121001-07	Girmatic / Parker
X 309.4	1	Sound absorber	6.0mm	KS	0671 06 00CH	121001-05	Parker
X 400	1	Compensator	1"	1.4404	Assiwel1 DN 25	106A40-61	Angst+Pfister
X 400.1	1	Compensator	1"	1.4404	Assiwel1 DN 25	106A40-61	Angst+Pfister
X 400.2	1	Expansion joint	3/4"	1.4404	FWH840.020.25	300 02805	Torgen GmbH
X 402	1	Pipe	ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 402.1	1	Pipe	ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 405	1	Orifice	zu ø 8.0mm	1.4435	D=10.0 d=0.6 s=1.5	350 07094	Marton / Zweifel
X 405.1	1	Pipe	ø 6.0x4.0mm	Polyamid	1100P0601 - Schwar	121001-09	Legris AG
X 500	1	Orifice	zu ø 15.0mm	1.4435	D=18.5 d=2.0 s=1.5	750 00049	Marton / Zweifel
X 600	1	Compensator	1"	1.4404	Assiwel1 DN 25	106A40-61	Angst+Pfister
X 705	1	Sound absorber	1"	Sinterbronze	06700034-1"	121001-25	Legris
Y 700	1	Safety valve	3/4"-1"	RG	0041EDA-40 / Stamp	133RG000020-32	Spence/USA
Y 702	1	Safety valve	3/4"-1"	RG	0041EDA-40 / Stamp	133RG000020-32	Spence/USA

32 Bill of Material to P & I Diagram Electric Steam Generator

page: 1 / 1

LIST OF SPARE PARTS / PROCESS TECHNOLOGY

DESIGN.	EACH	TERM	DIMENSION	MATERIAL	TYPE	PART NO.	MANUFACTURER/SUPPLIER
B 122	1	Level switch	3/8"	1.4435	L=120mm	60-300.15856	GEMS/Warrick
B 125	1	Level switch	3/8"	1.4435	L=135mm	60-300.15857	GEMS/Warrick
B 805	1	Pressure transducer	1/4"	1.4571	MIDAS S05/0-6bar/4	5510JUM-19	Jumo
F 608	1	Strainer	3/4" zu ø 1/2"	1.4404	03-000.5691	03-000.5691	Sauter Sul. /Künzli AG
H 100	1	Manual flow balancing valve	1"	316SS	054-44-6-6-TT-BW4-	123VA0W0025	Avintos AG
H 258	1	Manual flow balancing valve	1/2"	1.4401	RB-7310-015	131A4000015-02	Trigrress
H 259	1	Manual flow balancing valve	1/2"	316SS	05-44-6-6-TT-BW4-B	123VA0W0015-01	Avintos AG
H 608	1	Manual flow balancing valve	1/2"	MS verni.	B-RB-1700-015	131MS000015-05	Trigrress/Economic
H 608.1	1	Manual flow balancing valve	3/8"	MS verni.	B-RB-1700-010	131MS000010-03	Trigrress/Economic
H 609	1	Manual flow balancing valve	1/2"	316SS	05-44-6-6-TT-BW4-B	123VA0W0015-01	Avintos AG
H 808	1	Manual flow balancing valve	1/4"	MS verni.	B-RB-1720-008	131MS000008	Trigrress/Economic
P8	1	Pump		1.4408/1.4581	Y-2951 UL 0.12 kW	5140-53	Speck
PI 100	1	Mano-vacuummeter	1/4"	1.4301	MGS18.1/A-63 -1..+	111000D063-57	Moesch K. AG
PI 258	1	Manometer	1/4"	1.4304	MGS18.3/A / 0..+16	111000D063-40	Nuova Fima
R 102	1	Screwed heating unit	2"	1.4404	L=690mm/3x480V/60H	5350-32	Electrolux Professional
R 103	1	Screwed heating unit	2"	1.4404	L=690mm/3x480V/60H	5350-32	Electrolux Professional
R 259	1	Spring-loaded check valve	1/2"	1.4401	Small 0,02-0,03 ba	142A40015-01	Hamberger
S 112	1	Temperature limiter	1/2"	CU	EMF-50/U +160-200°	5410JUM-09	JUMO
S 113	1	Manometric switch	1/4" NPT	VA	K35 manual reset 6	5215-36	Danfoss
S 115	1	Manometric switch	1/4" NPT	VA	K35 auto. reset 6-	5215-37	Danfoss
S 344	1	Manometric switch	1/4"	1.4401	901.51.117R1	5215-26	Pewatron
V 631	1	Slanted seat valve	1/2"	1.4408	2000-A-013, 0-E-VAG	122A40E0015-01	Bürkert
V 631.1	1	Solenoid valve 3-way	-	KS	PS1-E21162BQ	121000-04	Girmatic / Parker
V 632	1	Solenoid valve 2-way	1/4"	1.4401	21L2K1T30-2/2weg-1	121A0SA008	ODE
V 633	1	Solenoid valve 2-way	1/4"	1.4401	21L2K1T30-2/2weg-1	121A0SA008	ODE
X 100.1	1	Tubing piece	1"-1"	PTFE/VA-Umfle.	FT800006.8519	350 07090	Angst + Pfister
X 238	1	Orifice	zu ø 8.0mm	1.4435	D=10.0 d=1.0 s=1.5	350 00682	Marton / Zweifel
X 259	1	Orifice	zu ø 12.0mm	1.4435	D=14.5 d=2.0 s=1.5	750 05293	Marton / Zweifel
Y 708	1	Safety valve	1/2"-3/4"	Bronze	0041DCA-55	133RG000015-24	Spence/USA

Graphic Symbols for Fittings

Symbol / Character	Term
	V Bellows valve piston operated
	V Bellows valve hand operated
	V Bellows valve piston operated (angle type)
	V Bellows valve piston operated (angle type)

Graphic Symbols for Fittings acc. to DIN 2429 part 2

Symbol / Character	Term
	A Steam Trap
	F Strainer
	F Gas filter / Air filter
	G Inspection glass
	H Ball valve hand operated
	H Ball valve motor driven
	H Ball valve pneumatic driven
	K Swing type check valve
	K Shutting clack hand operated
	K Shutting clack piston operated
	K Shutting clack pneumatic driven
	R Spring-loaded check valve
	V Ventilation
	V Shut-off valve hand operated
	V Diaphragm valve hand operated
	V Shut-off valve with float cock
	V Pressure reducing valve
	V Shut-off valve diaphragm driven
	V Diaphragm valve diaphragm driven
	V Shut-off valve piston operated
	V Diaphragm valve piston operated
	V Control valve with linear actuator
	V Control valve piston operated
	V 3/2 way valve hand operated
	V 3/2 way valve piston operated
	V Solenoid valve 2-way
	V Solenoid valve 3/2-way
	X Spraying nozzle
	X Orifice, flow reducer
	X Sound absorber

	X Pipe disconnecter
	X T-distribution
	Y Safety valve
	X Ball

Graphic symbols for equipment acc. to DIN 28004 part 3

Symbol / Character	Term
	B Sterilizer / Disinfection Equipment
	B Tank / Container
	D Steam Generator
	F Filter
	P Pump
	P Vacuum pump
	V Ventilator
	W Plate heat exchanger
	W Gilled radiator heat exchanger
	W Tubular heat exchanger
	X Gas radiator

Graphic symbols for pipes acc. to DIN 2429 part 2

Symbol / Character	Term
	X Tube
	X Expansion joint
	Flexible connection
	U-tube
	Pipe
	Control line
	Crossing with connection
	Crossing without connection
	Crossing without connection
	Branch-off

Graphic symbols for pneumatic acc. to DIN ISO 1219

Symbol / Character	Term
	R Flow control
	V Solenoid 5/2-way
	X Lubrication unit
	Z Pneumatic cylinder

Graphical symbols and identifying letters for measurement, open-loop control and closed-loop control in process engineering to DIN 19227 1


○ □ measuring location symbol
(Measurement/open-loop control/closed-loop control)

ⓓ ⓓ Process calculator symbol (Process computer)

The function of a measurement and control point is entered in the top half of the measurement and control point circuit using identifying letters shown in the table below. The measurement and control point number is entered in the bottom half of the measurement and control point circuit.

The measurement and control point number is formed from the identifying letters of the measurement and control functions and a consecutive number in the case of mechanical components and from the electrical-engineering designation to DIN 40719/1.74 and the sheet number and current path number in the case of electrical devices.

Example:



In the process computer point circuit, the process computer functions are entered in the top half of the process computer point circuit in accordance with the identifying letters shown in the table below.

Table for recognition character

Character	Group 1		Group 2
	Measured variable or other input variable as first character	as second character	Processing as additional character
A			State display, alarm
B			
C			Controlling
D	Density	Difference	
E	All electrical variables		Sensing element
F	Flow rate	Ratio	
G	Gauging, position or length		
H	Hand, manually initiated		
I			Indicating
J			
K	Time or time program		
L	Level		
M	Moisture or humidity		
N	User's choice		User's Choice
O	User's choice		Test point connection
P	Pressure or vacuum		
Q	Quality, e. g. analysis, concentration, conductivity	Integrate or totalise	
R	Nuclear radiation		Recording
S	Speed or frequency		Switching
T	Temperature		Transmitting, multifunction unit
U	Multi-variable		
V	Viscosity		Valve, damper, louver, actuating element
W	Weight or force		
X	Unclassified variables		Unclassified functions, e. g. TV camera, cathode ray tube, radioactive source
Y	User's choice		Computing relay, relay
Z			Emergency or safety acting
+			Upper limit value
/			Intermediate value
-			Lower limit value

Manufacturer:

Belimed Sauter AG
Zelgstrasse 8
8583 Sulgen
Switzerland

Tel. +41 71 644 85 00
Fax +41 71 644 86 00
contact@belimed.ch

Service organisations:

BELGIUM

Belimed SA - NV
Rue de Clairvaux 8
1348 Louvain-la-Neuve
Tel. +32 10 42 02 40
Fax +32 10 42 02 49
info@belimed.be

FRANCE

Belimed SAS
Parc GVIO
330 Allée des Hêtres, Hall E
69760 Limonest
Tel. +33 4 37 41 83 71
Fax +33 4 37 41 63 04
info@belimed.fr

AUSTRIA

Belimed GmbH
Grüne Lagune 1
8350 Fehring
Tel. +43 3155 40 699 0
Fax +43 3155 40 699 10
info@belimed.at

USA

Belimed Inc.
2284 Clements Ferry Road
Charleston, SC 29492
Tel. +1 843 216 7424
Fax +1 843 216 7707
info@belimed.us

CHINA

Belimed Medical Equipment
(Shanghai) Co. Ltd.
CaiLun Road #780
5th floor, Room H
ZhangJiang Hi-Tech Park
201203 Pudong, Shanghai
Tel. +86 40 082 066 05
Fax +86 21 513 709 96
info@belimed.cn

UNITED KINGDOM

Belimed Limited
Unit 4 Newbuildings Place
Dragons Green Road
Shipley
West Sussex, RH13 8GQ
Tel. +44 1403 738 811
Fax +44 1403 730 830
admin@belimed.co.uk

SWITZERLAND

Belimed Sauter AG
Zelgstrasse 8
8583 Sulgen
Tel. 0848 55 88 00
Fax +41 71 644 86 06
service.ch@belimed.com

OTHER COUNTRIES

Contact your local Belimed
distribution partner.

GERMANY

Belimed Deutschland GmbH
Edisonstrasse 7a
84453 Mühldorf am Inn
Tel. +49 8631 9896 500
Fax +49 8631 9896 560
service@belimed.de

THE NETHERLANDS

Belimed B.V.
Vlambloem 65
3068 JG Rotterdam
Tel. +31 10 286 17 50
Fax +31 10 456 56 97
service@belimed.nl

SLOVENIA

Belimed d.o.o.
Taborska cesta 38 E
1290 Grosuplje
Tel. +386 1 786 60 30
Fax +386 1 786 60 11
service@belimed.si

www.belimed.com

